Virtual®

CADbite Registration

Instructions for Use
Gebrauchsinformation
Mode d'emploi
Istruzioni d'uso
Instrucciones de Uso
Instruções de Uso
Bruksanvisning
Brugsanvisning

- Reflective vinylpolysiloxane (addition-type) bite registration material
- Reflektierendes Vinylpolysiloxan (additionsvernetzendes) Bissregistrierungsmaterial
- Matériau de type vinyl-polysiloxane possédant des propriétés réfléchissantes pour l'enregistrement d'occlusion
- VinilPoliSilossano riflettente per la registrazione del morso
- Material de registro de mordida (tipo adición) de vinilpolisiloxano reflectante.
- Material de registo de oclusão (tipo adição) de vinilpolisiloxano reflectante
- Reflekerande bettregistreringsmaterial av vinylpolysiloxan (additionstyp)
- Reflekerende vinylpolysiloxan (additionspolymeriserende) bid-registreringsmateriale

Käyttöohjeet
Bruksanvisning
Gebruiksaanwijzing
Οδηγίες Χρήσεως
Kullanım kılavuzu
Инструкция по применению
Instrukcja stosowania

- Heijastava vinyylipolysiloksaani (A-silikoni-perusteinen) purennan rekisteröintimateriaali
- Reflekerende bittregistreringsmateriale av vinylpolysiloksan (addisjonstype)
- Reflecterend, additiereactief vinylpolysiloxaan-afdrukmetaal
- Ανακλαστικό υλικό αποτύπωσης σύγκλεισης από βινυλοπολυσιλοξάνη (τύπου προσθήκης)
- Reflektif vinilpolisiloksan (ilave tip) karanis kaydi maddesi
- Рефлектирующий материал на основе поливинилсилоксана для регистрации прикуса
- Odbijający światło materiał poliwinylsiloksanowy (silikon typu addycyjnego), przeznaczony do rejestracji zwarcia

For dental use only!

Caution: U.S. Federal law restricts this device to sale by or on the order of a licensed dentist.

608387/0811/WE3

ivoclar
vivadent
clinical

English

Instructions for Use

Description

Virtual CADbite Registration is an addition-reaction silicone developed specifically to record the occlusal bite for the fabrication of indirect restorations. The material offers a thixotropic viscosity, high precision, a fast setting time with a final hardness of 32 Shore D and is dimensionally stable.

Virtual CADbite Registration is reflective and may be used for the optical recording of data when fabricating restorations using CAD/CAM/CIM technology, e.g. the CEREC system (antagonist bite registration).

Colors and Working Time

See table "Technical Data"

Technical Data

Color	Blue
Mixing Ratio [Base:Catalyst]	1:1
Working time after mixing	max. 30 sec
Set Time in Mouth	min. 45 sec
Detail Reproduction	2 µm
Linear Dimensional Change [24h]	< 0.05 %
Hardness Shore after Polymerization	32 Shore D

Composition

Virtual CADbite Registration is an addition-reaction silicone containing vinylpolysiloxane, methylhydrogensiloxane, organoplatinic complex, silica, pigment and food dyes.

Indication

- Virtual CADbite Registration is recommended for use by dental professionals to assist them in creating proper articulation of indirect restorations.
- Optical recording of the antagonist situation for CAD/CAM/CIM systems

Contraindication

The use of Virtual CADbite Registration material is contraindicated if the patient is known to be allergic to any of its ingredients.

Side Effects

Not known

Interactions

Latex gloves may inhibit the setting of Virtual CADbite Registration.

Do not touch tooth surfaces, preparations or retraction cords with latex gloves. Materials such as rubber dam or certain substances may also inhibit the setting reaction of vinylpolysiloxane. If the operator suspects that the preparation has been contaminated, it is recommended that the preparation be rinsed and dried.

Application

- Loading the Cartridge Dispenser

Fig. 1

Fig. 2

Fig. 3

Fig. 4

1. Press the black release lever located below the plunger on the back of the dispenser and pull the plunger as far back as possible (Fig. 1).
 2. Lift the cartridge lock and insert the cartridge with the "V" shape on the cartridge base turned downwards. Lower the cartridge lock (Fig. 2).
 3. Remove cartridge cap by turning it 1/4 turn counter clockwise (Fig. 3). Discard the cap.
 4. It is VERY IMPORTANT to clear or bleed the cartridge PRIOR to applying the mixing tip! This will assure a proper mix and set time. Gently press on the dispenser lever until both components (base and catalyst) begin to flow out of the cartridge evenly onto a mixing pad.
 5. Insert the mixing tip in a clean motion, so that the "V" shaped notch on the coloured base is aligned with the "V" shaped notch on the housing of the cartridge. While gripping the coloured base, and not the mixing tip, turn the mixing tip base 1/4 turn clockwise (Fig. 4).
 6. Begin mixing by pressing gently on the lever of the dispenser.
- **Conventional bite registration**
 1. Dispense bite registration material directly onto the occlusal surface of the lower teeth.
 2. The patients should be asked to close into their normal occlusal relationship and hold the position until the material is completely set (approx. 45 seconds).
 - **Optical data recording**
 1. Dispense Virtual CADbite Registration directly onto the occlusal surface of the prepared teeth.
 2. The patient should be asked to close their teeth together without any functional movement and hold the position until the material is completely set (approx. 45 seconds).
 3. If required, the bite registration can be removed and trimmed to enable data recording.
 4. Then reseat the bite registration on the preparation and take an image of the antagonist situation.

– Removing the cartridge

1. Lift the release lever to pull the slide back. Lift the cartridge lock and remove the cartridge.
2. Leave the mixing tip on the cartridge. The material will polymerize within the mixing tip, functioning as a natural cap.

– Disinfection

Bite registrations made with Virtual impression materials can be immersed in disinfection solution (glutaraldehyde 0.5%, benzalkonium chloride 0.5%) for 10 minutes without distortion.

Special Note

Virtual CADbite Registration should be processed at room temperature (23 °C). Lower temperatures, e.g. when stored in the refrigerator, will prolong the working time and the time the impression material should remain in the mouth, while higher temperatures will reduce it.

Vinylpolysiloxanes are chemically resistant.

Unpolymerized material may stain clothing.

Warnings

If uncured material comes into contact with the eyes, rinse immediately with plenty of water. If irritation persists, consult a physician. If the material comes into contact with the skin, rinse with copious amounts of water.

Storage

- Storage temperature: 2 - 28 °C / 36 - 82 °F
- Expiration date: see label on primary packaging
- Keep away from direct heat sources!

Keep out of the reach of children!

For use in dentistry only!

Date information prepared: 08/2011; Rev. 2

Ivoclar Vivadent AG, FL-9494 Schaan/Liechtenstein

These materials have been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes.

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Ivoclar Vivadent Marketing Ltd.
Ivoclar Vivadent AB
Ivoclar Vivadent Limited
Ivoclar Vivadent, Inc.