Instructions for Use

Page 2

Gebrauchsinformation

Seite 6

Mode d'emploi

Page 10

Istruzioni d'uso

Pagina 14

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Pagina 18

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Página 22

Bruksanvisning

Sidan 26

Brugsanvisning

Side 30

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Sivu 34

Bruksanvisning

Side 38

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Pagina 42

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ivoclar
vivadent
technical

Indication

The phosphate-bonded IPS PressVEST investment material (Type 1, Class 2) can be used with the following press ceramics in the Ivoclar Vivadent Programat press furnaces.

- IPS e.max® Press
- IPS e.max® ZirPress
- IPS Empress® Esthetic
- IPS Empress® Cosmo
- IPS InLine® PoM (Press-on Metal-Ceramics)

Contraindication

All other applications not listed in the indications.

Delivery form

- IPS PressVEST Powder, 25 bags, 100 g each, including 1 measuring cup
- IPS PressVEST Liquid, 0.5 litre
- IPS PressVEST Powder, 50 bags, 100 g each, 1 measuring cup
- IPS PressVEST Liquid, 1.0 litre

Storage

- Storage at 12–28 °C / 54–82 °F
- Store powder in a dry place.
- The liquid must not be exposed to temperatures below +5 °C / 41 °F (frost susceptible!). **If frozen once, the liquid is unusable (crystallization) and must not be used anymore.**

Processing temperatures

- Room temperature 18 °C – max. 23 °C / 64 °F – max. 73 °F
- Any other processing temperature decisively affects the setting behaviour.
- A higher processing temperature accelerates the setting behaviour.
- A lower processing temperature delays the setting behaviour.

Investment ring system and press plunger

- For IPS e.max and IPS InLine PoM ingots, it is imperative to use the IPS Investment Ring System (100 g, 200 g and 300 g).
- For IPS e.max and IPS InLine PoM ingots, it is imperative to use the IPS AloX plunger.
- For IPS Empress ingots, it is imperative to use the IPS Empress Investment Ring System (100 g, 200 g).
- For the IPS Empress Esthetic ingots, the IPS Empress AloX plunger or the IPS One-Way-Plunger is used.
- Use the IPS Silicone Ring for investment. The IPS Silicone Ring is used for both, the IPS Investment Ring System and the IPS Empress Investment Ring System.

Mixing ratio 100 g powder : 22 ml diluted liquid

IPS e.max® Press				
Indication	Concentration approx.	100-g investment ring Liquid : dist. water	200-g investment ring Liquid : dist. water	300-g investment ring Liquid : dist. water
(Thin) veneers, partial crowns, single crowns	60–70%	13 ml : 9 ml 15.5 ml : 6.5 ml	26 ml : 18 ml 31 ml : 13 ml	---
Hybrid abutments, hybrid abutment crowns	70%	16 ml : 6 ml	32 ml : 12 ml	---
Inlays	50%	11 ml : 11ml	22 ml : 22 ml	---
3-unit bridges in the anterior and premolar region	60–70%	---	26 ml : 18 ml 31 ml : 13 ml	---
Pressing over electroplated frameworks	60–70%	13 ml : 9 ml 15.5 ml : 6.5 ml	26 ml : 18 ml 31 ml : 13 ml	---

IPS e.max® ZirPress				
Indication	Concentration approx.	100-g investment ring Liquid : dist. water	200-g investment ring Liquid : dist. water	300-g investment ring Liquid : dist. water
Pressing over ZrO ₂ single tooth constructions	60–70%	13 ml : 9 ml 15.5 ml : 6.5 ml	26 ml : 18 ml 31 ml : 13 ml	39 ml : 27 ml
Pressing over ZrO ₂ frameworks bridge constructions	60–70%	---	26 ml : 18 ml 31 ml : 13 ml	39 ml : 27 ml

IPS Empress® Esthetic				
Indication	Concentration approx.	100-g investment ring Liquid : dist. water	200-g investment ring Liquid : dist. water	300-g investment ring Liquid : dist. water
Crowns, veneers	60–70%	13 ml : 9 ml 15.5 ml : 6.5 ml	26 ml : 18 ml 31 ml : 13 ml	---
Inlay MOD	50–60%	11 ml : 11 ml 13 ml : 9 ml	22 ml : 22 ml 26 ml : 18 ml	---
Inlay MOD + OD	50–60%	11 ml : 11 ml 13 ml : 9 ml	22 ml : 22 ml 26 ml : 18 ml	---

IPS Empress® Cosmo				
Indication	Concentration approx.	100-g investment ring Liquid : dist. water	200-g investment ring Liquid : dist. water	300-g investment ring Liquid : dist. water
Preprosthetic buildup in combination with CosmoPost	50–60%	11 ml : 11 ml 13 ml : 9 ml	22 ml : 22 ml 26 ml : 18 ml	---

IPS InLine® Pom – Press-on-metal ceramic				
Indication	Concentration approx.	100-g investment ring Liquid : dist. water	200-g investment ring Liquid : dist. water	300-g investment ring Liquid : dist. water
Pressing over opaquerized metal frameworks	60%	13 ml : 9 ml	26 ml : 18 ml	39 ml : 27 ml

The total quantity of liquid (liquid + dist. water) must not be altered. The indications regarding the liquid concentration are approximates and may be adjusted, if required. Do not dilute the liquid to below 50%.

Processing

Excess of separating agents may cause reactions with the investment material and affect the press result. Disperse excess material with compressed air before investment.

- Fill liquid into the mixing cup.
- Add powder.
- Thoroughly mix investment material with a spatula for 20 seconds until even wetting of the powder is achieved.

Vacuum mixer

Vacuum mixing time (at room temperature)	1 min
Speed	approx. 350 rpm
Processing time	6 min
Setting time	60 min

The processing time depends on the material temperature, mixing quantity, mixing time and mixing intensity. A higher material temperature and/or a longer mixing time shortens the processing time.

Investment

- Mix max. 400 g investment material in one go.
- Closely observe mixing ratio!
- Investment is carried out on a shaker under slight vibration. Allow the investment ring to set without manipulating it.

Some aspects influencing the setting expansion of the investment material

- Temperature of the investment material and the liquid
- Residual water in the mixing cup, atmospheric humidity

Preheating / holding times

Starting temperature	Room temperature
Temperature increase	5 °C / 41 °F per min
Holding time	30 min at 250 °C / 482 °F
Holding time after reaching the final temperature (850 °C / 1562 °F)	100-g investment ring – min. 45 min
	200-g investment ring – min. 60 min
	300-g investment ring – min. 90 min

If more than 2 investment rings are preheated at the same time, the holding time at the final temperature has to be prolonged by 15 minutes per additional investment ring.

Preheating of ingots and Alox plungers

IPS e.max® Press / IPS e.max® ZirPress / IPS InLine® PoM	
Ingot	No preheating
Alox plunger	No preheating
IPS Empress® Esthetic	
Ingot	Preheating
Alox plunger	Preheating
IPS One-Way-Plunger	No preheating
IPS Empress® Cosmo	
Ingot	No preheating
Alox plunger	Preheating

Important information

- Do not use a debubbler (tenside). Reactions with the investment material may occur.
- Consistent results can only be achieved under consistent conditions and if the Instructions for Use are observed.
- Dilute liquid exclusively with distilled or de-ionized water.
- Liquid concentration: The more distilled or de-ionized water is added to the mixing liquid, the lower the setting expansion. Do not dilute the liquid to below 50%.
- Use only dry, clean instruments to process the investment material. Do not use plaster mixing cups.
- Carefully cut open the powder bags and remove the entire contents of the bags.
- The stipulated minimum setting and holding times must be closely observed.
- Do not invest under pressure since this changes the product properties.
- Do not invest over the weekend (risk of cracks and press flashes).

To ensure thorough heating of the investment rings in the preheating furnace, the following points have to be observed:

- Always place the investment rings into the cool preheating furnace.
- Load the preheating furnace to maximally half of the available utility space.
- Always place the investment rings in the rear part of the preheating furnace.

- Always place the investment rings in the preheating furnace with the opening facing down and at an angle of 45°, if possible.
- Do not place the investment rings in direct contact with each other. Allow air circulation.
- Always remove pressed investment rings from the press furnace immediately after the end of the program and place them on a grid for cooling.
- When divesting, ensure to only sandblast (4 bar pressure) the excess of the investment material and not the ceramic object.


IMPORTANT! Safety notes

- Do not open the furnace during the first 20 minutes of preheating: Burn hazard caused by wax fumes that might ignite in the air!
- This material contains quartz, which may cause lung diseases (silicosis or cancer) when inhaled. Prevent the inhalation of dust and wear a suitable dust mask.

Physical data

(according to ISO 15912:2006)

- | | | |
|-----------------------------|--------------|--------------------------|
| – Flow characteristics: | 152 mm | Liquid concentration 60% |
| – Initial set: | 8.15 minutes | Liquid concentration 60% |
| – Compressive strength: | 3.4 MPa | Liquid concentration 60% |
| – Linear thermal expansion: | 0.97 % | Liquid concentration 60% |

This material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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Ivoclar Vivadent Limited
Ivoclar Vivadent, Inc.