ProBase[®] Cold

Instructions for Use
Verarbeitungsanleitung
Mode d'emploi
Istruzioni d'uso
Instrucciones de uso
Instruções de Uso
Bruksanvisning
Brugsanvisning
Käyttöohjeet
Bruksanvisning
Productinformatie
Οδηγίες Χρήσεως
Kullanma Talimatı
Инструкция по применению
Instrukcja stosowania

CE 0123

Complies with / entspricht:
ISO 20795-1; EN ISO 20795-1

For dental use only.
Rx only

Manufacturer:
Ivoclar Vivadent AG, FL-9494 Schaan/Liechtenstein
www.ivoclarvivadent.com

ivoclar
vivadent[®]
technical

english

Description

ProBase® Cold is a self-curing denture base material. It demonstrates excellent flow and modelling properties. It is easy and reliable to use with the pouring or packing technique, even when two or more saddles are present. The material is available in a variety of shades. As the powder and liquid can be dosed as desired within the usual limits, users can vary the consistency and working time of ProBase Cold.

Composition

Powder

Polymethyl methacrylate, softening agent, benzoyl peroxide, catalyst, pigments

Liquid

Methyl methacrylate, dimethacrylate, catalyst

Indication

- Partial dentures
- Combination dentures
- Relining
- Repairs
- Complete dentures

Contraindication

- Direct intraoral contact of unpolymerized material.
- If the patient is known to be allergic to any of the ingredients in ProBase Cold

Side effects

In individual cases, local allergic reactions to polymethyl methacrylate materials have been reported.

Application

Pouring technique

Preparation

Isolate boiled-out, well-wetted plaster surfaces with two layers of Ivoclar Vivadent Separating Fluid. Allow to dry, well roughen the teeth, apply mechanical retention, and wet with monomer to ensure an adequate bond with the denture base.

Dosage

- **Ideal mixing ratio**
15 g polymer (powder) : 10 ml monomer (liquid)

Mixing

Thoroughly mix polymer and monomer with the spatula. Subsequently, allow the mixture to rest for 15 seconds to permit any trapped air to rise.

Flow phase

At room temperature (23 °C / 73 °F), the flow phase is approx. 2.5 to 3 minutes. Pour the material into the saddle within this time span.

Modelling phase

The material is set after a transition period of approx. 4 minutes. It can be modelled during an additional 3 minutes.

High room temperatures shorten the flow and modelling phase.

Polymerization

Polymerization is carried out in a pressure device (e.g. Ivomat) at 40 °C / 104 °F and at 2 to 6 bar pressure for 15 minutes.

Finishing

Remove the precast, finish and polish in the usual manner.

Packing technique

Preparation

Isolate boiled-out, well-wetted plaster surfaces with two layers of Separating Fluid and allow to dry completely. Well roughen the teeth and wet with monomer to ensure an adequate bond with the denture base.

Dosage

Ideal mixing ratio for one denture: 20.5 g polymer (= 1st graduation on measuring cylinder): 10 ml monomer

Integrated dosage system

This system ensures an ideal mixing ratio and, therefore, minimum polymerization shrinkage of ProBase Cold. The measuring cylinder for the polymer indicates the quantity of material required for one or two medium-sized dentures. The graduation on the measuring cylinder for the monomer is in millilitres.

Mixing

Thoroughly mix polymer and monomer in the given mixing ratio with a spatula. Mix thoroughly. Allow the dough to mature in a closed mixing cup for approx. 3 to 4 minutes. Subsequently, work the dough within approx. 2 minutes.

High room temperatures shorten the working time.

Pressing

Place a sufficient quantity of the resin dough in the hand-warm, wetted and isolated flask halves. Carefully close flask and load with 80 bar pressure. Fix with a clamp.

Polymerization

Polymerization is carried out by means of the clamp or a pressing device under constant pressure for 30 minutes (if the room temperature is 23 °C / 73 °F).

Deflasking and finishing

Open the flask and remove plaster. Check occlusion of the denture. Subsequently, finish and polish in the usual manner.

Repair and correction possibilities with ProBase Cold

- Corrections and repairs of ProBase Hot, ProBase Cold and SR Ivocap® may be carried out with ProBase Cold by using the pouring technique. Thoroughly roughen the corresponding surfaces and wet with monomer.
- The residual monomer content after polymerizing the material according to the method described is < 4,5%.

Warnings

- The monomer contains methyl methacrylate (MMA).
- Methyl methacrylate is easily flammable and irritating (flash point +10 °C / 50 °F).
- Irritating to eyes, skin, and respiratory system.
- May cause sensitization by skin contact.
- Avoid contact of the skin with monomer or uncured material. Commercial medical gloves do not provide protection against the sensitizing effect of methacrylates.
- Do not breathe vapour.
- Keep away from sources of ignition – no smoking.
- Do not empty into drains.
- Take precautionary measures against static discharges.

Storage

- Store material in a cool, dark, well-ventilated place. Storage temperature: 2–28 °C (36–82 °F).
- Do not use the materials after the indicated date of expiration.
- Keep out of the reach of children.

Date information prepared: 08/2012

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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Pulver:

Polymethylmethacrylat, Weichmacher, Benzoylperoxid, Katalysator, Pigmente

Flüssigkeit:

Methylmethacrylat, Dimethacrylat, Katalysator

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Ivoclar Vivadent Limited
Ivoclar Vivadent, Inc.