SR Triplex® Hot

Instructions for Use
Verarbeitungsanleitung
Mode d'emploi
Istruzioni d'uso
Instrucciones de uso
Instruções de Uso

Bruksanvisning
Brugsanvisning
Käyttöohjeet
Bruksanvisning
Productinformatie
Οδηγίες Χρήσεως

CE 0123

For dental use only.
Caution: U.S. Federal Law restricts this device to
sale by or on the order of a licensed dentist.

Made in Liechtenstein
Ivoclar Vivadent AG, FL-9494 Schaan / Liechtenstein

The logo features a series of small dots forming a semi-circular arc above the text. The text 'ivoclar' is in a lowercase, sans-serif font, and 'vivadent' is in a larger, bold, lowercase, sans-serif font. Below this, the word 'technical' is written in a smaller, italicized, lowercase, sans-serif font.
ivoclar
vivadent®
technical

Instructions for Use

Description

SR Triplex Hot is a heat-curing, PMMA-based denture base material for the condensing technique. It offers you ease of handling and your patients comfort of fit. Additionally, a veined version of the material is available.

Composition

Powder

Polymethyl methacrylate, Catalyst, Pigments

Liquid

Methyl methacrylate stab., dimethacrylate

Indication

For use in the condensing technique for

- Complete dentures
- Partial dentures
- Combination dentures
- Relining

Contraindication

- Direct intraoral contact of unpolymerized material.
- If the patient is known to be allergic to any of the ingredients in SR Triplex Hot

Side effects

In individual cases, local allergic reactions to polymethyl methacrylate materials have been reported.

Application procedure

Isolate boiled-out, well-wetted plaster surfaces with two layers of Ivoclar Vivadent Separating Fluid and allow to dry. To ensure an adequate bond with the denture base, roughen the teeth thoroughly and wet them with monomer.

Dosage

- **Ideal mixing ratio for one denture**
23.4 g polymer (powder) : 10 ml monomer (liquid)
- **With dosage system**
 - 1 denture**
Polymer: 2nd graduation on measuring cylinder
Monomer: 10 ml
 - 2 dentures**
Polymer: 3rd graduation on measuring cylinder
Monomer: 20 ml

Integrated dosage system

The integrated dosage system ensures an ideal mixing ratio and minimum polymerization shrinkage of SR Triplex Hot. The measuring cylinder for the polymer indicates the

quantity of material required for either one or two medium-sized dentures. The graduation of the monomer cylinder is in millilitre.

For ideal dosage, observe the indicated quantities. Too high a monomer content alters the shade and handling properties of the material. Additionally, polymerization is impaired. This may result in irritation of the mucous membrane.

Mixing

Thoroughly mix polymer and monomer in the given ratio with a spatula. Leave to mature in the closed mixing cup at room temperature (23 °C / 73 °F) for approx. 10 minutes.

Working time

As soon as the material has matured sufficiently and is no longer sticky, it can be worked for approx. 20 minutes at 23 °C / 73 °F.

Dough time and working time depend on the room temperature.

Pressing

Place a sufficient quantity of resin dough in the hand-warm (approx. 40 °C / 104 °F), isolated flask halves. Carefully close flask, load with 80 bar pressure and fix with a clamp. **Maintain pressure.**

Polymerization

Heat-polymerization can be carried out in two ways:

- **Standard procedure (recommended method)**
Place closed flask in cold water, heat up to 100 °C / 212 °F and let boil for 45 minutes.
- **Alternative procedure**
Place flask in boiling water. Bring the water to the boil again and then let boil for 20 minutes.

The alternative procedure is suitable only for small to medium-sized dentures. Layer thicknesses must not exceed 1 cm.

The residual monomer content can be reduced by increasing the polymerization temperature and extending the polymerization time. We recommend using the standard procedure to keep residual monomer at minimum levels. If the standard procedure is used, the residual monomer content is < 2.2%.

Cooling

Let the flask cool at room temperature for 30 minutes. Subsequently, completely cool the flask with cold water.

Deflasking and finishing

Open the completely cooled flask and remove plaster.

Check occlusion of the denture. Subsequently, finish and polish in the usual manner.

Repairs and correction possibilities of SR Triplex Hot

Corrections and repairs can be carried out with cold-curing SR Triplex Cold or ProBase Cold. Thoroughly roughen the corresponding surfaces with a bur and wet with monomer. Process the repair material according to the pouring technique.

Warnings

- The monomer contains methyl methacrylate (MMA).
- MMA is easily flammable and irritating (flash point +10 °C / 50 °F).
- Irritating to eyes, skin, and respiratory system.
- May cause sensitization by skin contact.
- Avoid contact of the skin with monomer or uncured material. Commercial medical gloves do not provide protection against the sensitizing effect of methacrylates.
- Do not breathe vapour.
- Keep away from sources of ignition – no smoking.
- Do not empty into drains.
- Take precautionary measures against static discharges.

Storage

- Store the material in a cool, dark, well-ventilated place. Storage temperature: 2 to 28 °C / 36 to 82 °F.
- Do not use the material after the date of expiration.
- Keep out of the reach of children.

Date information prepared: 11/2010

Manufacturer:

Ivoclar Vivadent AG, FL-9494 Schaan / Liechtenstein

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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Pulver:

Polymethylmethacrylat, Katalysator, Pigmente

Flüssigkeit:

Methylmethacrylat stab., Dimethacrylat

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Ivoclar Vivadent Limited
Ivoclar Vivadent, Inc.