Indications Recommended for inlays, onlays, 3/4 crowns, telescope crowns, posts, short and long span bridges, PFM crowns.

Table with 10 columns: Au, Ag, Cu, Sn, Zn, In, Ni, Pd, Pt, Ni. Values: 84.0, 7.1, 5.7, 1.0, 0.8, 0.8, 0.8, 0.8, 0.8, 0.8.

Composition

Cast Temperature: 1260-1320C/2300-2410F

METAL PREPARATION Carefully divest and clean the object with Al2O3. Do not use a hammer for divesting the object to prevent deformation.

OXIDATION Blast surface with 50-100 micron Al2O3 at max. 4.5 bar/65 psi pressure before oxidation.

HEAT TREATMENT Hardening: 510C/950F for 30 minutes; air cool.

SOLDERS AND FLUXES Design the soldering patty as small as possible and preheat it in the furnace at approximately 600C/1112F.

POLISHING After glazing remove oxide and flux residue and finish/polish the framework with rubber finishers/polishers.

INDICATIONS Recommended for inlays, onlays, 3/4 crowns, crowns, telescope crowns, conus crowns, posts, short and long span bridges, PFM crowns.

CONTRAINDICATIONS For patients with known allergy/sensitivity to any major or minor constituents of this alloy, consult with a physician is recommended.

SIDE EFFECTS In individual cases, sensitivity or allergies to components of this alloy may occur.

INTERACTIONS Galvanic effects may occur between different alloys in the same oral environment.

For additional information look into the alloy property chart.

FRANCE Ivoclar Vivadent SAS, B.P. 118, 74410 Saint-Jorioz, Tel. +33 450 88 64 00

GERMANY Ivoclar Vivadent GmbH, Dr. Adolf-Schneider-Str. 2, 73479 Ellwangen, Jagst, Tel. +49 71 61 6322 1657

INDIA Ivoclar Vivadent Marketing (India) Pvt. Ltd., 502/504, Rajhaja Plaza, 15 B Shah Industrial Estate, VERA DESAI ROAD, ANDHERI (WEST), Mumbai, 400 053, Tel. +91 22 2673 0302

CHINA Ivoclar Vivadent (Shanghai) Trading Co., Ltd., 2/F Building 1, 881 Wuding Road, Jing An District, Shanghai 200040, Tel. +86 21 6032 1657

COLOMBIA Ivoclar Vivadent Marketing Ltd., Calle 134 No. 7-8-83, Of 520 Bogotá, Tel. +57 1 613 3365

ITALY Ivoclar Vivadent s.r.l., Via Isozno 67/69, 40033 Casalecchio di Reno (BO), Tel. +39 051 6113555

JAPAN Ivoclar Vivadent K.K., 1-28-24-4F Hongo, Bunkyo-ku, Tokyo 113-0033, Tel. +81 3 6903 3535

KOREA Ivoclar Vivadent Ltd., 12 Omega Street, Rosedale, PO Box 303011 North Harbour, Auckland 0751, Tel. +64 9 914 9999

NETHERLANDS Ivoclar Vivadent B.V., De Warendy Business Park, Leicesters LE19 4GE, Tel. +44 116 284 78 80

RUSSIA Ivoclar Vivadent Marketing Ltd., Prospekt Andropova 18, Korpus 6, Office 10-06, 115422 Moscow, Tel. +7 499 418 0300

USA Ivoclar Vivadent Inc., 175 Pineview Drive, Amherst, NY 14228, Tel. +1 800 533 6825

UK Ivoclar Vivadent Limited, Compass Building, Filderspar Close, Warens Business Park, Endebey, Leicesters LE19 4GE, Tel. +44 116 284 78 80

ISRAEL Ivoclar Vivadent Ltd., 116-284-78-80, Tel. +972 3 623 54 89

USA Ivoclar Vivadent Inc., 175 Pineview Drive, Amherst, NY 14228, Tel. +1 800 533 6825

EN INSTRUCTIONS FOR USE

MODELLATION Design the framework in a reduced anatomic shape taking the planned veneer into consideration. Single crowns require a thickness of minimum 0.3 mm; abutment crowns thickness a minimum of 0.5 mm.

SPRUNG Provide the modeled bridge framework or coping with sprues of a suitable size. Use the direct or indirect technique being sure that the reservoir is positioned in the heat center.

INVESTMENT Weigh the wax pattern including the sprue to determine the quantity of the alloy to be used. (See wax conversion sheet/formula: weight x density = gr. of alloy).

BURN-OUT The suggested burnout temperature: 750-820C/1380-1510F

MELTING AND CASTING Use a separate carbon/ceramic crucible for each alloy and preheat the crucible (ceramic) in the burnout furnace. Used and new alloys must be in a ratio of 1:1.

METAL PREPARATION Carefully divest and clean the object with Al2O3. Do not use a hammer for divesting the object to prevent deformation.

OXIDATION Blast surface with 50-100 micron Al2O3 at max. 4.5 bar/65 psi pressure before oxidation. Subsequently, clean with ultrasonic or steam.

HEAT TREATMENT Hardening: 510C/950F for 30 minutes; air cool.

SOLDERS AND FLUXES Design the soldering patty as small as possible and preheat it in the furnace at approximately 600C/1112F.

POLISHING After glazing remove oxide and flux residue and finish/polish the framework with rubber finishers/polishers.

INDICATIONS Recommended for inlays, onlays, 3/4 crowns, crowns, telescope crowns, conus crowns, posts, short and long span bridges, PFM crowns.

CONTRAINDICATIONS For patients with known allergy/sensitivity to any major or minor constituents of this alloy, consult with a physician is recommended.

SIDE EFFECTS In individual cases, sensitivity or allergies to components of this alloy may occur.

INTERACTIONS Galvanic effects may occur between different alloys in the same oral environment.

For additional information look into the alloy property chart.

FRANCE Ivoclar Vivadent SAS, B.P. 118, 74410 Saint-Jorioz, Tel. +33 450 88 64 00

GERMANY Ivoclar Vivadent GmbH, Dr. Adolf-Schneider-Str. 2, 73479 Ellwangen, Jagst, Tel. +49 71 61 6322 1657

INDIA Ivoclar Vivadent Marketing (India) Pvt. Ltd., 502/504, Rajhaja Plaza, 15 B Shah Industrial Estate, VERA DESAI ROAD, ANDHERI (WEST), Mumbai, 400 053, Tel. +91 22 2673 0302

CHINA Ivoclar Vivadent (Shanghai) Trading Co., Ltd., 2/F Building 1, 881 Wuding Road, Jing An District, Shanghai 200040, Tel. +86 21 6032 1657

COLOMBIA Ivoclar Vivadent Marketing Ltd., Calle 134 No. 7-8-83, Of 520 Bogotá, Tel. +57 1 613 3365

ITALY Ivoclar Vivadent s.r.l., Via Isozno 67/69, 40033 Casalecchio di Reno (BO), Tel. +39 051 6113555

JAPAN Ivoclar Vivadent K.K., 1-28-24-4F Hongo, Bunkyo-ku, Tokyo 113-0033, Tel. +81 3 6903 3535

KOREA Ivoclar Vivadent Ltd., 12 Omega Street, Rosedale, PO Box 303011 North Harbour, Auckland 0751, Tel. +64 9 914 9999

NETHERLANDS Ivoclar Vivadent B.V., De Warendy Business Park, Leicesters LE19 4GE, Tel. +44 116 284 78 80

RUSSIA Ivoclar Vivadent Marketing Ltd., Prospekt Andropova 18, Korpus 6, Office 10-06, 115422 Moscow, Tel. +7 499 418 0300

USA Ivoclar Vivadent Inc., 175 Pineview Drive, Amherst, NY 14228, Tel. +1 800 533 6825

UK Ivoclar Vivadent Limited, Compass Building, Filderspar Close, Warens Business Park, Endebey, Leicesters LE19 4GE, Tel. +44 116 284 78 80

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