ivoclar vivadent

ISO 13485
Quality Management System Certified
Made in U.S.A.

www.ivoclarvivadent.com

Rev. Date: /_17
Rev. 7

Rx ONLY
For dental use only!

passion plus innovation

Indications: Recommended for: onlays, 3/4 crowns, short span bridges, telescope crowns, conus crowns, posts, long span bridges, PFM crowns, implant superstructure.

Composition: Au 59.4, Zr 2.0, Pd 9.5, Ag 25.5, Cu <1.0, Ni <1.0, Hf <1.0, Pt 1.0, Other <1.0

EN INSTRUCTIONS FOR USE

MODELLATION
Design the framework in a reduced anatomic shape taking the planned veneer into consideration. For non ceramic veneering material use mechanical retentions. Single crowns require a thickness of minimum 0.3 mm; abutment crowns thickness a minimum of 0.5 mm. Make sure the framework demonstrates adequate stability of shape. Avoid sharp angles. Design the connector areas to be adequate for the position and alloy being used.

SPRUNG
Provide the modeled bridge framework or coping with sprues of a suitable size. Use the direct or indirect technique being sure that the reservoir is positioned in the heat center. The connection sprues between the reservoir and the coping should be 2.5-3.0 mm in length and width.

INVESTMENT
Weigh the wax pattern including the sprue to determine the quantity of the alloy to be used. (See wax conversion sheet/formula: weight x density = gr. of alloy). Use investment following the manufacturer's instruction.

BURN-OUT
The suggested burnout temperature: 650-760C/1200-1400F

MELTING AND CASTING
Use a separate carbon/ceramic crucible for each alloy and preheat the crucible (ceramic) in the burnout furnace. Used and new alloys must be in a ratio of 1:1. Depending on the type of casting machine, follow the manufacturer's instructions for use. If you use Ivoclar Vivadent Magic Wand, set propane to 0.35 bar/5 psi and oxygen to 0.7 bar/10 psi pressure. Keep the alloy in the reducing atmosphere of the flame between the inner and outer cones. Do not use flux. After casting bench cool to room temperature.

CASTING Temperature: 1155-1215C/2110-2220F

METAL PREPARATION
Carefully divest and clean the object with AL₂O₃. Do not use a hammer for divesting the object to prevent deformation. Finish the framework with carbide burs or with ceramic-bonded grinding instruments. For a veneering with a non ceramic material the framework must be finished and polished. Treat the application surface according to manufacturers instructions. Avoid inhalation of dust during grinding!

OXIDATION
Blast surface with 50-100 micron AL₂O₃ at max. 4.5 bar/65 psi pressure before oxidation. Subsequently, clean with ultrasonic or steam. Place the object on the firing tray and provide adequate support. Place the tray in a porcelain furnace set at a low temperature of 400C/750F and increase the furnace temperature to 700C/1290F under full vacuum with 3 min. hold time at top temperature.

HEAT TREATMENT
Hardening: 450C/840F for 30 minutes; air cool.

SOLDERS AND FLUXES
Design the soldering patty as small as possible and preheat it in the furnace at approximately 600C/1121F. The soldering gap should be the same thickness as the soldering strip. Allow the object to cool slowly after soldering.

Pre-Solder: Universal Solder PKF
Post-Solder: 585 Fine Solder
Laser weld material: Laser C&B Yellow

POLISHING
After glazing remove oxide and flux residue and finish/polish the framework with rubber finishers/polishers.

INDICATIONS
Recommended for: onlays, 3/4 crowns, crowns, short span bridges, telescope crowns, conus crowns, posts, long span bridges, PFM crowns, implant superstructure.

CONTRAINDICATIONS
For patients with known allergy/sensitivity to any major or minor constituents of this alloy, consultation with a physician is recommended.

SIDE EFFECTS
In individual cases, sensitivity or allergies to components of this alloy may occur.

INTERACTIONS
Galvanic effects may occur between different alloys in the same oral environment.
For additional information look into the alloy property chart.

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