Onlays, 3/4 crowns, telescope crowns, conus crowns, short and long span bridges, posts, implant superstructure.

Table with 10 columns: Au, Pt, Pd, Ag, Cu, Ni, Zn, In, Ga, Fe. Values: 50.0, 1.0, 1.0, 1.0, 1.0, 1.0, 1.0, 1.0, 1.0, 1.0.

Composition

Instructions for Use, Gebrauchsinformation, Istruzioni d'uso, Instrucciones de Uso, Brüksanvisning, Brugsanvisning, Käyttöohjeet.

Yellow, gold-based dental casting alloy, Type 4

Ivoclar Vivadent Worldwide

- Australia, Austria, Brazil, Canada, China, Colombia, France, Germany, India, Italy, Japan, Mexico, New Zealand, Poland, Russia, Spain, Sweden, Switzerland, Taiwan, Turkey, USA, UK, Vietnam.

ISO 13485 Quality Management System Certified Made in USA. www.ivoclarvivadent.com

EN INSTRUCTIONS FOR USE

MODELLATION Wax to full contour for crown and bridge frame design. For composite, build up the framework in a reduced anatomic shape taking the planned veneer into consideration.

SPRUNG Provide the modeled bridge framework or coping with sprues of a suitable size. Use the direct or indirect technique being sure that the reservoir is positioned in the heat center.

INVESTMENT Weigh the wax pattern including the sprue to determine the quantity of the alloy to be used. (See wax conversion sheet/formula: weight x density = gr. of alloy).

BURN-OUT The suggested burnout temperature: 650-760C/1200-1400F High heat temperature investment: 480-540C/900-1000F

MELTING AND CASTING Use a separate carbon/ceramic crucible for each alloy. Use and new alloy must be in a ratio of 1:1. Depending on the type of casting machine, follow the manufacturer's instructions for use.

METAL PREPARATION Carefully divest and clean the object with AL2O3 glass beads, or a pickling agent (such as Prevox®). Do not use a hammer for divesting the object to prevent deformation.

HEAT TREATMENT Annealing: 760C/1400F for 15 minutes; quench immediately (water) Hardening: 455C/850F for 15 minutes; air cool.

SOLDERS AND FLUXES Design the soldering patty as small as possible and preheat it in the furnace at approximately 600C/112F. The soldering gap should be the same thickness as the soldering strip.

POLISHING After soldering or heat treatment, remove oxide and flux residue and finish and polish the framework with rubber finishers and polishers.

INDICATIONS Recommended for onlays, 3/4 crowns, crowns, telescope crowns, conus crowns, short and long span bridges, posts, implant superstructure.

CONTRAINDICATIONS For patients with known allergy/sensitivity to any major or minor constituents of this alloy, consultation with a physician is recommended.

SIDE EFFECTS In individual cases, sensitivity or allergies to components of this alloy may occur.

INTERACTIONS Galvanic effects may occur between different alloys in the same oral environment.

For additional information look into the alloy property chart.


French Federation Ivoclar Vivadent Marketing Ltd. Prospekt Andropova 18, Korpus 6 Office 10-06 115432 Moscow Tel: +7 499 418 0300 Fax: +7 499 418 0310 www.ivoclarvivadent.ru

SPRUNG METHOD

RECOMMENDATIONS DIRECT: single crowns, inlays and onlays INDIRECT: multiple units and multiple single crowns

INSTRUCTIONS: 1. Select a sprue with a reservoir equal to or larger than the thickest cross-section of the restoration.

2. Maintain the reservoir(s) within the heat center of the investment; keep the restoration(s) approximately 5 mm from the end of the investment and no closer than 5 mm from the sides. 3. Make sure that the sprue is connected to the thickest part of the restoration.

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