**ivoclar vivadent**

Recommended for onlays, 3/4 crowns, crowns, telescope crowns, conus crowns, short and long span bridges, posts, implant superstructure.

**Indications**

70.7	Au	3.6	Pt	13.7	Ag	10.0	Cu	1.0	In	0.8	Sn	1.0	Zn	1.0	Fe	Other	<1.0
------	----	-----	----	------	----	------	----	-----	----	-----	----	-----	----	-----	----	-------	------

**Composition**

Yellow, gold-based dental casting alloy, Type 4

**Academy Gold™ XH**

**Instructions for Use**

**Instructions d'usage**

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**Brüksanvisning**

**Mode d'emploi**

**Brüksanvisning**

**Brüksanvisning**

**Brüksanvisning**

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**ivoclar Vivadent WorldWide**

<p><b>Australia</b></p> <p>ivoclar Vivadent Pty. Ltd. 1 - 5 Overseas Drive P.O. Box 367 Noble Park, Vic. 3174 Tel. +61 3 979 595 99 Fax +61 3 979 596 45 www.ivoclarvivadent.com.au</p> <p><b>Brazil</b></p> <p>ivoclar Vivadent Ltda. Alameda Caiapós, 723 Centro Empresarial Tamboaré CEP 06460-110 Barueri - SP Tel. +55 11 2424 7400 Fax +55 11 2424 7440 www.ivoclarvivadent.com.br</p> <p><b>Canada</b></p> <p>ivoclar Vivadent Inc. 1-6600 Dixie Road Mississauga, Ontario L5T 2Y2 Tel. +1 905 670 8459 Fax +1 905 670 3102 www.ivoclarvivadent.ca</p> <p><b>China</b></p> <p>ivoclar Vivadent (Shanghai) Trading Co., Ltd. 21F Building 1, 881 Wuding Road Jing An District Shanghai 200040 Tel. +86 21 6032 1657 Fax +86 21 6176 0968 www.ivoclarvivadent.cn</p> <p><b>Colombia</b></p> <p>ivoclar Vivadent Marketing Ltd. Calla 13a No. 7-B-83, Of. 520 Bogotá Tel. +57 1 673 33 99 Fax +57 1 633 16 63 www.ivoclarvivadent.co</p> <p><b>France</b></p> <p>ivoclar Vivadent SAS Barren - SP Tel. +33 1 2424 7400 Fax +33 1 2424 7440 www.ivoclarvivadent.com.fr</p> <p><b>Germany</b></p> <p>ivoclar Vivadent GmbH Dr. Adolf-Schneider-Str. 2 73479 Ellwangen, Jagst Tel. +49 79 61 8 89 9 Fax +49 79 61 83 26 www.ivoclarvivadent.de</p> <p><b>India</b></p> <p>ivoclar Vivadent Marketing (India) Pvt. Ltd. 503/504 RahaJa Plaza 15 B Shah Industrial Estate Veera Desai Road, Andheri (West) Mumbai, 400 053 Tel. +91 22 2673 0302 Fax +91 22 2673 0301 www.ivoclarvivadent.in</p> <p><b>Italy</b></p> <p>ivoclar Vivadent s.r.l. Via Isozno 67/69 40033 Casalecchio di Reno (BO) Tel. +39 051 6113555 Fax +39 051 6113565 www.ivoclarvivadent.it</p> <p><b>Japan</b></p> <p>ivoclar Vivadent K.K. 1-28-24-4F Hongo Bunkyo-ku Tokyo 113-0033 Tel. +81 3 6903 3535 Fax +81 3 5844 3657 www.ivoclarvivadent.jp</p> <p><b>Mexico</b></p> <p>ivoclar Vivadent S.A. de C.V. Av. Insurgentes Sur No. 863 Piso 14, Col. Napoléon 03810 México, D.F. Tel. +52 55 5062 1000 Fax +52 55 5062 1029 www.ivoclarvivadent.com.mx</p> <p><b>New Zealand</b></p> <p>ivoclar Vivadent Ltd. 12 Omega Street, Rosedale PO Box 303011 North Harbour Auckland 0751 Tel. +64 9 914 9999 Fax +64 9 914 9990 www.ivoclarvivadent.co.nz</p> <p><b>Palau</b></p> <p>ivoclar Vivadent Sp. z o.o. ul. Jana Pawla II 78 00175 Warszawa Tel. +48 22 635 54 96 Fax +48 22 635 54 69 www.ivoclarvivadent.pl</p> <p><b>Russian Federation</b></p> <p>ivoclar Vivadent Marketing Ltd. Prospekt Andropova 18, Korpus 6 Office 10-06 115432 Moscow Tel. +7 499 418 0300 Fax +7 499 418 0310 www.ivoclarvivadent.ru</p> <p><b>Spain</b></p> <p>ivoclar Vivadent S.L.U. Cta. Fuencarral, 24-Portal I Bajo 28108 Alcobendas (Madrid) Tel. +34 913 75 78 20 Fax +34 913 75 78 38 www.ivoclarvivadent.es</p> <p><b>Sweden</b></p> <p>ivoclar Vivadent AB Dalsvägen 14 169 56 Solna Tel. +46 8 514 93 930 Fax +46 8 514 93 940 www.ivoclarvivadent.se</p> <p><b>Turkey</b></p> <p>ivoclar Vivadent Liaison Office Tesvikkiye Mahallesi Sakayık Sokak Nisantasi Plaza No: 40 Kat: 6 Daire: 31-32 34021 Sisi Istanbul Tel. +90 212 343 0802 Fax +90 212 343 0842 www.ivoclarvivadent.com.tr</p> <p><b>UK</b></p> <p>ivoclar Vivadent Limited Compass Building Feldspar Close Warrens Business Park Enderby Leicester LE19 4SE Tel. +44 116 284 78 80 Fax +44 116 284 78 81 www.ivoclarvivadent.com.uk</p>	<p><b>EC</b> <b>REP</b></p> <p><b>Manufacturer:</b> ivoclar Vivadent Inc. 175 Pineview Drive Amherst, NY 14228 USA Tel. +1 800 533 6825 Fax +1 716 691 2285 www.ivoclarvivadent.us</p>
--	--

**EN INSTRUCTIONS FOR USE**

**MODELLATION**  
Wax to full contour for crown and bridge frame design. For composite, build up the framework in a reduced anatomic shape taking the planned veneer into consideration. Single crowns require a thickness of minimum 0.3 mm; abutment crowns thickness a minimum of 0.5 mm. Make sure the framework demonstrates adequate stability of shape. Avoid sharp angles. Design the connector areas to be adequate for the position and alloy being used. If a composite or resin veneer is required, mechanical retention is recommended.

**SPRUNG**  
Provide the modeled bridge framework or coping with sprues of a suitable size. Use the direct or indirect technique being sure that the reservoir is positioned in the heat center. The connection sprues between the reservoir and the coping should be 2.5-3.0 mm in length and width.

**INVESTMENT**  
Weigh the wax pattern including the sprue to determine the quantity of the alloy to be used. (See wax conversion sheet/formula: weight x density = gr. of alloy). Use investment following the manufacturer's instruction.

**BURN-OUT**  
The suggested burnout temperature:  
High heat temperature investment: 650-760C/1200-1400F  
Low heat temperature investment: 480-540C/900-1000F

**MELTING AND CASTING**  
Use a separate carbon/ceramic crucible for each alloy. Used and new alloy must be in a ratio of 1:1. Depending on the type of casting machine, follow the manufacturers instructions for use. Ideally a compressed air and natural gas torch should be used to melt C&B alloys because propane and oxygen is much too hot and can easily overheat these alloys. If you are using propane and oxygen the pressure should be a 0.15 bar/2 psi for propane and 0.35 bar/5 psi for oxygen. Keep the alloy in the reducing atmosphere (if the flame is between the inner and outer cones. Use casting flux if needed. After casting bench cool to room temperature.

**Casting Temperature:** 980-1040C/1800-1905F

**METAL PREPARATION**  
Carefully divest and clean the object with Al<sub>2</sub>O<sub>3</sub> glass beads, or a pickling agent (such as Prevox®). Do not use a hammer for divesting the object to prevent deformation. Finish the framework with carbide burs or with ceramic-bonded grinding instruments. Avoid inhalation of dust during grinding!

**HEAT TREATMENT**  
**Annealing:** 760C/1400F for 15 minutes; quench immediately (water)  
**Hardening:** 345C/650F for 30 minutes; air cool.

**SOLDERS AND FLUXES**  
Design the soldering patty as small as possible and preheat it in the furnace at approximately 600C/112F. The soldering gap should be the same thickness as the soldering strip. Allow the object to cool slowly after soldering.  
**Solder:** .615, .585 Fine Solder  
**Flux:** Bondal Flux  
**Laser weld material:** Laser C&B Yellow

**POLISHING**  
After soldering or heat treatment, remove oxide and flux residue and finish and polish the framework with rubber finishers and polishers.

**INDICATIONS**  
Recommended for onlays, 3/4 crowns, crowns, telescope crowns, conus crowns, short and long span bridges, posts, implant superstructure.

**CONTRAINDICATIONS**  
For patients with known allergy/sensitivity to any major or minor constituents of this alloy, consultation with a physician is recommended.

**SIDE EFFECTS**  
In individual cases, sensitivity or allergies to components of this alloy may occur.

**INTERACTIONS**  
Galvanic effects may occur between different alloys in the same oral environment.  
**For additional information look into the alloy property chart.**

**SPRUNG METHOD**

**RECOMMENDATIONS**

**DIRECT:** single crowns, inlays and onlays

**INDIRECT:** multiple units and multiple single crowns

**INSTRUCTIONS:**

- Select a sprue with a reservoir equal to or larger than the thickest cross-section of the restoration.
- Maintain the reservoir(s) within the heat center of the investment; keep the restoration(s) approximately 5 mm from the end of the investment and no closer than 5 mm from the sides.
- Make sure that the sprue is connected to the thickest part of the restoration.
- The connection between the sprue and the restoration should be flared ("trumpet" configuration); eliminating turbulence of the alloy (causing investment erosion) and facilitating an undisturbed flow of alloy during casting and solidification.
- Maintain proper alloy weight in order to eliminate the negative effect of a too large burden during solidification. The rule of thumb for proper alloy weight determination is: total wax weight x specific gravity of the alloy = proper alloy weight.

**SUGGESTIONS:**

- Thickness and design of wax pattern: follow your manufacturers instructions.
- Use chill-vents when casting heavy and/or large restorations.

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