EN INSTRUCTIONS FOR USE
MODELLATION
Wax to full contour for crown and bridge frame design. For composite, build up the framework in a reduced anatomic shape taking the planned veneer into consideration. Single crowns require a thickness of minimum 0.3 mm; abutment crowns thickness a minimum of 0.5 mm. Make sure the framework demonstrates adequate stability of shape. Avoid sharp angles. Design the connector areas to be adequate for the position and alloy being used. If a composite or resin veneer is required, mechanical retention is recommended.
SPRUNG
Provide the modeled bridge framework or coping with sprues of a suitable size. Use the direct or indirect technique being sure that the reservoir is positioned in the heat center. The connection sprues between the reservoir and the coping should be 2.5-3.0 mm in length and width.
INVESTMENT
Weigh the wax pattern including the sprue to determine the quantity of the alloy to be used. (See wax conversion sheet/formula: weight x density = gr. of alloy). Use investment following the manufacturer's instruction.
BURN-OUT
The suggested burnout temperature:
High heat temperature investment: 650-760C/1200-1400F
Low heat temperature investment: 480-540C/900-1000F
MELTING AND CASTING
Use a separate carbon/ceramic crucible for each alloy. Used and new alloy must be in a ratio of 1:1. Depending on the type of casting machine, follow the manufacturers instructions for use. Ideally a compressed air and natural gas torch should be used to melt C&B alloys because propane and oxygen is much too hot and can easily overheat these alloys. If you are using propane and oxygen the pressure should be a 0.15 bar/2 psi for propane and 0.35 bar/5 psi for oxygen. Keep the alloy in the reducing atmosphere until the flame between the inner and outer cones. Use casting flux if needed. After casting bench cool to room temperature.
Casting Temperature: 1000-1080C/1830-1980F
METAL PREPARATION
Carefully divest and clean the object with Al₂O₃ glass beads, or a pickling agent (such as Prevox®). Do not use a hammer for divesting the object to prevent deformation. Finish the framework with carbide burs or with ceramic-bonded grinding instruments. Avoid inhalation of dust during grinding!
HEAT TREATMENT
Annealing: 760C/1400F for 15 minutes; quench immediately (water)
Hardening: 230C/450F for 15 minutes; air cool.
SOLDERS AND FLUXES
Design the soldering patty as small as possible and preheat it in the furnace at approximately 600C/112F. The soldering gap should be the same thickness as the soldering strip. Allow the object to cool slowly after soldering.
Solder: .650, .615, .585 Fine Solder
Flux: Bondal Flux
Laser weld material: Laser C&B Yellow
POLISHING
After soldering or heat treatment, remove oxide and flux residue and finish and polish the framework with rubber finishers and polishers.
INDICATIONS
Recommended for inlays, onlays, 3/4 crowns, crowns, telescope crowns, conus crowns, short and long span bridges, posts, implant superstructure, partial dentures.
CONTRAINDICATIONS
For patients with known allergy/sensitivity to any major or minor constituents of this alloy, consultation with a physician is recommended.
SIDE EFFECTS
In individual cases, sensitivity or allergies to components of this alloy may occur.
INTERACTIONS
Galvanic effects may occur between different alloys in the same oral environment.
For additional information look into the alloy property chart.

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