ivoclar vivadent

Indications

Recommended for Inlays, Onlays, 3/4 Crowns, Crowns

Composition

Au 76.8, Pt 1.0, Pd 1.0, Ag 12.8, Cu 8.3, Sn 2.9, Zn 1.0, In 1.0, Ga 1.0, Ni 1.0, Fe 0.5, Other <1.0

Brüksanvisning

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Mode d'emploi

Gebruichsinformatiön

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Instructions for Use

Rich yellow, gold-based dental casting alloy, Type 2

Harmony Medium

Ivoclar Vivadent Worldwide

<p>Australia Ivoclar Vivadent Pty. Ltd. 1 - 5 Overseas Drive P.O. Box 367 Noble Park, Vic. 3174 Tel. +61 3 979 595 99 Fax +61 3 979 596 45 www.ivoclarvivadent.com.au</p> <p>Brazil Ivoclar Vivadent Ltda. Alameda Caiapós, 723 Centro Empresarial Tamboré CEP 06460-110 Barueri - SP Tel. +55 11 2424 7400 Fax +55 11 2424 7440 www.ivoclarvivadent.com.br</p> <p>Canada Ivoclar Vivadent Inc. 1-6600 Dixie Road Mississauga, Ontario L5T 2Y2 Tel. +1 905 670 8459 Fax +1 905 670 3102 www.ivoclarvivadent.com</p> <p>China Ivoclar Vivadent (Shanghai) Trading Co., Ltd. 21F Building 1, 881 Wuding Road Jing An District Shanghai 200040 Tel. +86 21 6032 1657 Fax +86 21 6176 0968 www.ivoclarvivadent.cn</p> <p>Colombia Ivoclar Vivadent Marketing Ltd. Calla 13A No. 7-B-83, Of. 520 Bogotá Tel. +57 1 671 63 99 Fax +57 1 633 16 63 www.ivoclarvivadent.com</p>	<p>France Ivoclar Vivadent SAS B.P. 118 74410 Saint-Jorioz Tel. +33 450 88 64 00 Fax +33 450 88 91 52 www.ivoclarvivadent.fr</p> <p>Germany Ivoclar Vivadent GmbH Dr. Adolf-Schneider-Str. 2 73479 Ellwangen, Jagst Tel. +49 79 61 8 89 9 Fax +49 79 61 83 26 www.ivoclarvivadent.de</p> <p>India Ivoclar Vivadent Marketing (India) Pvt. Ltd. 503/504 Rahaqa Plaza 15 B Shah Industrial Estate Veera Desai Road, Andheri (West) Mumbai, 400 053 Tel. +91 22 2673 0302 Fax +91 22 2673 0301 www.ivoclarvivadent.in</p> <p>Italy Ivoclar Vivadent s.r.l. Via Isonzo 67/69 40033 Casalecchio di Reno (BO) Tel. +39 051 6113555 Fax +39 051 6113565 www.ivoclarvivadent.it</p>	<p>Russian Federation Ivoclar Vivadent Marketing Ltd. Prospekt Andropova 18, Korpus 6 Office 10-06 115432 Moscow Tel. +7 499 418 0300 Fax +7 499 418 0310 www.ivoclarvivadent.ru</p> <p>Spain Ivoclar Vivadent S.L.U. C/ta. Fuencarral, 24-Portal 1 Bajó 28108 Alcobendas (Madrid) Tel. +34 913 75 78 20 Fax +34 913 75 78 38 www.ivoclarvivadent.es</p> <p>Sweden Ivoclar Vivadent AB Dalsvägen 14 169 56 Solna Tel. +46 8 514 93 930 Fax +46 8 514 93 940 www.ivoclarvivadent.se</p> <p>Turkey Ivoclar Vivadent Liason Office Tesvikkiye Mahallesi Sakayik Sokak Nisantasi Plaza No: 40 Kat: 6 Daire: 31-32 34021 Sisi Istanbul Tel. +90 212 343 0802 Fax +90 212 343 0842 www.ivoclarvivadent.com.tr</p> <p>UK Ivoclar Vivadent Limited Compass Building Feldspar Close Warrens Business Park Enderby Leicester LE19 4SE Tel. +44 116 284 78 89 Fax +44 116 284 78 81 www.ivoclarvivadent.co.uk</p>	<p>Japan Ivoclar Vivadent K.K. 1-28-24-4F Hongo Bunkyo-ku Tokyo 113-0033 Tel. +81 3 6903 3535 Fax +81 3 5844 3657 www.ivoclarvivadent.jp</p> <p>Mexico Ivoclar Vivadent S.A. de C.V. Av. Insurgentes Sur No. 863 Piso 14, Col. Napoléon 03810 México, D.F. Tel. +52 55 5062 1000 Fax +52 55 5062 1029 www.ivoclarvivadent.com.mx</p> <p>New Zealand Ivoclar Vivadent Ltd. 12 Omega Street, Rosedale PO Box 303011 North Harbour Auckland 0751 Tel. +64 9 914 9999 Fax +64 9 914 9990 www.ivoclarvivadent.co.nz</p> <p>Poland Ivoclar Vivadent Sp. z o.o. ul. Jana Pawla II 78 00175 Warszawa Tel. +48 22 635 54 96 Fax +48 22 635 54 69 www.ivoclarvivadent.pl</p>
--	---	---	--

EN INSTRUCTIONS FOR USE

MODELLATION
Wax to full contour for crown and bridge frame design. For composite, build up the framework in a reduced anatomic shape taking the planned veneer into consideration. Single crowns require a thickness of minimum 0.3 mm; abutment crowns thickness a minimum of 0.5 mm. Make sure the framework demonstrates adequate stability of shape. Avoid sharp angles. Design the connector areas to be adequate for the position and alloy being used. If a composite or resin veneer is required, mechanical retention is recommended.

SPRUNG
Provide the modeled bridge framework or coping with sprues of a suitable size. Use the direct or indirect technique being sure that the reservoir is positioned in the heat center. The connection sprues between the reservoir and the coping should be 2.5-3.0 mm in length and width.

INVESTMENT
Weigh the wax pattern including the sprue to determine the quantity of the alloy to be used. (See wax conversion sheet/formula: weight x density = gr. of alloy). Use investment following the manufacturer's instruction.

BURN-OUT
The suggested burnout temperature:
High heat temperature investment: 650-760C/1200-1400F
Low heat temperature investment: 480-540C/900-1000F

MELTING AND CASTING
Use a separate carbon/ceramic crucible for each alloy. Used and new alloy must be in a ratio of 1:1. Depending on the type of casting machine, follow the manufacturers instructions for use. Ideally a compressed air and natural gas torch should be used to melt C&B alloys because propane and oxygen is much too hot and can easily overheat these alloys. If you are using propane and oxygen the pressure should be a 0.15 bar/2 psi for propane and 0.35 bar/5 psi for oxygen. Keep the alloy in the reding atmosphere if the flame between the inner and outer cones. Use casting flux if needed. After casting bench cool to room temperature.

Casting Temperature: 1000-1060C/1830-1940F

METAL PREPARATION
Carefully divest and clean the object with Al₂O₃ glass beads, or a pickling agent (such as Prevox®). Do not use a hammer for divesting the object to prevent deformation. Finish the framework with carbide burs or with ceramic-bonded grinding instruments. Avoid inhalation of dust during grinding!

HEAT TREATMENT
Annealing: 705C/1300F for 30 minutes; quench immediately in water.
Hardening: 205C/400F for 30 minutes; air cool.

SOLDERS AND FLUXES
Design the soldering patty as small as possible and preheat it in the furnace at approximately 600C/1112F. The soldering gap should be the same thickness as the soldering strip. Allow the object to cool slowly after soldering.
Solder: .650, .615, .585 Fine Solder
Flux: Bondal Flux
Laser weld material: Laser C&B Yellow

POLISHING
After soldering or heat treatment, remove oxide and flux residue and finish and polish the framework with rubber finishers and polishers.

INDICATIONS
Recommended for inlays, onlays, 3/4 crowns, crowns.

CONTRAINDICATIONS
For patients with known allergy/sensitivity to any major or minor constituents of this alloy, consultation with a physician is recommended.

SIDE EFFECTS
In individual cases, sensitivity or allergies to components of this alloy may occur.

INTERACTIONS
Galvanic effects may occur between different alloys in the same oral environment.
For additional information look into the alloy property chart.

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