Recommended for onlays, 3/4 crowns, telescopes, crowns, conus crowns, posts, short and long span bridges, PFM crowns, implant superstructures, partial dentures.

Indications

Table with 2 rows and 12 columns showing composition percentages for different materials.

Composition

Instructions for Use, Gebrauchsinformation, Modus d'emploi, Bruksanvisning, Instrucciones de Uso, Käyttöohjeet, and Brugsanvisning.

White, cobalt based dental metal ceramic casting alloy, Type 5

Worldwide distribution list for Ivoclar Vivadent, including Australia, Brazil, Canada, China, Colombia, Germany, India, Japan, Mexico, New Zealand, Poland, Russia, Sweden, Switzerland, Taiwan, Turkey, UK, USA, and Venezuela.

ISO 13485 Quality Management System Certified, Made in U.S.A., and Ivoclar Vivadent logo.

EN INSTRUCTIONS FOR USE

MODELLATION Design the framework in a reduced anatomic shape taking the planned veneer into consideration. Single crowns require a thickness of minimum 0.3 mm; abutment crowns thickness a minimum of 0.5 mm.

SPRUNG Provide the modelled bridge framework or coping with sprues of a suitable size. Use the direct or indirect technique being sure that the reservoir is positioned in the heat center. The connection sprues between the reservoir and the coping should be 2.5-3.0 mm in length and width.

MELTING AND CASTING Use a separate ceramic crucible for each alloy and preheat the crucible in the burnout furnace. Use and new alloys must be in a ratio of 1:1. Depending on the type of casting machine, follow the manufacturer's instructions for use.

OXIDATION Not required. Optional Blast surface with 50 micron Al2O3 at max. 4.5 bar/65 psi pressure before oxidation.

SOLDERS AND FLUXES Design the soldering patty as small as possible and preheat it in the furnace at approximately 600C/1112F. The soldering gap should be the same thickness as the soldering strip.

CONTRAINDICATIONS For patients with known allergy/sensitivity to any major or minor constituents of this alloy, consultation with a physician is recommended.

SIDE EFFECTS In individual cases, sensitivity or allergies to components of this alloy may occur. Ivoclar Vivadent makes no claims regarding the MRI-compatibility of its dental alloys.

INTERACTIONS Galvanic effects may occur between different alloys in the same oral environment.

SPRUNG METHOD

RECOMMENDATIONS: DIRECT: single crowns, inlays and onlays. INDIRECT: multiple units and multiple single crowns.

INSTRUCTIONS: 1. Select a sprue with a reservoir equal to or larger than the thickest cross-section of the restoration.

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