IPS Ceramic Etching Gel

English

Instructions for Use

Description

The use of IPS Ceramic Etching Gel generates retentive bonding surfaces on glass-ceramic restorations as a preparation for cementation. The etch pattern results in enhanced bonding strength between the cementation material and the glass-ceramic.

IPS Ceramic Etching Gel is for professional use only and must not be applied in the oral cavity.

Composition
< 5 % hydrofluoric acid HF
Indication
Etching of the bonding surfaces to increase the micro-retentions on glass-ceramic restorations (leucite, lithium disilicate, Fluorapatite ceramics) as a preparation for incorporation.
Contraindication
– Any application in the oral cavity
– Etching of oxide ceramics (e.g. IPS e.max

Processing
Note
Preferably etch the reconstruction immediately before incorporation of the restoration. This prevents contamination during transport or try-in.

terial to be etched. Prolonging the reaction time does not additionally enhance the bonding strength between th estora ion and the cementation material.

– IPS Empress CAD	60 sec.
– IPS Classic®	60 sec.
– IPS InLine®	120 sec.
– IPS d.SIGN®	120 sec.
– IPS d.SIGN®	60 sec.

- Apply approximately 1 ml of IPS Ceramic Etching Gel from the restoration under running water into a cup (poly thyl-3. Apply approximately 1 ml of IPS Ceramic Etching Gel from the surfaces and apply the bonding agent (e.g. Monobond-S) according to the Instructions for Use in conjunction with the bonding agent.
- Before cementation, rinse the thoroughly dry the ceramic bonding

Neutralization of the IPS Ceramic Etching Gel

- To neutralize the diluted solution of etching gel and water, add the neutralizing powder and let it react for 5 minutes.
- 1 measuring spoon of IPS Ceramic Neutralizing Powder is suitable to neutralize the quantity of IPS Ceramic Etching Gel required for a large MOD inlay.

Note

The hydrofluoric acid in IPS Ceramic Etching Gel attacks quartz, silicate and borate glasses, as well as sanitary ceramics and various metals and alloys (e.g. high-grade steel), Nickel, copper, polyethylene, PVC, and Teflon are resistant to hydrofluoric acid.

Shelf life
See label / packaging
Special safety notes for the use of IPS Ceramic Etching Gel:
– Hydrofluoric acid is highly toxic. It is strongly corrosive and does not cause any warning pain on the surface of skin and mucous membranes, but causes subsequent, painful in-depth effect.
– Contact with skin, eyes, and clothing must be prevented at all costs, since the material is extremely toxic and corrosive. In case of accidental skin contact, immediately rinse with copious amount of running water and immediately contact a physician!
– Wear acid resistant protective gloves, protective clothing, and goggles when working with IPS Ceramic Etching Gel.
– Close package immediately after use.
– Neutralize spilled IPS Ceramic Etching Gel with IPS Ceramic Neutralizing Powder and rinse with copious amounts of water.
– Store out of the reach of children!

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