Viteo® Base Ti

Viteo® Base Ti

EN Instructions for Use

DE Gebrauchsinformation

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SV Bruksanvisning

DA Brugsanvisning

Viteo® Base Ti System:

Viteo® Base Ti
SD

Viteo® Base Ti
MD

Viteo® Screw

Viteo® Base
Trimmer SD+MD

Viteo® Base
Press Sleeve SD

Viteo® Base
Press Sleeve MD

Viteo® Channel
Pin

Viteo® Holder

Product Description

The Viteo Base Ti, with its blasted bonding surface and internal anti-rotation protection, is suitable for PRESS and CAD/CAM-produced restoration materials. The Viteo Base Ti, which can be shortened, is compatible with various implant systems.

The Viteo Base Ti is a prefabricated, prosthetic component for endosseous dental implants to support prosthetic restorations in partially or fully edentulous jaws. The Viteo Base Ti is secured to the dental implant with a screw.

Indications

- Implant-retained single-tooth restorations

Contraindications

- For restorations with angulation correction of more than 20° to the implant axis
- Modification of the implant interface
- Reduction of the marginal step width
- Modification of the abutment geometry, with the exception of shortening the abutment height
- Shortening the abutment height to less than 4 mm
- Shortening using unsuitable instruments
- Multiple use of the same Viteo Base Ti
- Bruxism
- Multi-unit restorations
- As a cast-on abutment coping
- Other applications not listed in the indications above
- Applications beyond the indication for use of the implant manufacturer
- The use of third party components

Caution:

Small diameter implants and angled abutments are not recommended for use in the posterior region of the mouth. Please note warnings, cautions and pre-cautions of the implant manufacturer.

Side Effects

Viteo Base Ti should not be used if the patient is known to be allergic to any of the components. The Viteo Base Ti has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artefact in the MR environment. The safety of Viteo Base Ti in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

Composition (constituents)

- Viteo Base Ti: Ti-6Al-4V ELI
- Viteo Screw: Ti-6Al-4V ELI

Directions for Use**Preparing the Viteo Base Ti**

Place the Viteo Base Ti selected for the implant system onto the relevant laboratory analog and secure with the appropriate screw to max. 5 Ncm. Align the anti-rotation protection (vertical groove) in a mesial or distal direction. Ensure there is sufficient space between the Viteo Base Ti and the antagonist teeth. The recommended thickness of the restoration material must be observed (see chapter Manufacturer's Information).

If it is necessary to shorten the abutment height, cut to the required height using a thin separating disc held at a 90° angle. The laser marking shows the minimum required abutment height of 4 mm! First deburr the cut outer edges before placing the Viteo Base Trimmer guidance sleeve on top. Then reduce using the Viteo Base Trimmer (SD or MD). Place the guidance sleeve onto the Viteo Base Ti and push down into the end position. Insert the Viteo Base Trimmer vertically into the guidance sleeve and then grind the surface with gentle pressure at max. 6000 rpm. Clean away the grindings from time to time with compressed air. Remove the Viteo Base Trimmer from the guidance sleeve.

If the Viteo Base Trimmer has contact with the inner side of the guidance sleeve, then the minimum abutment height of 4 mm has been reached. Smooth over the prepared surface of the screw channel with a rotating instrument and smooth the front surface using a silicone polisher. If unsuitable rotating instruments are used or if the surface is ground incorrectly, the area can become hot and cause damage to the Viteo Base Ti.

Note:

When using Viteo Base Trimmer, please wear protective glasses.

Viteo
Base Ti
6 mm

Viteo
Base Ti
4 mm

Application Viteo Base Trimmer

If the blasted bonding surface has been damaged or ground, then re-blast the area. Place the Viteo Base Ti into the ultrasonic bath and clean using a steam cleaner, dry with compressed air and then screw onto a laboratory analog. Apply hard modelling wax onto the emergence profile and into the screw channel in order to protect these areas when blasting. Blast the bonding surface carefully using Al₂O₃, 50 µm, and max. 2 bar pressure until the surface is evenly matt. Avoid blasting the area below the platform with Al₂O₃, as this may

LOT

Rx ONLY

For dental use only

CE 0123

Manufacturer:

Ivoclar Vivadent AG

FL-9494 Schaan/Liechtenstein

www.ivoclarvivadent.com

Date information prepared:

2019-03-08/Rev. 2

689827/WE1+SV+DA

ivoclar
digital

ivoclar
digital

cause instability and damage to the implant interface. Remove all residues of the modelling wax with an instrument and clean using a steam cleaner.

Note:

Check every package and individual product to ensure of the entirety and validity of all contents. The supplied screw is suitable for use in the laboratory and also for clinical application. If used for clinical application, only new/unused screws may be used. Viteo Base Ti with less than 4 mm abutment height must not be used. The Viteo Base Trimmer is designed for multiple use. The guidance sleeve must be exchanged regularly. By blasting the bonding surfaces with Al_2O_3 , the adhesion of the cemented restoration is improved. Prerequisites for this are that the bonding surfaces are cleaned thoroughly with a steam cleaner before cementation and that the surface conditioning is carried out according to the manufacturer's instructions.

Fabrication of the restoration for the PRESS technique

Select the appropriate Viteo Base Press Sleeve (SD/MD) and place onto the Viteo Base Ti. Check the fitting accuracy at the screw channel and marginal edge. Apply a separating agent to the Viteo Base Ti in the transition area to the Viteo Base Press Sleeve before fabricating the restoration in wax or with a burn-out acrylic resin. In order to lengthen the screw channel (if necessary), apply a separating agent to the Viteo Screw Channel Pin and insert it into the screw channel. Lengthen using either modelling wax or modelling acrylic. If the abutment height is shortened, the Viteo Base Press Sleeve must also be shortened accordingly and the screw channel and Viteo Base Ti front surface must be adapted so that the restoration sits accurately in these areas. The restoration is fabricated according to the clinical situation. The design and finishing requirements must be observed. The modelled restoration is converted into ceramic using the press process (see chapter Manufacturer's Information).

Note:

It is not advisable to apply modelling wax or modelling acrylic directly to the blasted surface of the Viteo Base Ti! Do not treat the Viteo Base Press Sleeves with hot water or hot steam, due to the risk of deformation!

Scanning and design of the restoration for the CAD/CAM procedure

Select the relevant Scan Abutment which is suited to the type of implant used (e.g. ELOS Accurate Intra Oral Scan Abutments) and screw onto the laboratory analog using the supplied Viteo Scan Abutment screw driver. Determine the required abutment height (4 mm or 6 mm) before creating the order in the CAD software. Then, when placing the order in the CAD software, select the relevant implant system, implant diameter and Viteo Base Ti with the determined abutment height and platform diameter from the Viteo Base Ti library. Carry out the scan process and create the desired restoration design. Next, transfer the finished design into the production process (observe the manufacturer's instructions). The current list of available scan abutments / software can be downloaded from the website (www.ivoclarvivadent.com).

Temporary fixation of the restoration on the Viteo Base Ti (preparation for clinical try-in)

A clinical try-in can be carried out before the ceramic restoration is permanently cemented onto the Viteo Base Ti. For this, the two components are temporarily fixed to one another using low-viscosity impression material (e.g. Virtual® Extra Light Body Fast Set). Screw the Viteo Base Ti onto a laboratory analog. Place the restoration onto the Viteo Base Ti and check the position and relationship by marking with a permanent waterproof pen. Clean the Viteo Base Ti and the restoration using a steam cleaner; then blow dry with compressed air. Seal the screw channel with the Viteo Screw Channel Pin. The length of the Viteo Screw Channel Pin can be shortened with a scalpel. Apply the silicone directly onto the bonding surface of the Viteo Base Ti and the restoration, and then bring both parts together. Align both parts according to the markings made earlier. Fix the objects to each other in the correct position until the silicone has completely set (check manufacturer's information). Remove excess impression material with a suitable instrument (e.g. scalpel). Carefully carry out the clinical try-in. The clinical try-in serves for the control of the emergence profile and contact points only.

Permanent cementation of the restoration on the Viteo Base Ti

Conditioning the Viteo Base Ti

Clean the Viteo Base Ti thoroughly in an ultrasonic bath and with a steam cleaner, and then allow to dry. The surface must have an even matt colouring. Before the definitive cementation, screw the Viteo Base Ti onto a laboratory analog. Place the restoration onto the Viteo Base Ti and mark the position and relationship with a waterproof permanent pen. After cleaning with a steam cleaner, take care not to contaminate the bonding surface as this will have a detrimental effect on the cementation. For the metal bond apply Monobond® Plus and allow it to take effect for 60 sec. After the exposure time has elapsed, blow away any excess and dry with water/oil-free compressed air (see chapter Manufacturer's Information).

Conditioning the bonding surface of the restoration

The pre-conditioning process required on the bonding surface of the restoration depends on the restoration material (see chapter Manufacturer's Information).

Restoration material	Pre-conditioning the bonding surface	Conditioning
Lithium disilicate glass-ceramic (IPS e.max® Press/IPS e.max® CAD)	20 s with IPS® Ceramic etching gel	60 s with Monobond® Plus
	20 s rubbing / 40 s exposure with Monobond Etch & Prime®	
PMMA (Telio® CAD)	30 s SR Connect, 40 s light polymerization	

Cementation with Multilink® Hybrid Abutment

Lay out the cleaned and conditioned components to be adhesively cemented (restoration / Viteo Base Ti fixed to laboratory analog). Before the application of Multilink® Hybrid Abutment, screw the Viteo Screw Channel Pin into the screw channel. The length of the Viteo Screw Channel Pin can be shortened using a scalpel. Alternatively, to protect the screw, insert a sponge pellet into the screw channel. Apply Multilink Hybrid Abutment directly from the syringe in a thin even layer onto the bonding surface of the Viteo Base Ti and to the bonding surface of the restoration. Place the restoration onto the Viteo Base Ti until the previously marked guidelines are aligned. Bring the parts together by exerting gentle pressure (restoration / Viteo Base Ti) and check the correct end position. In the end position, hold the parts together firmly for 5 sec. Remove the excess cement during the initial curing phase (2 – 3 min) with a suitable dental technical instrument (e.g. LeCron). Remove the Viteo Screw Channel Pin after the initial curing phase has elapsed (2 – 3 min) by turning and pulling. Fix the parts by applying gentle pressure (e.g. diamond tweezers). Apply glycerine gel to the cement joint to prevent the formation of an inhibition layer. The entire auto polymerization is complete after 7 min. Allow the cement to cure completely.

Rinse the glycerine gel off with water once the polymerization is finished. Then place the laboratory analog with the Viteo Base Ti and the restoration into the Viteo Holder. To open the grip, turn the clasp in an anticlockwise direction and press upwards. Place the laboratory analog into the open grip and tighten. Polish the adhesive joint with a rubber polisher at a low rotational speed (< 5000 rpm). Any cement deposits left in the screw channel can be removed with a suitable rotating instrument (see chapter Manufacturer's Information).

Note:

Do not move the objects; keep them in a fixed position until Multilink Hybrid Abutment has fully cured (e.g. diamond tweezers). When removing the Viteo Screw Channel Pin ensure that the restoration and the Viteo Base Ti remain in a fixed position and do not separate.

Clinical insertion

The Viteo Base Ti and the Viteo Screw are not delivered in a sterile condition. The Viteo Base Ti must be screwed in place using a brand new screw. Keep the original instruments required for screwing in position at the ready. In preparation for definitive insertion, remove the temporary restoration / healing caps, clean the implant lumen (e.g. Cervitec® Liquid) and check the peri-implant gingival conditions (emergence profile). Avoid gingival contusion between implant and restoration when inserting the restoration. Screw in the Viteo Screw to the level of torque recommended by the implant manufacturer (see combination table).

For further information on insertion, see chapter Manufacturer's Information.

Note:

The Viteo Base Ti is for single use only.

Sterilization

Steam sterilization takes place using the fractionated vacuum process. The steam sterilization process may only be carried out by devices compliant with national or local standards. The equipment and devices must be properly maintained and serviced at regular intervals.

Components	Sterilization parameters*
Viteo Base Ti / Viteo Screw	3 min / 132 °C/270 °F

*The local requirements and hygiene standards for dental practices must be observed.

Do not store the Viteo Base Ti after sterilization; instead it must be inserted immediately. Sterilization of the entire restoration is only possible in certain cases (Viteo Base Ti / restoration material) (see chapter Manufacturer's Information). The creator of the restoration is responsible for informing the dentist about the sterilization possibilities.

Manufacturer's information

In the following instruction manuals (IfU), there is additional information on the fabrication and insertion of implant-retained prosthetic restorations.

- IfU IPS e.max Press Abutment Solutions
- IfU IPS e.max CAD Abutment Solutions
- IfU Telio CAD Abutment Solutions
- IfU Multilink Hybrid Abutment

For further information on using Ivoclar Vivadent products, please contact our specialist consultants or visit our website www.ivoclarvivadent.com.

Keep out of reach of children!

Only for dental medical and dental technical use!

Some of the products and/or indications have not yet been released in all countries. Please contact your national Ivoclar Vivadent organization regarding

this situation.

This material was developed for use in the dental field and must be used according to the instructions provided. Liability cannot be accepted for damages resulting from failure to observe the instructions of the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the instructions. These regulations also apply if the materials are mixed or used in conjunction with products of other manufacturers.

Viteo®, IPS e.max®, Multilink®, Monobond®, Telio®, Virtual® are registered trademarks of Ivoclar Vivadent AG, Schaan

Nobel Biocare®, NobelActive™, Biomet 3i™, OSSEOTITE® Certain®, OsseoSpeed™, ELOS ACCURATE® are not registered trademarks of Ivoclar Vivadent AG, Schaan.

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- IFU IPS e.max Press Abutment Solutions
- IFU IPS e.max CAD Abutment Solutions
- IFU Telio CAD Abutment Solutions
- IFU Multilink Hybrid Abutment

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
Ivoclar Vivadent Pty. Ltd.

1 – 5 Overseas Drive, P.O. Box 367, Noble Park, Vic. 3174, Australia
Tel. +61 3 9795 9599, Fax +61 3 9795 9645
www.ivoclarvivadent.com.au

Ivoclar Vivadent GmbH

Tech Gate Vienna, Donau-City-Strasse 1, 1220 Wien, Austria
Tel. +43 1 263 191 10, Fax +43 1 263 191 111
www.ivoclarvivadent.at

Ivoclar Vivadent Ltda.

Alameda Caiapós, 723, Centro Empresarial Tamboré
CEP 06460-110 Barueri – SP, Brazil
Tel. +55 11 2424 7400
www.ivoclarvivadent.com.br

Ivoclar Vivadent Inc.

1-6600 Dixie Road, Mississauga, Ontario, L5T 2Y2, Canada
Tel. +1 905 670 8499, Fax +1 905 670 3102
www.ivoclarvivadent.us

Ivoclar Vivadent Shanghai Trading Co., Ltd.

2/F Building 1, 881 Wuding Road, Jing An District, 200040 Shanghai, China
Tel. +86 21 6032 1657, Fax +86 21 6176 0968
www.ivoclarvivadent.com

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Ivoclar Vivadent Limited

Compass Building, Feldspar Close, Warrens Business Park
Enderby, Leicester LE19 4SD, United Kingdom
Tel. +44 116 284 7880, Fax +44 116 284 7881
www.ivoclarvivadent.co.uk

Ivoclar Vivadent, Inc.

175 Pineview Drive, Amherst, N.Y. 14228, USA
Tel. +1 800 533 6825, Fax +1 716 691 2285
www.ivoclarvivadent.us

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Ivoclar Vivadent Marketing Ltd.
Ivoclar Vivadent AB