## Instructions for Use

### ■ PRODUCT DESCRIPTION

Au-based dental metal-ceramic alloy, Type 4

### ■ INDICATIONS\*

Crowns, Telescope Crowns, Conus Crowds, Wide Bridges, Cast Posts/-Cores, Bars, Attachments, Implant Superstructures, Partial Dentures

### ■ WAXING/MODELLATION

Design the framework in a reduced anatomical shape taking the planned veneer into consideration. Single crowns require a minimum thickness of 0.3 mm. Abutment crowns require a minimum thickness of 0.5 mm. The waxed framework design provides adequate support for the veneering material. Avoid sharp angles. Connectors must have the required dimensions to provide resistance to deformation. Create large surface areas for planned soldering, with a gap of 0.05–0.2 mm.

### ■ SPRING

Provide the modeled single-tooth restoration or bridge framework with springs of a suitable size. In general the reservoir, sprue leads, and connector spurs, whether pear shaped or traditional, must be sized according to the specific technique used. When using the direct or indirect technique be sure that the reservoir is positioned directly behind the center of the connector between the reservoir and the casting should be a maximum of 2.5–3.0 mm in length and width. The wax pattern including the spurs must be weighed in grams in order to determine the needed amount of alloy. Wax conversion formula: wax weight (gram) x alloy density = grams of alloy required.

### ■ INVESTING

Use a phosphate-bonded investment material. Follow the manufacturer's instructions.

### ■ PREHEATING / BURN-OUT

Recommended burn-out temperature: 900 °C

### ■ MELTING AND CASTING

Flame: Propan 0.35 bar; Oxygen 0.7 bar  
Other specifics may be required by the type of casting machine. It is recommended to use a separate and clean ceramic crucible for each alloy. Preheat the ceramic crucible in the burnout furnace. The recommended ratio of used material to new material is 1:1. Do not use flux.  
Casting Temperature: 1400 °C

### ■ FRAMEWORK FINISHING

After bench cooling, carefully clean and clean the casting with aluminum oxide (Al<sub>2</sub>O<sub>3</sub>). Do not use a hammer for divesting. Finish the casting with carbide burs and/or with ceramic-bonded grinding instruments. Blast the surface with 50–110 micron aluminum oxide (Al<sub>2</sub>O<sub>3</sub>) at 2.0 bar. Subsequently, steam clean or ultrasonic clean with distilled water or ethanol and dry the framework.

### ■ OXIDATION

Place the framework on the firing tray providing adequate support. To achieve a uniform result follow the oxidation cycle.  
Temperature: 950 °C; Holding time: 5 min; Vacuum: No  
If the oxide layer is stained, grind and blast the surface again. Repeat the oxide firing. Use the appropriate ceramic veneering material, following the manufacturer's instructions.  
Use a ceramic investment with the recommended firing temperature: 1000 °C

### ■ HEAT TREATMENT

Hardening: 595 °C for 20 min; bench cool

### ■ SOLDERING AND LASER WELDING

The soldering gap should not be wider than the thickness of the soldering material. Use the soldered casting to test the joint under fire as sparingly.  
Pre Solder: Porta Lot 1090 W Flux: High Fusing Bondal Flux  
Post Solder: Porta OP Lot W-2 Flux: Bondal Flux  
Laser Welding Wire: Porta P6  
Laser-Schweiß-Draht: Porta P6

### ■ POLISHING

Carefully remove any oxide and flux residue. Smooth the metal surfaces with rubber polishers. Polish to a high gloss finish using polishing paste. Subsequently, clean with ultrasonic cleaning equipment or careful steam cleaning.

## ADDITIONAL SAFETY CONCERNS AND INSTRUCTIONS

### ■ CONTRAINDICATIONS

Pateients with known allergy/sensitivity to any major or minor elements of this alloy, consultation with a physician is recommended. Alloy is not to be used for any application not included within the indications.

### ■ SIDE EFFECTS

In individual cases, sensitivity or allergies to elements of this alloy may occur. Ivoclar Vivadent makes no claims regarding the MRI-compatibility of its dental alloys. It is recommended that the patient be made aware of the potential for dental alloys to affect MRI results and to disclose the presence of dental alloys to the MRI Technician prior to conducting a test.

### ■ INTERACTIONS

Galvanic effects may occur between different or dissimilar alloys in the same oral environment.

### ■ CAUTION

Metal vapors and metal dust are harmful if inhaled. Therefore, the use of extraction equipment and / or suitable protective masks is advised!

### ■ STORAGE CONDITION

Store in a dry environment at room temperature.

### ■ DISCLAIMER

This material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the products under suitable and safe conditions for any purpose not explicitly stated in the Instructions. These regulations also apply if the materials are used in conjunction with products of other manufacturers.

PROCESSING DATA		
Investment material:	phosphate-bonded	
Preheating/Burn-out Temperature:	900 °C	
Crucible:	Ceramic crucible	
Casting Temperature:	1400 °C	
Oxidation:	Temperature: 950 °C; Holding time: 5 min; Vacuum: No	
CTE:	(25–500 °C): 14.1 x 10 <sup>-6</sup> /K	
Refractive/Ceramic Material:	IP5 Style <sup>®</sup> , IP5 InLine <sup>®</sup>	
Hardening:	595 °C, for 20 min; bench cool	
Pre Solder / Flux:	Porta Lot 1090 W High Fusing Bondal Flux	
Post Solder / Flux:	Porta OP Lot W-2 Bondal Flux	
Laser Welding Wire:	Porta P6	

TECHNICAL DATA (ISO 22674:2016 & ISO 9693-1:2012)		
Type / Color:	4 White	4
Density (g/cm <sup>3</sup> ):	12.8	
Smelting Interval (Solidus/Liquidus):	1120–1260 °C	
Elastic Modulus (GPa):	125	
	Porcelain Fired	Hardened
Tickers Hardness:	290	260
Vickers Hardness (MPa):	800	840
Tensile Strength (MPa):	800	840
0.2% Proof Stress (MPa):	540	550
Elongation (%):	15	9

\* See TYPE CLASSIFICATION DUE TO PHYSICAL PROPERTIES

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