Instructions for Use

EN

■ PRODUCT DESCRIPTION

Pd-based dental metal-ceramic alloy, Type 4

■ INDICATIONS*

3/4 Crowns, Crowns, Telescopical Crowns, Conus Crowns, Bridges, Wide Bridges, Cast Posts /Cores, Bars, Attachments, Implant-Suprastructures, Partial Dentures

■ WAXING/MODELLATION

Design the framework in a reduced anatomical shape taking the planned veneer into consideration. Single crowns require a minimum thickness of 0.3 mm. Abutment crowns require a minimum thickness of 0.5 mm. The design should be such that the casting provides adequate support for the veneering material. Avoid sharp angles. Connectors must have the required dimensions to provide resistance to deformation. Create large surface areas for planned soldering, with a gap of 0.05–0.2 mm.

■ SPRING

Provide the modelled single-tooth restoration or bridge framework with springs of a suitable size. In general the reservoir, spring leads, and connector spurs, whether pear shaped or traditional, must be sized according to the specific technique. When using the direct or indirect technique be sure that the reservoir is positioned directly behind the connector between the reservoir and the casting should be a maximum of 1.1 mm. Springs 2.5–3.0 mm in length and width. The wax pattern including the spurs must be weighed in grams in order to determine the needed amount of alloy. Wax conversion formula: wax weight (gram) x alloy density = grams of alloy required

■ INVESTING

Use a phosphate-bonded investment material. Follow the manufacturer's instructions.

■ PREHEATING / BURN-OUT

Recommended burn-out temperature: 900 °C

■ MELTING AND CASTING

Flame: Propano 0.35 bar; Oxygen 0.7 bar

Other specifics may be required by the type of casting machine. It is recommended to use a separate and clean ceramic crucible for each alloy. Preheat the ceramic crucible in the burnout furnace. The recommended ratio of used material to new material is 1:1. Do not use flux.

■ FRAMEWORK FINISHING

After being cooled, carefully dress and clean the casting with aluminum oxide (Al₂O₃). Do not use a hammer for distressing. Finish the casting with carbide burs and/or with ceramic-bonded grinding instruments. Blast the surface with 50–110 micron aluminum oxide (Al₂O₃) at 2.0 bar. Subsequently, steam clean or ultrasonic clean with distilled water or ethanol and dry the framework.

■ OXIDATION

Place the framework on the firing tray providing adequate support. To achieve a uniform result follow the oxidation cycle.

Temperature: 930 °C; Holding time: 5 min; Vacuum: No

If the oxide layer is stained, grind and blast the surface again. Repeat the oxide firing. Use the appropriate ceramic veneering material, following the manufacturer's instructions.

Use a ceramic investment for firing temperature: 980 °C

■ HEAT TREATMENT

Hardening: 600 °C for 15 min; bench cool

■ SOLDERING AND LASER WELDING

The soldering gap should not be wider than the thickness of the soldering material. Allow the soldered casting to cool slowly to avoid cracking.

Pre Solder: Porta IP Lot V-2

Post Solder: Porta OP Lot W-2

Laser Soldering Wire: Porta P6

Flux: High Fusing Bondal Flux

Flux: Bondal Flux

■ POLISHING

Carefully remove any oxide and flux residue. Smooth the metal surfaces with rubber polishers. Polish to a high gloss finish using polishing paste. Subsequently, clean using ultrasonic cleaning equipment or careful steam cleaning.

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ADDITIONAL SAFETY CONCERNS AND INSTRUCTIONS

■ CONTRAINDICATIONS

For patients with known allergy/sensitivity to any major or minor elements of this alloy, consultation with a physician is recommended. Alloy is not to be used for any application not included within the indications.

■ SIDE EFFECTS

In individual cases, sensitivity or allergies to elements of this alloy may occur. Ivoclar Vivadent makes no claims regarding the MRI-compatibility of its dental alloys. It is recommended that the patient be made aware of the possibility of allergic reactions to dental alloys to alert MRI results and to disclose the presence of dental alloys to the MRI Technician prior to conducting a test.

■ INTERACTIONS

Galvanic effects may occur between different or dissimilar alloys in the same oral environment.

■ CAUTION

Metal vapors and metal dust are harmful if inhaled. Therefore, the use of extraction equipment and/or suitable protective masks is advised!

■ LAGERBEDINGUNGEN

Store in a dry environment at room temperature.

■ DISCLAIMER

This material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the products and their suitability for any purpose not explicitly stated in the Instructions. These regulations also apply if the materials are used in conjunction with products of other manufacturers.

PROCESSING DATA	phosphate-bonded
Investment Material:	phosphate-bonded
Preheating/Burn-out Temperature:	900 °C
Crucible:	Ceramic crucible
Casting Temperature:	1350 °C
Oxidation:	Temperature: 930 °C; Holding time: 5 min; Vacuum: No
CTE:	(25–500 °C): 14.8 x 10 ⁻⁶ K ⁻¹
Anteil keramische Material:	IP5 Style [®] ; IP5 InLine [®]
Hardening:	600 °C for 15 min; bench cool
Pre Solder / Flux:	Porta IP Lot V-2 High Fusing Bondal Flux
Post-soldadura / Flux:	Porta OP Lot W-2 Bondal Flux
Laser Welding Wire:	Porta P6
Technische Daten (ISO 22674:2016 & ISO 9693-1:2012)	
Type / Color:	4 White
Density (g/cm ³):	10.6
Schmelzintervall (Solidus/Liquidus):	1080–1200 °C
Elastic Modulus (GPa):	126
Material:	Porcelain Fired
Vickers Hardness:	270 280
Tensile Strength (MPa):	870 880
0.2% Proof Stress (MPa):	560 600
Elongation (%):	9 6

* See TYPE CLASSIFICATION DUE TO PHYSICAL PROPERTIES

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