CRT® buffer

Instructions for Use	Käyttöohjeet
Gebrauchsinformation	Bruksanvisning
Mode d'emploi	Productinformatie
Istruzioni d'uso	Οδηγίες Χρήσεως
Instrucciones de uso	Kullanma Talimatı
Instruções de Uso	Инструкция
Bruksanvisning	Instrukcja stosowania
Brugsanvisning	

- | | |
|---|--|
| - For determining the buffer capacity of saliva | - Sylijen puskurikapasiteetin määrittelyyn |
| - Zur Bestimmung der Pufferkapazität des Speichels | - For bestemmelse av bufferkapasiteten i spytt |
| - Pour la détermination du pouvoir tampon de la salive | - Voor het bepalen van de buffer capaciteit in het speeksel |
| - Per la determinazione del potere tampone della saliva | - Για τον προσδιορισμό της ρυθμιστικής ικανότητας του σάλιου |
| - Para la determinación de la capacidad buffer de la saliva | - Tükürüğün tampon kapasitesini belirlemek içinidir |
| - Para determinar a capacidade buffer da saliva | - для определения буферной способности слюны |
| - För bedömning av salivens buffertkapacitet | - Do okrešlania pojemnošci buforowej śliny |
| - Til bestemmelse af bufferkapaciteten i saliva | |

For dental use only!

Caution: U.S. Federal Law restricts this device to sale by or on the order of a licensed dentist.

Made in Liechtenstein

Ivoclar Vivadent AG, FL-9494 Schaan/Liechtenstein

553517/1011/WE3/REV1/G

ivoclar
vivadent®
clinical

English

Description

CRT buffer is used to determine the buffer capacity of saliva by means of a test strip featuring a special indicator system.

Indication

- In vitro diagnostics
- Determination of the buffer capacity of saliva

Application

The test is professionally conducted by the dentist or skilled personnel.

For at least one hour before the test is conducted, patients should

- neither eat nor drink anything,
- not chew any chewing gums,
- not smoke,
- not brush their teeth,
- not use any mouthwashes.

Step-by-step procedure

1. Have your patient sit upright.
2. Stimulate salivation by having your patient chew a paraffin pellet.
3. Collect the saliva in a calibrated container.

Tip: Collect the saliva over a defined period of 5 minutes and determine the salivation rate (ml saliva/min) at the same time.

Findings for adults:

Normal: ≥ 1 ml saliva/min.

Too low: < 0.7 ml saliva/min.

4. Carefully remove the CRT buffer test strip from the packaging without touching the yellow test field.
5. Place test strip on a stable, absorbent paper with the yellow test field facing upwards.
6. Wet the entire yellow test field with saliva using a pipette.
 - Allow saliva to drop off. Do not scratch the test field with the pipet.
 - Prevent the formation of bubbles in the saliva.

7. To determine the buffer capacity of saliva, compare the colour of the test field with the colour samples (pictures) after exactly 5 minutes of reaction time.

Evaluation

Blue indicates a high, yellow a low buffer capacity of saliva.

Note

- If the test field is blotchy, rate the buffer capacity according to the most unfavourable of the colours present.
- If in doubt, repeat the test. Use a new test strip for this purpose.

Disposal

Dispose of used test strips with normal practice waste.

Storage

- Storage temperature: 2–28 °C / 36–82 °F
- Shelf life: see label / packaging
- Do not use CRT buffer after the indicated date of expiration

Keep out of the reach of children.

For use in dentistry only.

Date information prepared: 10/2011, Rev 1

Ivoclar Vivadent AG

FL-9494 Schaan/Liechtenstein

This product has been developed solely for use in dentistry. Processing and application should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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Ivoclar Vivadent Limited
Ivoclar Vivadent, Inc.