**IPS Classic®**

**Instructions for Use**

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**CE 0123**

**ivoclar**  
**vivadent**  
*technical*

# Table of contents | Índice | Índice

3

Material | *El Material* | O Material

5

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
Step-by-Step

# Material

## A sound concept is the key to success

The fully developed IPS Classic metal-ceramic that has clinically proven itself for many years offers a high degree of individuality and creativity. Controlled particle size distribution permits excellent build-up and maintenance of shape of the ceramic material even after several firing cycles.

Generally, IPS Classic is suitable for use with alloys with a CTE of approx 13.8 to 15.2 x 10<sup>-6</sup> /K at 25–500°C. Depending on the type and the composition of the alloy, the firing temperature and the cooling phase in particular have to be taken into con-

sideration. Alloys with a CTE in the lower range (e.g. approx 13.8–14.5 at 25–500°C) can be processed with standard cooling.

IPS Opaque Dentin, IPS Classic Dentin, IPS Classic Incisal:

CTE (25–500°C) [10 <sup>-6</sup> /K] <sup>1)</sup>	12.9 ± 0.5
Flexural strength (biaxial) [MPa] <sup>1)2)</sup>	≥50
Chem. solubility [µg/cm <sup>2</sup> ] <sup>1)</sup>	≤100
Firing temperature [°C]	910–920

<sup>1)</sup> according to ISO 6872:2008

<sup>2)</sup> typical mean value for the flexural strength is 80 MPa

Classification: Dental ceramics Type I / Class 1

## Indication

Conventional multi-layer veneering ceramic for the most popular dental alloys in the CTE range of 13.8 – 15.2 x 10<sup>-6</sup>/K (25 – 500°C)

## Contraindication

- If patients are known to be allergic to any of the ingredients, the material should not be used
- Bruxism
- Veneering of titanium and zirconium oxide frameworks
- Any other use not listed in the indications

## Important processing restrictions

- Exceeding or falling short of the stipulated veneering layer thicknesses
- Failure to observe the layer thickness ratio between the framework and layering ceramic
- Mixing with and processing in conjunction with other dental ceramics
- Veneering of dental alloys not within the stipulated CTE range
- Failure to observe the necessary minimum connector and framework thicknesses
- Use of IPS Shade, IPS Shade V and IPS Stains-P during the layering procedure (e.g. between the dentin and incisal layer)
- Mixing of paste materials with materials in powder form

## Side effects

If patients are known to be allergic to any of the components in the materials, IPS Classic restorations should not be used.

## Important information

Ceramic furnaces of other manufacturers often feature opening mechanisms different from that of Ivoclar Vivadent furnaces. Therefore, the firing conditions may also differ. Make sure that these varying firing conditions are taken into account when working with IPS Classic.

## Warning

Finishing ceramic restorations results in grinding dust. Avoid inhalation. Use suction equipment or protective masks.

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### – Opaquer pastes and glazing materials

contain ceramic material and 25–40 % glycols

### – IPS Classic Powder Opaquer Liquid

Componentes: Water, glycol, acetic acid, additive

### – IPS Classic Build-Up Liquid N, S, L

Componentes: Water, butandiol and chloride

### – IPS Classic Margin Build-Up Liquid

Componentes: Water and cellulose derivative

### – IPS Classic Glaze and Stain Liquid

Componentes: Ethylene glycol

### – IPS Model Sealer

Componentes: Ethyl acetate, nitro-cellulose, softener

### – IPS Ceramic Separating Liquid

Componentes: Paraffin oil

### – IPS Margin Liquid

Componentes: Wax dissolved in hexane

### – IPS Margin Sealer

Componentes: Ethyl acetate, nitro-cellulose, softener

### Warning

- Hexane is highly flammable and detrimental to health. Avoid contact of the material with skin and eyes. Do not inhale vapours. Keep away from sources of ignition.

### Storage

- Protect material in powder form from moisture
- Store material in paste and liquid form at 2–28°C/36–82°F

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**Note: Protect liquids from sunlight and heat!**

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The compatibility of the following alloys and IPS Classic have been thoroughly tested.

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\* Cooling to | Enfriamiento | Esfriamento 800 °C / 1472 °F

\*\*Cooling to | Enfriamiento | Esfriamento 700 °C / 1292 °F

**The product range may vary from country to country.**

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Long-term cooling  
*Enfriamiento lento*  
Esfriamiento lento

### Important

- If the minimum requirements cannot be observed, cooling to \*800 °C, or \*\*700 °C (depending on the alloy type), is required in conjunction with all main firings and glaze firings.
- With ceramic layer thicknesses of over 1.5 mm up to max. 2.5 mm, as well as with voluminous restorations (e.g. implant-retained reconstructions) in combination with high gold and base metal alloys, cooling to \*800 °C or \*\*700 °C must be conducted. This also applies to soldered restorations.

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


# Instructions for Use

## Starting situation

Fabricate a master model or a model with detachable segments on the basis of the impression in the usual manner. It is advisable to apply a sealer to harden the surface and to protect the die. The application of a sealer must not cause any changes in the dimensions of the die. After that, a spacer may be applied in the usual manner.

## Fabricating the framework

When fabricating the framework, make sure that the minimum wall thickness after finishing is 0.3 mm for single crowns and 0.5 mm for bridge abutments. These dimensions are the prerequisite for the stability of the metal framework and the durable bond between the metal and ceramic material. If the stipulated framework and connector dimensions are not observed, the resulting stress within the material will lead to delamination of the ceramic material and distortion of the framework.

## Recommended procedure

1. Complete modelling of the anatomical tooth shape
2. Reduce to make room for the veneering materials
3. Waxing of the contact points and connector areas
4. Checking the occlusal and proximal contact points

## Contouring

The framework reflects the reduced anatomical tooth shape (functional support of the veneering ceramic). In this way, the ceramic material may be applied in even layers and is adequately supported. For that purpose, the properties of the different alloys (e.g. firing stability) must be taken into consideration.

- Undersized metal frameworks result in increased shrinkage of the veneering ceramic and require additional corrective firings.
- If the metal framework is too small, the veneering ceramic is not adequately supported, which may lead to cracks and delamination, particularly in very thick ceramic layers.

## Finishing the metal framework

The cast metal framework is finished using tungsten carbide metal burs or ceramic-bonded grinding instruments. The marginal area of the framework is reduced up to the inner edge of the chamfer or shoulder preparation (labial or circular) to make room for the ceramic shoulder.

- Use reduced pressure when working with softer alloys.
- Work in one direction only to avoid overlapping and inclusions in the metal surface.
- Do not use diamond grinding instruments. Diamond particles may be trapped in the alloy and cause bubbles in the ceramic material during firing.


## Oxide firing

After grinding, carefully blast the framework with  $\text{Al}_2\text{O}_3$  50–100  $\mu\text{m}$ . The required pressure is 1.5–2.0 bar.

- Use only disposable, pure  $\text{Al}_2\text{O}_3$  to blast the alloy surface.
- Observe the Instructions for Use of the alloy manufacturer.

Blasting improves the mechanical bond. It results in the object surface being roughened and considerably enlarged. In order to prevent inclusions of blasting medium residue in the ceramic, we recommend blasting the alloys with the indicated pressure while keeping the nozzle at a flat angle to the object surface. A contaminated metal surface may result in the formation of bubbles during ceramic firing. When finishing the frameworks, the instructions of the alloy manufacturer must be observed at all times. Oxidation is carried out according to the instructions of the alloy manufacturer.

## Schematic diagram of the blasting direction

## Correct angle

for blasting the alloy surface

Before the oxidation firing, clean the metal framework using a brush under running water. Then, thoroughly clean it with steam or in the ultrasonic cleaner. Allow the framework to dry after cleaning. Oxidize the framework according to the instructions given by the alloy manufacturer. Provide ample support for the framework on the firing tray. This is particularly important for long-span bridges. After oxidation, carefully check the framework for porosity or irregular oxide layer. Refinish and reoxidize if necessary.

Please refer to the “Framework Design Guidelines for Metal-Ceramic Restorations” for additional information on framework design. They can be ordered from your Ivoclar Vivadent contact address.

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Apply the 2<sup>nd</sup> opaquer layer so that it covers the entire framework. The fired opaquer should have a silky-mat appearance (egg-shell gloss).

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Mix the powder opaquer thoroughly with the Powder Opaquer Liquid and apply the first opaquer layer (wash) thinly using a brush.

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Apply the 2<sup>nd</sup> opaquer layer evenly in such a way that the metal framework is entirely covered with opaquer. The fired opaquer should have a silky-mat appearance (egg-shell gloss).

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The IPS Classic V Powder Opaquer and Powder Opaquer Liquid are ideally suitable for the application with conventional spray-on techniques. Mix the powder opaquer to a thin consistency, depending on the spray-on system used. Observe the instructions of the manufacturer of the spray-on systems.

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#### Important

- Use distilled water to rewet the mixed or the already applied powder opaquer.
- The firing tray with the opaquerized metal framework should only be placed in the firing chamber and removed from it once the furnace head is completely open and the beeper has sounded.

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#### Important

##### IPS Classic V Powder Opaquer

Alloys (CTE of approx.  $13.8$  to  $15.2 \times 10^{-6}/K$  at  $25-500$  °C) with a solidus point of  $\geq 1080$  °C are suitable for opaquerizing with the powder opaquer at a firing temperature of  $960$  °C.

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Layered IPS Classic Dentin material, where the mamelon shape is only outlined then, Incisal material is applied in small portions and the anatomical shape is completed by means of over-contouring.

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Once the contact points have been provided, the bridge restorations are separated in the interdental spaces down to the opaquer using a sharp instrument (e.g. thin scalpel, razor blade).

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The proximal areas are supplemented with the same materials used for the 1<sup>st</sup> dentin and incisal firing and the restoration is adjusted using Incisal material.

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Providing a true-to-nature surface texture including growth lines and convex/concave areas.

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## Glaze firing / Cocción de glaseado / Queima de Glazeamento

Apply glazing material in the usual manner using a brush.

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Slightly wet the restoration using IPS Classic Glaze and Stain liquid and apply shade adjustments and/or individualized characterizations.

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Remove IPS Classic Glazing material from its container and mix thoroughly. After that, apply glazing material in the usual manner using a brush. Avoid pooling or applying too thick of a layer.

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Completely fired IPS Classic bridge on the model.

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# A-D shades / Cores A-D

Assortment	Materials / Shade groups	A	B	C	D	Special materials
IPS Classic V	Paste + Powder Opaquer (O)	A1 A2 A3 A3.5 A4	B1 B2 B3 B4	C1 C2 C3 C4	D2 D3 D4	GO-Gingiva
	Intensive Paste Opaquer + Powder Opaquer (IO)	IO-A	IO-B	IO-C	IO-A/IO-B	IO-white, IO-violet
	Dentin (D)	A1 A2 A3 A3.5 A4	B1 B2 B3 B4	C1 C2 C3 C4	D2 D3 D4	
	Incisal (S)	S1 S2 S2 S4 S4	S1 S2 S3 S4	S2 S2 S3 S4	S1 S2 S3	
	Transparent (T)	clear, neutral				
IPS Opaque Dentin V	Opaque Dentin (Op.D.)	A1 A2 A3 A3.5 A4	B1 B2 B3 B4	C1 C2 C3 C4	D2 D3 D4	
	Opaque Dentin (Op.D.)	yellow, orange, brown				
IPS Margin V	Margin Material (M)	A1 A2 A3 A3.5 A4	B1 B2 B3 B4	C1 C2 C3 C4	D2 D3 D4	
	Intensive (M)	yellow, orange, brown				
IPS Impulse	Occlusal Dentin (Oc.D.)	orange	yellow	brown	orange, yellow	
	Mamelon Mat. (MM)	MM1, MM2, MM3, MM4, MM orange				
	Incisal (S)	yellow-grey, grey				
	Transparent (T)	yellow-grey, grey, blue				
	Molar Incisal (MS)	MS				
	Incisal Edge Mat. (S)	light yellow, yellow				
IPS Effect	Effect Material	E1 super opal, E2 opal, E3 whitish opal, E4 white opal, E5 red-brown opal				
IPS Gingiva	Gingiva Material	G1, G2, G3, G4, G5 Modifier GM1, GM2, GM3, GM4				
IPS Shade V	Dentin Stains	A1 A2 A3 A3.5 A4	B1 B2 B3 B4	C1 C2 C3 C4	D2 D3 D4	
IPS Stains P	Characterization Stains	white, orange, bamboo, caramel brown, copper brown, cork brown, mahogany, azure, black				basic yellow, basic red, basic blue

# Ivoclar Vivadent – worldwide

## Ivoclar Vivadent AG

Bendererstrasse 2  
9494 Schaan  
Liechtenstein  
Tel. +423 235 35 35  
Fax +423 235 33 60  
[www.ivoclarvivadent.com](http://www.ivoclarvivadent.com)

## Ivoclar Vivadent Pty. Ltd.

1 – 5 Overseas Drive  
P.O. Box 367  
Noble Park, Vic. 3174  
Australia  
Tel. +61 3 9795 9599  
Fax +61 3 9795 9645  
[www.ivoclarvivadent.com.au](http://www.ivoclarvivadent.com.au)

## Ivoclar Vivadent GmbH

Tech Gate Vienna  
Donau-City-Strasse 1  
1220 Wien  
Austria  
Tel. +43 1 263 191 10  
Fax: +43 1 263 191 111  
[www.ivoclarvivadent.at](http://www.ivoclarvivadent.at)

## Ivoclar Vivadent Ltda.

Alameda Caiapós, 723  
Centro Empresarial Tamboré  
CEP 06460-110 Barueri – SP  
Brazil  
Tel. +55 11 2424 7400  
Fax +55 11 3466 0840  
[www.ivoclarvivadent.com.br](http://www.ivoclarvivadent.com.br)

## Ivoclar Vivadent Inc.

1-6600 Dixie Road  
Mississauga, Ontario  
L5T 2Y2  
Canada  
Tel. +1 905 670 8499  
Fax +1 905 670 3102  
[www.ivoclarvivadent.us](http://www.ivoclarvivadent.us)

## Ivoclar Vivadent Shanghai

Trading Co., Ltd.  
2/F Building 1, 881 Wuding Road,  
Jing An District  
200040 Shanghai  
China  
Tel. +86 21 6032 1657  
Fax +86 21 6176 0968  
[www.ivoclarvivadent.com](http://www.ivoclarvivadent.com)

## Ivoclar Vivadent Marketing Ltd.

Calle 134 No. 7-B-83, Of. 520  
Bogotá  
Colombia  
Tel. +57 1 627 3399  
Fax +57 1 633 1663  
[www.ivoclarvivadent.co](http://www.ivoclarvivadent.co)

## Ivoclar Vivadent SAS

B.P. 118  
F-74410 Saint-Jorioz  
France  
Tel. +33 4 50 88 64 00  
Fax +33 4 50 68 91 52  
[www.ivoclarvivadent.fr](http://www.ivoclarvivadent.fr)

## Ivoclar Vivadent GmbH

Dr. Adolf-Schneider-Str. 2  
D-73479 Ellwangen, Jagst  
Germany  
Tel. +49 7961 889 0  
Fax +49 7961 6326  
[www.ivoclarvivadent.de](http://www.ivoclarvivadent.de)

## Wieland Dental + Technik GmbH & Co. KG

Lindenstrasse 2  
75175 Pforzheim  
Germany  
Tel. +49 7231 3705 0  
Fax +49 7231 3579 59  
[www.wieland-dental.com](http://www.wieland-dental.com)

## Ivoclar Vivadent Marketing (India) Pvt. Ltd.

503/504 Raheja Plaza  
15 B Shah Industrial Estate  
Veera Desai Road, Andheri (West)  
Mumbai, 400 053  
India  
Tel. +91 22 2673 0302  
Fax +91 22 2673 0301  
[www.ivoclarvivadent.in](http://www.ivoclarvivadent.in)

## Ivoclar Vivadent s.r.l.

Via Isonzo 67/69  
40033 Casalecchio di Reno (BO)  
Italy  
Tel. +39 051 6113555  
Fax +39 051 6113565  
[www.ivoclarvivadent.it](http://www.ivoclarvivadent.it)

## Ivoclar Vivadent K.K.

1-28-24-4F Hongo  
Bunkyo-ku  
Tokyo 113-0033  
Japan  
Tel. +81 3 6903 3535  
Fax +81 3 5844 3657  
[www.ivoclarvivadent.jp](http://www.ivoclarvivadent.jp)

## Ivoclar Vivadent Ltd.

12F W-Tower, 1303-37  
Seochong-dong, Seocho-gu,  
Seoul 137-855  
Republic of Korea  
Tel. +82 2 536 0714  
Fax +82 2 596 0155  
[www.ivoclarvivadent.co.kr](http://www.ivoclarvivadent.co.kr)

## Ivoclar Vivadent S.A. de C.V.

Av. Insurgentes Sur No. 863,  
Piso 14, Col. Napoles  
03810 México, D.F.  
México  
Tel. +52 55 5062 1000  
Fax +52 55 5062 1029  
[www.ivoclarvivadent.com.mx](http://www.ivoclarvivadent.com.mx)

## Ivoclar Vivadent BV

De Fruittuinien 32  
2132 NZ Hoofddorp  
Netherlands  
Tel. +31 23 529 3791  
Fax +31 23 555 4504  
[www.ivoclarvivadent.com](http://www.ivoclarvivadent.com)

## Ivoclar Vivadent Ltd.

12 Omega St, Rosedale  
PO Box 303011 North Harbour  
Auckland 0751  
New Zealand  
Tel. +64 9 914 9999  
Fax +64 9 914 9990  
[www.ivoclarvivadent.co.nz](http://www.ivoclarvivadent.co.nz)

## Ivoclar Vivadent Polska Sp. z o.o.

Al. Jana Pawła II 78  
00-175 Warszawa  
Poland  
Tel. +48 22 635 5496  
Fax +48 22 635 5469  
[www.ivoclarvivadent.pl](http://www.ivoclarvivadent.pl)

## Ivoclar Vivadent Marketing Ltd.

Prospekt Andropova 18 korp. 6/  
office 10-06  
115432 Moscow  
Russia  
Tel. +7 499 418 0300  
Fax +7 499 418 0310  
[www.ivoclarvivadent.ru](http://www.ivoclarvivadent.ru)

## Ivoclar Vivadent Marketing Ltd.

Qlaya Main St.  
Siricon Building No.14, 2<sup>nd</sup> Floor  
Office No. 204  
P.O. Box 300146  
Riyadh 11372  
Saudi Arabia  
Tel. +966 11 293 8345  
Fax +966 11 293 8344  
[www.ivoclarvivadent.com](http://www.ivoclarvivadent.com)

## Ivoclar Vivadent S.L.U.

Carretera de Fuencarral n°24  
Portal 1 – Planta Baja  
28108-Alcobendas (Madrid)  
Spain  
Telf. +34 91 375 78 20  
Fax: +34 91 375 78 38  
[www.ivoclarvivadent.es](http://www.ivoclarvivadent.es)

## Ivoclar Vivadent AB

Dalvägen 14  
S-169 56 Solna  
Sweden  
Tel. +46 8 514 939 30  
Fax +46 8 514 939 40  
[www.ivoclarvivadent.se](http://www.ivoclarvivadent.se)

## Ivoclar Vivadent Liaison Office

: Tesvikkiye Mahallesi  
Sakayik Sokak  
Nisantas' Plaza No:38/2  
Kat:5 Daire:24  
34021 Sisli – Istanbul  
Turkey  
Tel. +90 212 343 0802  
Fax +90 212 343 0842  
[www.ivoclarvivadent.com](http://www.ivoclarvivadent.com)

## Ivoclar Vivadent Limited

Ground Floor Compass Building  
Feldspar Close  
Warrens Business Park  
Enderby  
Leicester LE19 4SE  
United Kingdom  
Tel. +44 116 284 7880  
Fax +44 116 284 7881  
[www.ivoclarvivadent.co.uk](http://www.ivoclarvivadent.co.uk)

## Ivoclar Vivadent, Inc.

175 Pineview Drive  
Amherst, N.Y. 14228  
USA  
Tel. +1 800 533 6825  
Fax +1 716 691 2285  
[www.ivoclarvivadent.us](http://www.ivoclarvivadent.us)

Manufacturer:

Ivoclar Vivadent AG, 9494 Schaan/Liechtenstein  
[www.ivoclarvivadent.com](http://www.ivoclarvivadent.com)

## Rx only

Date information prepared: 2014-10-28, Rev. 1

These materials have been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes.