Instructions for Use

Gebrauchsinformation

Mode d'emploi

Istruzioni d'uso

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Bruksanvisning

Bruksanvisning

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Bruksanvisning

Productinformatie

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For dental use only!

Caution: U.S. Federal Law restricts this device to sale by or on the order of a licensed dentist.

Made in Liechtenstein

Ivoclar Vivadent AG, FL-9494 Schaan

Complies with/entspricht:

ISO 4049 / EN ISO 4049

CE 0123

ivoclar
vivadent®
clinical

Instructions for Use

Description

Heliomolar HB is a light-curing, radiopaque microfilled composite specifically designed for the restoration of posterior teeth. Heliomolar HB cures with light in the wavelength range of 400–500 nm (blue light). The heavy-body consistency of Heliomolar HB facilitates the establishment of properly shaped proximal contact areas.

Shades

Heliomolar HB is available in the following 11 shades:

A–D

A1

A2

A3

B1

B2

C2

D2

WE

GE

TW*

TG*

* Not radiopaque due to the high translucency

Composition

The monomer matrix consists of Bis-GMA, urethane dimethacrylate and decandiol dimeth-acrylate (22.6 wt%; (Shades TW, TG: 21.0 wt%)). The fillers are composed of highly dispersed silicon dioxide, ytterbium tri-fluoride and copolymers (77.0 wt%; (Shades TW, TG: Fillers are composed of highly dispersed silicon dioxide and copolymers (78.6 wt%)). Additional contents are catalysts, stabilizers, and pigments (0.4 wt%).

The total content of inorganic fillers is 66.7 wt% (Shades TW, TG: 65.4 wt%) or 46.0 vol% (Shades TW, TG: 44.9 vol%). The particle size ranges from 0.04 to 0.2 µm.

Indication

- Restorations in posterior teeth (Classes I and II)
- Restorations in primary teeth

Contraindication

The placement of Heliomolar HB restorations is contraindicated

- if a sufficiently dry operating field cannot be established or if the recommended application technique cannot be applied;
- if the patient is known to be allergic to any of the ingredients in Heliomolar HB.

Side effects

In rare cases, components of Heliomolar HB may cause a sensitizing reaction in patients with a hypersensitivity to any of the material's ingredients. In these cases, the material should not be applied.

To avoid irritation of the pulp, protect deep portions of the cavity preparation in close proximity to the pulp with a suitable pulp/dentin protector (place a calcium hydroxide liner in areas close to the pulp).

Interactions

Materials containing eugenol or clove oil may inhibit the polymerization of composites. Such materials should therefore not be used in conjunction with Heliomolar HB. Cationic mouth rinses and plaque disclosing agents containing chlorhexidine may cause discoloration.

Application procedure

1. Shade selection

Clean the teeth prior to shade selection. The teeth should be moist for shade selection

2. Isolation

Appropriate isolation is essential. The use of a rubber dam is recommended.

3. Cavity preparation

The cavity is prepared according to the principles of the

adhesive technique, i.e. preserving the tooth structure. Do not prepare any sharp internal edges or additional undercuts in caries-free areas. Remove the caries and/or the old restoration, which dictates the proportions of the cavity.

4. Pulp protection / Base / Liner

Do not place a base/liner if not required.

If the cavity is very deep, cover the areas close to the pulp with a calcium hydroxide material and then cover this layer with a pressure resistant cement (e.g. a glass ionomer such as Vivaglass® Liner or a zinc phosphate cement, e.g. PhosphaCEM® PL). Leave the remaining surface of the cavity free for bonding with a dentin-enamel adhesive such as Syntac®, Excite® or AdheSE®.

5. Placement of matrix / interdental wedge

Place a matrix in cavities extending into the proximal area and secure it with a wedge.

6. Conditioning / Application of the bonding agent

Conditioning and application of the bonding agent according to the Instructions for Use of the product in use. Given the coordination of Ivoclar Vivadent materials, we recommend using Syntac (with phosphoric acid etching) or Excite (with phosphoric acid etching) or the AdheSE self-etching adhesive.

7. Application of Heliomolar HB

To improve the marginal adaptation of the high-viscosity composite material, we recommend coating the cavity first with a layer of Helio-molar Flow that should be maximum 1 mm thick.

Apply Heliomolar HB in increments of max. 2.0 mm and adapt the material to the cavity walls with a suitable instrument (e.g. P1 instrument). Light cure each increment with a polymerization unit $> 500 \text{ mW/cm}^2$ (e.g. bluephase®) for 20 seconds*. Hold the light emission window of the curing light as close to the restorative as possible. If a metal matrix has been applied, light cure from the bucco-proximal and/or linguo-/palato-proximal aspect after having removed the matrix.

The shade range of Heliomolar HB includes two transparent enamel shades which may be used as top layer to further enhance the esthetic appearance of the final restoration.

* If the exposure time is 40 s, increments of 2.5 mm (A2, A3, B2, C2, GE) and 3.0 mm (A1, B1, D2, WE, TW, TG,) may be applied.

8. Finishing / Checking the occlusion / Polishing

Remove excess material with suitable finishers (e.g. Astropol® F) or fine diamonds after polymerization. Remove proximal excess with diamond or tungsten carbide finishers or finishing strips. Check occlusion and articulation and make appropriate grinding adjustments to prevent premature contacts or undesired occlusal deflections on the surface of the restoration. Polish the restoration to a high gloss by means of silicone polishers (e.g. Astropol P/ Astropol HP, Astrobrush®), polishing discs and polishing strips.

Additional information

1. If necessary, additional Heliomolar HB can be directly applied to polymerized material. If the Heliomolar HB restoration has already been polished, it must first be roughened and wetted with Heliobond before fresh Heliomolar HB is applied.
2. Heliomolar HB should be applied at room temperature. Cool material may be difficult to dispense.

Warning

- Avoid contact of unpolymerized Heliomolar HB with the skin / mucous membrane or eyes. Unpolymerized Heliomolar HB may cause a slight irritating effect and may lead to a sensitization against methacrylates.
- Commercial medical gloves do not provide protection against the sensitizing effect of methacrylates.

Storage

- Do not use Heliomolar HB after the indicated date of expiration.
- Storage temperature: 2–28 °C (36–82 °F)
- Close syringes and Cavifils immediately after use. Exposure to light causes premature polymerization.
- Shelf life: see note on label and packaging

**Keep out of the reach of children.
For use in dentistry only.**

Date information prepared:
02/2005

Manufacturer:
Ivoclar Vivadent AG
9494 Schaan / Liechtenstein

This material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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Ivoclar Vivadent AB
Ivoclar Vivadent UK Limited
Ivoclar Vivadent, Inc.