## Instructions for Use

### Bruksanvisning Brugsvejledning Käyttöohjeet Bruksanvisning Productinformatie Οδηγίες Χρήσεως

- Dual-curing/light-curing resin based dental luting material
- Dual- og ljusthærdende resinbaseret tandfyllnadsmaterial
- Dualhærdende/lyshærdende dental plastretentionscement
- Kaksoiskovetteinen/valokovetteinen yhdistelmämuovi-sementti hammaslääketieteelliseen käyttöön
- Dualherdende / lysherdende odontologisk sementeringsmateriale
- Dualuithardend/lichtuithardend tandheelkundig bevestigingscomposiet
- Ρητινώδης κοιλία διπλού πολυμερισμού/ φωτοπολυμεριζόμενη

For dental use only.

Caution: U.S. Federal Law restricts this device to sale by or on the order of a licensed dentist.

Made in Liechtenstein

Ivoclar Vivadent AG, FL-9494 Schaan/Liechtenstein

CE 0123

ivoclar  
vivadent®

## Instructions for Use

### Description

Variolink II is a dual-curing (light- and self-curing) luting composite system for the adhesive luting of ceramic, and composite restorations.

Variolink II may also be applied in the light-curing technique only (e.g. luting of veneers). For this purpose, only the Variolink II Base is to be used.

Variolink II is based on the Advanced Composite Technology of Tetric®Ceram. The special filler composition of this technology gives Variolink II excellent physical properties, high abrasion resistance, excellent radiopacity, very good optical properties, as well as continuous fluoride release. Variolink II is available in 6 base shades and 2 catalyst shades as well as in 3 degrees of viscosity. Thus, an appropriate version of Variolink II is available for the various clinical requirements.

### Shades

Base	Chromascop	A-D
bleach XL	010	-
transparent	T	-
white	110	A1
white opaque	-	-
yellow (universal)	210	A3
brown	340	A4
<b>Catalyst high/low viscosity:</b>		
yellow (universal)	210	A3
transparent	T	-
<b>Catalyst ultra viscosity</b>		
yellow (universal)	210	A3

### Working time

About 3.5 min. at 37 °C/99 °F

## Mixing ratio

Mix base paste and catalyst in a 1:1 ratio.

## Composition

The monomer matrix is composed of Bis-GMA, urethane dimethacrylate, and triethylene glycol dimethacrylate. The inorganic fillers are barium glass, ytterbium trifluoride, Ba-Al-fluorosilicate glass, and spheroid mixed oxide.

Additional contents: catalysts, stabilizers, and pigments.

The particle size is 0.04–3.0  $\mu\text{m}$ . The mean particle size is 0.7  $\mu\text{m}$ .

	Monomer %wt.	Filler %wt.	Filler %vol.
Base	26.3	73.4	46.7
Catalyst low viscosity	27.9	71.2	43.6
Catalyst high viscosity	22.0	77.2	52.0

## Indication

Adhesive luting of ceramic, and composite restorations (inlays, onlays, crowns, metal-free adhesive bridges, veneers (SR Adoro<sup>®</sup>, IPS Empress<sup>®</sup>), core build-ups), made in the laboratory or chairside.

## Contraindication

The application of Variolink II is contraindicated:

- if a dry working field cannot be established, or if the stipulated technique cannot be applied
- if the patient is known to be allergic to any of the ingredients of the luting composite


## Side effects

Systemic side effects not known to date. In individual cases, allergic reactions have been reported. Irritation of the pulp may occur in very deep or medium cavities, without appropriate pulp/dentin protection or base/liner.

## Interactions

Phenolic substances (e.g. eugenol) inhibit polymerization. Consequently, application of materials containing these substances is to be avoided.

## Flowchart for adhesive luting with Variolink II

## Try-in and adhesive placement of tooth-coloured restorations

### Step by Step

For detailed information also refer to the Instructions for Use of the respective products.

#### 1. Removal of the temporary restoration

Remove possible residue of the temporary luting composite from the cavity or the prepared die with a polishing brush and an oil- and fluoride-free cleaning paste. Rinse with water spray.

#### Note:

Use only eugenol-free temporary luting agents (e.g. Provilink®), as eugenol may inhibit the polymerization of Variolink II. When working on inlays and using Systemp®.inlay/Systemp®.onlay as a temporary restorative material, the application of a temporary luting composite is generally not necessary.

#### 2. Try-in of the restoration

To aid in achieving optimum final restoration results, try-in of the restoration with a try-in paste (e.g. Variolink II Try-In). To avoid fracture of the restoration, do not check the occlusion at this stage. If necessary, make corrections with a fine diamond at medium speed and with slight pressure. Polish ground surfaces.

#### 3. Isolation

When applying the adhesive luting technique with composites, isolation with a rubber dam is recommended.

#### Restoration – pre-treatment

#### Note:

The following steps “Restoration 4a and 5a” may also be done by the dental assistant in order to relieve the dentist.

#### 4a. Cleaning

Rinse off the glycerine gel with water and dry the restoration with water- and oil-free air.

#### 5a. Pre-treatment of the restoration

According to the instructions of the manufacturer.

#### Cavity / Preparation – pre-treatment

#### 4b. Cleaning

Clean cavity/preparation with water spray and dry.

#### 5b. Pre-treatment of the cavity

- Apply phosphoric acid gel (37% e.g. Total Etch) to the prepared enamel and then flow the etchant onto the prepared dentin. Use a disposable brush or syringe tip to push the etchant into any restricted areas of the preparation, leave for 15-30 seconds on enamel and for 10–15 seconds on dentin.
- Remove all etchant gel with av vigorous water spray for at least 5 seconds. Excess water is removed with either:
  - a) the application of 1–2 bursts of clean, dry air or
  - b) high-volume evacuation tip directly over the preparation surface for 1–2 seconds, or
  - c) blotting with a clear applicator or lint-free absorbent material.
- **Do not overdry dentin**, moist preparation surfaces should be apparent.
- The etched enamel surface should have a chalky white appearance. If not, or if the etched enamel surface has been contaminated, repeat etching.
- Subsequently, apply the adhesive (see instructions for use Excite DSC, Syntac)

#### Note:

Using a matrix band in the proximal regions inhibits accidental etching of the adjacent teeth and facilitates the subsequent removal of excess Variolink II. Hold matrix band in place with a suitable interdental wedge.

#### 6. Variolink II

Mix Variolink II in a 1:1 ratio on a mixing pad for 10 s (mix carefully with a spatula). The working time of the mixed Variolink II is about 3.5 min., at a temperature of 37 °C/99 °F.

## Attention

- Variolink II is a light-/dual-curing material and thus sensitive to blue light (operatory light, ambient light). Therefore, Variolink II should be mixed immediately before application. Avoid intensive light during application.
- The result of mixing the base pastes (bleach XL 010, white (110/A1), white opaque, brown (340/A4) with the catalyst pastes is a mixed colour.
- When mixing is too hard, air may be trapped in the paste. This may delay polymerization and compromise mechanical properties.
- Make sure to remove excess material in time, especially in areas that are difficult to reach (proximal or gingival margins). Given the hardness of fully polymerized luting composites like Variolink II, the removal of excess with rotary or oscillating instruments (finishing diamonds, discs, Proxoshape, etc.) is extremely difficult and time consuming.

**We recommend** the base with the low-viscosity catalyst paste for the placement of crowns and partial crowns, the base with the high-viscosity catalyst paste for inlays/onlays and veneers. The dentist, however, is free to choose the proper viscosity for each indication.

Should the tooth substance be discoloured, we recommend using shade white opaque, since this opaque version of Variolink II has better masking properties.

As veneers demonstrate only limited thickness and thus high translucency, Variolink II can also be applied with **light-curing only**. In that case, only the base paste is used to place the restoration. The catalyst paste is un-necessary. Consequently, a selection between different consistencies is not possible in this case. For controlled direct placement of Variolink II base utilize a Variolink II application tip.

## 7. Placement

### a) Inlays/Onlays:

Apply mixed Variolink II with a brush or a spatula in the cavity and/or, if necessary (with concave shapes or to avoid trapping air), on the inner surface of the restoration.

### b) Crowns:

Apply mixed Variolink II on the inner surface of the restoration and/or, if necessary (to avoid trapping air), on the preparation.

### c) Veneers:

Apply mixed Variolink II with a brush or a spatula on the preparation and/or, if necessary, on the inner surface of the restoration. For direct placement of base alone utilizing the Variolink II application tip.

## Attention

- First, place restoration in situ with slight pressure and remove excess with a brush or another suitable instrument.
- When working on large restorations, make sure to remove excess material in time, especially in areas that are difficult to reach (proximal or gingival margins).
- Increase pressure and maintain for several seconds. Remove the Variolink II excess thus forced out with a brush.
- To aid in the clean-up process maintain the pressure and tack the restoration on an isolated surface area with light-polymerization for 10 to 20 seconds.

## Important:

At this stage, do not polymerize any proximal or marginal regions.

- After tacking the restoration in place remove any additional marginal excess with a spatula or other suitable instrument. Make sure that Variolink II excess is not brushed out of the margins.

## 8. Polymerization

Polymerize Variolink II step by step for at least 40 s per segment. Begin with the proximal margins (polymerize through transparent interdental wedges).

## Remarks regarding oxygen inhibition

As any other composite, Variolink II is subject to oxygen inhibition. This means that the surface layer (approx. 50 µm) does not polymerize during light-curing, as it comes in contact with atmospheric oxygen. This problem can be solved in two ways:

- 1) Do not remove the small amounts of excess Variolink II before polymerization. Thus, only the excess is affected by oxygen inhibition and may easily be removed after polymerization during finishing of the margins.
- 2) Cover restoration margins with a glycerine gel (e.g. Liquid Strip) after removal of excess, but before the polymerization. Thus, oxygen inhibition can be avoided.

### 9. Finishing/Polishing

- Remove polymerized excess with a finishing diamond (grain size below 25 µm), and flexible discs.
- Apply finishing and polishing strips in the proximal regions.
- Check occlusion and functions and make corrections if necessary.
- Polish restoration margins with silicone polishers (Politip-F/-P, Astropol®) or discs.

### 10. Topical fluoride application

Fluoride the tooth after completion of the restoration (e.g. with Fluor Protector).

#### Warning:

Avoid contact of Variolink II with skin, mucous membrane, and eyes. Unpolymerized Variolink II may cause slight irritations and may lead to a sensitization against methacrylates.

#### Storage

- Do not use Variolink II after expiry date.
- Cool storage (2–8 °C/36–46 °F) is indicated for the catalyst paste.
- **Note:** Syringes should be closed immediately after use. Exposure to light causes premature polymerization.
- Shelf life: see expiry date

**Keep material out of the reach of children!**

**For use in dentistry only!**

**Date information prepared:**

04/2003

#### Manufacturer:

Ivoclar Vivadent AG

FL-9494 Schaan / Liechtenstein

This material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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## **Ivoclar Vivadent UK Limited**
## **Ivoclar Vivadent, Inc.**
## **Dentigenix Inc.**