

Vivaglass CEM PL



Scientific Documentation



Table of contents

1	Introduction.....	3
2	Technical Data Sheet.....	5
3	Physical properties	6
4	Biocompatibility	8
5	Literature	8

1 Introduction

Glass ionomers were first described by Alan Wilson and Brian Kent in the late 1960s [1]. The first commercially manufactured dental cements that were based on their findings appeared on the market as early as in 1975 [2]. However, the technology employed was not as advanced as that of the glass ionomers used for conventional cementation in dentistry today.

Originally, glass ionomer cements were composed of an aqueous solution of polyacrylic acid, which reacted with calcium aluminium fluorosilicate glass powder. These first cements were very opaque, demonstrated a long setting time and were very sensitive to moisture. Due to their biocompatibility, these cements were also employed in non-dental fields, eg as components of bone or cartilage replacement material [3].




Figure 1: First step of the neutralization reaction

The properties of the various glass ionomers, such as handling properties and compressive strength, can be modified via the following parameters [4]:

- Molar weight of the employed polymer
- Concentration of the acid solution
- Powder/liquid ratio
- Use of complexing agents (e.g. tartaric acid)

When mixing the powder and liquid, the dissolved polymer acid reacts with the alkaline surface of the glass. This reaction is called a “neutralization reaction”, which is the reaction between an acid and a base that produces water and a salt. The acid first reacts with the calcium ions of the glass (see fig 1). The formation of salt with the aluminum ions takes place slightly more slowly. However, the bond that forms is much stronger and the resulting aluminum carboxylates are not soluble in water. Gradually, a thoroughly cross-linked inorganic structure forms. In the process, the density and viscosity of the cement progressively increases (see fig 2).




Figure 2: Increase of viscosity in the course of neutralization (according to [5])

For glass ionomer cements, aluminosilicate-type glasses are used exclusively. They are obtained by sintering the components at temperatures between 1200 °C and 1550 °C (2192 °F - 2822 °F). By cooling them down suddenly very porous structures form, which are milled into fine particles with a size smaller than 45 µm or even smaller than 15 µm. The following typical compositions are used: $\text{SiO}_2 - \text{Al}_2\text{O}_3 - \text{CaO}$ or $\text{SiO}_2 - \text{Al}_2\text{O}_3 - \text{CaF}_2$. Glass ionomers which contain the latter ingredients are capable of releasing fluoride and thus demonstrate a caries-inhibiting effect.

2 Technical Data Sheet

Standard – Composition (in weight %)

<u>Powder</u>	Ionomer glass	72.0
	Polyacrylic acid	28.0
	Pigments	< 0.1
<u>Liquid</u>	Water	85.0
	Tartaric acid	15.0
	Parabene	< 0.3

Physical properties

In accordance to ISO/FDIS 9917:2002 – Water-based cements

Mixing ratio: Powder/Liquid 3:1

Film thickness	15 ± 5 μm
Setting time (at 37 °C, 100 % rel. humidity)	6.0 – 7.0 Min.
Compressive strength	110 ± 5 MPa
Acidic solubility	0.16 ± 0.1 mm
Lead content	< 10 ppm

3 Physical properties

Vivaglass CEM PL stands out from among conventional glass ionomer cements because of its good radiopacity and aesthetic appearance (see fig 3). Its solubility in water and acid meets the requirements of valid standards. With a shear bond strength of 6 MPa, Vivaglass CEM PL exhibits an exceptionally high dentin bonding ability when compared to other glass ionomers.




Figure 3: Comparison of the transparency of several conventional cements




Figure 4:
Shear bond strength on dentin after 24-h immersion in H₂O




Figure 5:
Radiopacity of different glass ionomer cements




Figure 6: Compressive strength

A compressive strength of 115 MPa clearly surpasses the 75 MPa prescribed by EN ISO 9917.

This feature as well as the comparatively high dentin bonding value ensure reliable adhesion between the restoration and the dental hard tissues.

During the reaction of the carboxylic acid with calcium aluminofluorosilicate glass, fluoride anions are continuously released from the glass (Figure 7).




Figure 7: Measurement of the fluoride release in water of different cements

Summary

Vivaglass Cem PL is a traditional glass ionomer cement. Comparable products have been marketed for 30 years with excellent clinical success rates [6]. The properties of Vivaglass Cem compare favourably with those of standard products, and the product meets the standard requirements in all the important parameters (or shows even better values than those prescribed by ISO 9917).

Even though no extensive clinical studies have been carried out so far, a good clinical performance of the product is expected, judging from the data available and the initial clinical experiences gathered.

4 Biocompatibility

The biocompatibility of glass ionomer cements has been confirmed in scientific studies [6].

All substances used are typical components of glass ionomer cements that have been employed in dentistry for years. Carcinogenic, mutagenic or teratogenic effects are not known to date [7].

Compared with standard products, no increased toxicological risk is to be expected with Vivaglass Cem PL based on the present state of knowledge.

5 Literature

- [1] B. E. Kent, A. D. Wilson; Br. Dent. J. 135, 322 (1973)
- [2] G. J. Mount; A Colour Atlas of Glass-Ionomer Cements; 2nd ed.; London; Martin Dunitz; 1994
- [3] I. M. Brook, G. T. Craig, D. J. Lamb; Clin. Mater. 7, 295 (1991)
- [4] J. W. Nicholson; Biomaterials 19, 485 (1998)
- [5] E. A. Wasson; PhD thesis, Brunnel Univerisity, Uxbridge, 1992 (siehe [4])
- [6] J. W. McLean, A. D. Wilson; Glasionomerzemente; Quintessenz Verlags-GmbH 1988
- [7] Y. Li; J. Dent. Res. 69, 1188 (1990)

This documentation contains a survey of internal and external scientific data ("Information"). The documentation and Information have been prepared exclusively for use in-house by Ivoclar Vivadent and for external Ivoclar Vivadent partners. They are not intended to be used for any other purpose. While we believe the Information is current, we have not reviewed all of the Information, and we cannot and do not guarantee its accuracy, truthfulness, or reliability. We will not be liable for use of or reliance on any of the Information, even if we have been advised to the contrary. In particular, use of the information is at your sole risk. It is provided "as-is", "as available" and without any warranty express or implied, including (without limitation) of merchantability or fitness for a particular purpose.

The Information has been provided without cost to you and in no event will we or anyone associated with us be liable to you or any other person for any incidental, direct, indirect, consequential, special, or punitive damages (including, but not limited to, damages for lost data, loss of use, or any cost to procure substitute information) arising out of your or another's use of or inability to use the Information even if we or our agents know of the possibility of such damages.

Ivoclar Vivadent AG
Research & Development
Scientific Service
Bendererstrasse 2
9494 Schaan
Liechtenstein

Contents: Dr. Thomas Völkel
Edition: April 2003
