IPS Ivocolor[®]

EN Short Instructions

- Stains and glazes

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- Malfarben und Glasuren

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Date information prepared:
2020-03-17 / Rev.0
740009/WE1

The logo graphic features a series of seven dots of varying sizes arranged in a slight arc above the text.

ivoclar
vivadent[®]
passion vision innovation

Intended use

Intended purpose

Ceramic materials for dental restorations

Use

For dental use only!

Description

IPS Ivocolor® is a universal stains and glaze assortment for the individualized staining and characterization of ceramic materials. The range of products has been coordinated with the layering, press and CAD ceramics from Ivoclar Vivadent and can be used regardless of the CTE of the ceramic.

Technical data

Characteristics	Unit	Specifications	
Flexural strength	MPa	≥ 50*	
Chemical solubility	µg/cm ²	≤ 100*	
Coefficient of thermal expansion (25 – T _G °C)	10 ⁻⁶ K ⁻¹	Glaze	9.3 ± 0.5*
		Glaze FLUO / Essence / Shade pigment content ≤ 4%	9.2 ± 0.5*
		Essence / Shade pigment content > 4 to ≤ 12%	9.1 ± 0.5*
		Essence / Shade Pigment content > 4 to ≤ 12%	8.8 ± 0.5*
Glass transition temperature T _G	°C	460 ± 20*	
Radioactivity U ²³⁸	Bq/g	≤ 1*	

* The performance criteria defined in EN ISO 6872:2015 Dentistry – Ceramic materials (ISO 6872:2015) were met.
Classification: Dental ceramics Type II / Class I

Indication

Staining and/or characterization and glazing of ceramic materials from Ivoclar Vivadent:

- Veneering ceramics
- Press ceramics
- Milling ceramics
- Zirconium oxide (frameworks and full-contour restorations)

Contraindication

- Any other use not listed in the indications.
- IPS Ivocolor is not suitable for staining, characterizing and glazing of partially crystallized IPS e.max® CAD restorations (in their blue state) or for IPS e.max CAD Veneering Solutions.
- If the patient is known to be allergic to any of the constituents, IPS Ivocolor must not be used.

Side effects

There are no known side effects to date.

Hypersensitivity to the product or its constituents cannot be excluded in individual cases.

Composition

- **IPS Ivocolor Shade**
Components: Alkali aluminosilicate glass
Additional components: Solvent, pigments
- **IPS Ivocolor Essence**
Components: Alkali aluminosilicate glass
Additional components: Pigments
- **IPS Ivocolor Glaze Powder/FLUO**
Components: Alkali aluminosilicate glass
- **IPS Ivocolor Glaze Paste/FLUO**
Components: Alkali aluminosilicate glass
Additional components: Solvent
- **IPS Ivocolor Mixing Liquids/Essence Fluid**
Components: Solvent

Application

Processing




- **IPS Ivocolor Shade**
IPS Ivocolor Shade Dentin and Shade Incisal are intended for the surface staining of the restorations in the desired A-D tooth shade. The ready-to-use pastes feature a gel structure and have to mixed or diluted, if required, with the IPS Ivocolor Mixing Liquids. To facilitate the application, wetting the restoration surface with IPS Mixing Liquid before staining is recommended. IPS Ivocolor is then applied in a thin layer.

- **IPS Ivocolor Essence**
IPS Ivocolor Essence stains are intended for the individualized surface characterization of restorations as well as for flushing in layering ceramics.
The powders may be used either pure or they can be mixed with each other in any combination. IPS Ivocolor Mixing Liquids or IPS Ivocolor Essence Fluid can be used for mixing with IPS Ivocolor (when mixing for the first time). IPS Ivocolor is then applied in a thin layer.
- **IPS Ivocolor Glaze Powder/FLUO**
IPS Ivocolor Glaze Powder/FLUO are intended for glazing the restoration.
The powders are mixed with the IPS Ivocolor Mixing Liquids to the desired consistency. The IPS Ivocolor Mixing Liquid is added dropwise to the Glaze Powder and mixed to a homogeneous paste. If mixed correctly, the consistency of the glazing material must not be too thin/fluid in order to ensure the gloss effect of the fired restoration. The glaze material is then applied in a thin layer. The application of a thicker glaze layer is recommended for zirconium oxide restorations.
- **IPS Ivocolor Glaze Paste/FLUO**
IPS Ivocolor Glaze Paste/FLUO are intended for glazing the restoration.
The gel-type pastes are mixed with the IPS Ivocolor Mixing Liquids to the desired consistency. If mixed correctly, the consistency of the glazing material must not be too thin/fluid in order to ensure the gloss effect of the fired restoration. The glaze material is then applied in a thin layer. The application of a thicker glaze layer is recommended for zirconium oxide restorations.
- **IPS Ivocolor Mixing Liquid allround**
Liquid for mixing with all IPS Ivocolor powders and pastes.
- **IPS Ivocolor Mixing Liquid longlife**
Liquid for mixing with all IPS Ivocolor powders and pastes.
The viscous consistency and slow evaporation of the liquid results in a higher viscosity of the mixed powders or pastes and thus prolongs the processing time.
- **IPS Ivocolor Essence Fluid**
Liquid for mixing with IPS Ivocolor Essence powder to achieve a paste-like consistency.
IPS Ivocolor Essence Fluid is only suitable for the initial mixing of IPS Ivocolor Essence.

Note

The consistency of the material affects the application behaviour and therefore the firing results. If the glazing material is applied correctly, the restoration demonstrates an even gloss after firing and the desired surface structure is maintained. More intensive shades are achieved by several staining procedures and repeated firing, not by applying thicker layers.

Component Combination Table

	Initial Mixing			Refreshing of Dried Out Materials		
	IPS Ivocolor Mixing Liquid allround	IPS Ivocolor Mixing Liquid longlife	IPS Ivocolor Essence Fluid	IPS Ivocolor Mixing Liquid allround	IPS Ivocolor Mixing Liquid longlife	IPS Ivocolor Essence Fluid
Components						
 IPS Ivocolor Essence	✓	✓	✓	✓	✓	✗
 IPS Ivocolor Shade	✓	✓	✗	✓	✓	✗
 IPS Ivocolor Glaze Powder/FLUO	✓	✓	✗	✓	✓	✗
 IPS Ivocolor Glaze Paste/FLUO	✓	✓	✗	✓	✓	✗

Shade Combination Table IPS Ivocolor® Shade

Combination table for the staining of IPS veneering ceramics and all-ceramic materials of the IPS e.max® System*

	A-D Tooth Shade																			
	BL1	BL2	BL3	BL4	A1	A2	A3	A3.5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4
Shade 0																				
Shade 1																				
Shade 2																				
Shade 3																				
Shade 4																				
Shade 5																				
Shade 6																				
Shade 7																				
Shade 8																				

* Other shade combinations apply for IPS e.mx ZirCAD LT discs in LT group shades.

Firing Parameters

Stain and Glaze firing in the staining and cut-back technique (monolithic/partially veneered)

	Stand-by temperature B [°C/°F]	Closing time* S [min]	Heating rate t↗ [°C/°F/min]	Firing temperature T1 [°C/°F]	Holding time H [min]	Vacuum 1 [°C/°F]	Vacuum 2 [°C/°F]	Long-term cooling** L [°C/°F]	Cooling rate tl [°C/°F/min]
IPS e.max® Press	403/757	IRT 6	60/108	710/1310	1	450/842	709/1308	0	0
IPS e.max® CAD	403/757	IRT 6	60/108	710/1310	1	450/842	709/1308	0	0
IPS e.max® ZirPress	403/757	IRT 6	60/108	710/1310	1	450/842	709/1308	450/842	0
IPS Empress® Esthetic	403/757	IRT 6	60/108	710/1310	1	450/842	709/1308	0	0
IPS Empress® CAD	403/757	IRT 6	60/108	710/1310	1	450/842	709/1308	0	0
IPS InLine® PoM	403/757	IRT 6	60/108	710/1310	1	450/842	709/1308	0	0
IPS e.max® ZirCAD	403/757	IRT 6	45/81	710/1310	1	450/842	709/1308	450/842	0

*IRT standard mode ** Long-term cooling for the last firing

Note: Long-term cooling is required for the final Stain and/or Glaze firing of IPS e.max ZirCAD veneered with IPS e.max Ceram or pressed over with IPS e.max ZirPress.

If several IPS e.max ZirCAD restoration units (multi-unit bridges with solid bridge pontics or several full-contour restorations) are fired with the same firing cycle, the thorough heating of the objects to be fired may take longer. This can be counteracted by reducing the heating rate or prolonging the holding time at the final temperature.

Stain and Glaze firing in the layering technique (full-contour)

	Stand-by temperature B [°C/°F]	Closing time* S [min]	Heating rate t [°C/°F/min]	Firing temperature T1 [°C/°F]	Holding time H [min]	Vacuum 1 [°C/°F]	Vacuum 2 [°C/°F]	Long-term cooling** L [°C/°F]	Cooling rate tl [°C/°F/min]
IPS e.max® Ceram/ IPS e.max® ZirCAD	403/757	IRT 6	60/108	710/1310	1	450/842	709/1308	450/842	0
IPS e.max® Ceram/ IPS e.max® Press	403/757	IRT 6	60/108	710/1310	1	450/842	709/1308	0	0
IPS e.max® Ceram/ IPS e.max® CAD	403/757	IRT 6	60/108	710/1310	1	450/842	709/1308	0	0
IPS Style® Ceram	403/757	IRT 6	60/108	750/1382	1	450/842	749/1380	0	0
IPS d.SIGN®	403/757	IRT 6	60/108	800/1472	1	450/842	799/1470	0	0
IPS InLine®	403/757	IRT 6	60/108	830/1526	1	450/842	829/1524	0	0

*IRT standard mode ** Long-term cooling for the last firing

Note: If IPS e.max ZirCAD restorations are veneered with IPS e.max Ceram, long-term cooling to 450°C (842°F) is required for the final Stain and/or Glaze firing cycle.

Important processing restrictions

- Non-observance of the stipulated limitations of use and processing steps may lead to failure.
- **IPS Ivocolor has only been tested and approved in conjunction with Ivoclar Vivadent IPS ceramic materials.** The user is solely responsible for using IPS Ivocolor in conjunction with other ceramic materials.
-  The correct preparation of the surface depends on the restoration material. Procedure according to the Instructions for Use of the respective material.
- The restoration must be free of contaminations and grease residue before staining and glazing.
- A steam jet or ultrasonic bath is suitable for cleaning the restoration.
- The Combination Table is to be observed when processing the IPS Ivocolor components.
- Pooling should be avoided and the material must not be applied too thickly.
- The indicated firing parameters are guidance values. They are valid for the Ivoclar Vivadent furnaces Programat® P310, P510, P710, P300, P500, P700, EP 3010, EP 5010, EP 3000, EP 5000. For furnaces of an older generation (e.g. P80, P90, P95, P100, X1) these temperatures are also guidance values. However, depending on the age of the firing chamber, the values may deviate by +/- 10°C (18°F).
- If a non-Ivoclar Vivadent furnace is used, temperature corrections may be necessary.
- Regional differences in the power supply or the operation of several electronic devices by means of the same circuit may render adjustments of the temperatures necessary.
- To maintain the material properties, open product containers must be resealed immediately after use.

Safety information

- In the case of serious incidents related to the product, please contact Ivoclar Vivadent AG, Bendererstrasse 2, 9494 Schaan/Liechtenstein, website: www.ivoclarvivadent.com, and your responsible competent authority.
- The current Instructions for Use are available in the download section of the Ivoclar Vivadent AG website (www.ivoclarvivadent.com).

Disposal information

- Remaining stocks must be disposed of according to the corresponding national legal requirements.

Shelf life and storage

- Storage temperature:
Shade/Paste/Mixing Liquids/Essence Fluid: 2 – 28°C (36 – 82 °F)
Powder: No special storage conditions required.
- Do not use the product after the indicated date of expiration.
- Expiration date: see note on packaging.

Additional information

Keep material out of the reach of children!

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the product for its suitability and use for any purpose not explicitly stated in the Instructions.

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	Mélange initial			Réhumidification des matériaux séchés		
	IPS Ivocolor Mixing Liquid allround	IPS Ivocolor Mixing Liquid longlife	IPS Ivocolor Essence Fluid	IPS Ivocolor Mixing Liquid allround	IPS Ivocolor Mixing Liquid longlife	IPS Ivocolor Essence Fluid
Composants						
 IPS Ivocolor Essence	✓	✓	✓	✓	✓	✗
 IPS Ivocolor Shade	✓	✓	✗	✓	✓	✗
 IPS Ivocolor Glaze Powder/FLUO	✓	✓	✗	✓	✓	✗
 IPS Ivocolor Glaze Paste/FLUO	✓	✓	✗	✓	✓	✗

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	Mezcla inicial			Reconstitución de materiales secos		
	IPS Ivocolor Mixing Liquid allround	IPS Ivocolor Mixing Liquid longlife	IPS Ivocolor Essence Fluid	IPS Ivocolor Mixing Liquid allround	IPS Ivocolor Mixing Liquid longlife	IPS Ivocolor Essence Fluid
Componentes						
 IPS Ivocolor Essence	✓	✓	✓	✓	✓	✗
 IPS Ivocolor Shade	✓	✓	✗	✓	✓	✗
 IPS Ivocolor Glaze Powder/FLUO	✓	✓	✗	✓	✓	✗
 IPS Ivocolor Glaze Paste/FLUO	✓	✓	✗	✓	✓	✗

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	Mistura inicial			Renovação de materiais secos		
	IPS Ivocolor Mixing Liquid allround	IPS Ivocolor Mixing Liquid longlife	IPS Ivocolor Essence Fluid	IPS Ivocolor Mixing Liquid allround	IPS Ivocolor Mixing Liquid longlife	IPS Ivocolor Essence Fluid
Componentes						
 IPS Ivocolor Essence	✓	✓	✓	✓	✓	✗
 IPS Ivocolor Shade	✓	✓	✗	✓	✓	✗
 IPS Ivocolor Glaze Powder/FLUO	✓	✓	✗	✓	✓	✗
 IPS Ivocolor Glaze Paste/FLUO	✓	✓	✗	✓	✓	✗

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Ivoclar Vivadent Limited
Ivoclar Vivadent, Inc.