


Ivoclar Vivadent All-Ceramics | Vollkeramik



Ivoclar Vivadent All-Ceramics | Vollkeramik

Indications | Indikationen

	IPS e.max®					IPS Empress®	
	IPS e.max® Press	IPS e.max® CAD	IPS e.max® ZirCAD	IPS e.max® ZirPress	IPS e.max® Ceram	IPS Empress® Esthetic	IPS Empress® CAD
Thin Veneers I Dünne Veneers	✓	✓					
Veneers	✓	✓		✓	✓	✓	✓
Occlusal veneers I Okklusale Veneers	✓	✓					
Inlays / Onlays	✓	✓				✓	✓
Partial crowns I Teilkronen	✓	✓			✓ ²⁾	✓	✓
Minimally invasive crowns ⁴ I Minimal-invasive Kronen ⁴	✓	✓	✓				
Crowns I Kronen	✓	✓/✓ ³	✓	✓ ¹⁾	✓ ²⁾	✓	✓
3-unit anterior bridges I 3-gliedrige Frontzahnbrücken	✓	✓/✓ ³	✓	✓ ¹⁾	✓ ²⁾		
3-unit premolar bridges I 3-gliedrige Prämolarenbrücken	✓	✓/✓ ³	✓	✓ ¹⁾	✓ ²⁾		
3-unit posterior bridges I 3-gliedrige Seitenzahnbrücken		✓ ³⁾	✓	✓ ¹⁾	✓ ²⁾		
3- and multiple-unit bridges I 3- und vielgliedrige Brücken		✓ ³⁾	✓	✓ ¹⁾	✓ ²⁾		
Hybrid abutments I Hybrid-Abutments	✓	✓					
Hybrid abutment crowns I Hybrid-Abutment-Kronen	✓	✓					

1) In combination with IPS e.max ZirCAD I in Kombination mit IPS e.max ZirCAD

2) One layering ceramic for all IPS e.max components I eine Schichtkeramik für alle IPS e.max-Komponenten

3) IPS e.max Veneering Solutions (CAD-on technique) in combination with IPS e.max ZirCAD I IPS e.max Veneering Solutions (CAD-on-Technik) in Kombination mit IPS e.max ZirCAD

4) Material thicknesses vary depending on the material I Schichtstärken variieren in Abhängigkeit vom Material



Please read the corresponding Instructions for Use
Bitte entsprechende Gebrauchsinformation beachten

Date information prepared: 2017-05-30
627459