

CAD/CAM blocks

Overview of the range of indications

Indikationen	 IPS Empress® CAD	 IPS e.max® all ceramic all you need		 Telio® CAD
		IPS e.max® CAD	IPS e.max® ZirCAD	
Thin veneers		✓		
Occlusal veneers		✓		
Veneers	✓	✓		
Inlays	✓	✓		
Onlays	✓	✓		
Partial crowns	✓	✓		
Minimal-invasive crowns (≥1 mm)		✓		
Crowns	✓	✓	✓	
3-unit bridges (up to the second premolar as the terminal abutment)		✓		
3-unit bridges			✓	
Crowns and bridges (Veneering Solutions with ZrO ₂ framework)		✓	✓	
Hybrid abutment		✓		
Hybrid abutment crown		✓		
Crown frameworks		✓	✓	
Interlocked crown frameworks			✓	
Multi-unit bridge frameworks			✓	
Temporary crowns				✓
Temporary bridges (2 connected bridge pontics at most)				✓
Temporary hybrid abutment crowns				✓

CAD/CAM partners

The blocks have to be processed with an authorized CAD/CAM system.

Depending on the CAD/CAM system in use, the indications may be restricted.

More information is available on the Internet from www.ivoclarvivadent.com

Preparation date: 2017-04/Rev. 1

www.ivoclarvivadent.com

Ivoclar Vivadent AG

Bendererstr. 2 | 9494 Schaan | Liechtenstein | Tel. +423 235 35 35 | Fax +423 235 33 60