

IVOCLAR VIVADENT TAM SERAMİK ÜRÜNLERİ

 **IPS Empress®**

 **IPS e.max™**

all ceramic
all you need

PRESS
Teknolojisi



CAD / CAM
Teknolojisi



IPS Empress Esthetic Veneer
IPS Empress Universal Shade/Stains



IPS e.max Ceram



ivoclar vivadent®

IVOCLAR VIVADENT TAM SERAMİK ENDİKASYONLARI

	IPS Empress		IPS e.max				
	IPS Empress Esthetic	IPS Empress CAD	IPS e.max Press	IPS e.max ZirPress	IPS e.max ZirCAD	IPS e.max CAD	IPS e.max Ceram
İnce venerler			✓				
Venerler	✓	✓	✓	✓		✓	✓
Oklüzal venerler			✓				
İnleyler/ Onleyler	✓	✓	✓			✓	
Parçalı kuronlar	✓	✓	✓			✓	✓ ²⁾
Anterior kuronlar	✓	✓	✓	✓ ¹⁾	✓	✓	✓ ²⁾
Posterior kuronlar	✓	✓	✓	✓ ¹⁾	✓	✓	✓ ²⁾
3-üyel anterior köprüler			✓	✓ ¹⁾	✓		✓ ²⁾
3-üyel premolar köprüler			✓	✓ ¹⁾	✓		✓ ²⁾
3-üyel posterior köprüler				✓ ¹⁾	✓		✓ ²⁾
4 - 6 üyel köprüler				✓ ¹⁾	✓		✓ ²⁾
7 - 12 üyel köprüler				✓ ¹⁾	✓		✓ ²⁾
İnley köprüler				✓ ¹⁾	✓		

¹⁾ IPS e.max ZirCAD ile birlikte kullanıldığında.

²⁾ Tüm IPS e.max bileşenleri için tek yığma seramiği.