

Programat® P510

Software Revision Log

Version	Changes from previous version	Status
V5.0 28.04.2017	<p>Software optimized</p> <ul style="list-style-type: none"> • Now with DSA-function (Digital Shade Assistant) only available starting from serial number 4010005000 (200V-240V / 50-60Hz) and 4011005000 (110-120V / 50-60Hz) • New program structure with new program parameters for IPS Style, IPS e.max Press Layering Technique, IPS e.max ZirPress Layering Technique, IPS e.max ZirCAD Staining Technique, IPS e.max CAD Crystall/Glaze, IPS e.max CAD Layering Technique, IPS e.max CAD Layering Technique, IPS d.SIGN, IPS e.max ZirCAD Layering Technique • Optimized Vacuum test program • Various software optimizations 	Currently delivered version
V4.0 08.04.2016	<p>Software optimized</p> <ul style="list-style-type: none"> • New programs for IPS Style and IPS Ivocolor • In-process dehumidification: Function that protects the vacuum system from contamination and its consequential damages • Various software optimizations 	No longer available!
V3.1 09.04.2015	<p>Software optimized</p> <ul style="list-style-type: none"> • Various software optimizations 	No longer available!
V3.0 01.10.2014	<p>Software optimized</p> <ul style="list-style-type: none"> • Optimized Programat Infrared Technology: <ul style="list-style-type: none"> ○ Enhanced object detection – particularly of smaller objects ○ IRT Plus mode: especially suitable for opaquer firing cycles and complex restorations ○ New IRT menu with fast access directly from the program screen • New firing programs for IPS InLine Powder Opaquer (P21, P221, P241) • Various software optimizations 	No longer available!

Programat® P510

V2.1 06.06.2014	Software optimized <ul style="list-style-type: none"> The software was attuned to meet the requirements of the Programat App (for iPhones) and various software optimization measures were implemented The response time of the touch screen was sped up Various optimization measures regarding the Programat Infrared technology were integrated 	No longer available!
V2.0 18.11.2013	Software optimized <ul style="list-style-type: none"> New Speed Crystallization program for IPS e.max CAD (Important: The speed program must only be used in combination with the specially designed IPS e.max CAD Speed Crystallization Tray) Additional software improvements 	No longer available!
V1.1 20.09.2013	Basic version	No longer available!

Version 28.04.2017