

Programat® CS2

Software Revision Log

Version	Changes from previous version	Status
V4.1 01.02.2018	Software optimized <ul style="list-style-type: none"> • Push function can be tested • Various languages have been improved • Various software optimizations 	Currently delivered version
V4.0 28.04.2017	Software optimized <ul style="list-style-type: none"> • New program structure with new program parameters • Various software optimizations 	No longer available!
V3.0 08.04.2016	Software optimized <ul style="list-style-type: none"> • In-process dehumidification: Function that protects the vacuum system from contamination and its consequential damages • Various software optimizations 	No longer available!
V2.1 09.04.2015	Software optimized <ul style="list-style-type: none"> • Various software optimizations 	No longer available!
V2.0 01.10.2014	Software optimized <ul style="list-style-type: none"> • Various software optimizations 	No longer available!
V1.1 06.06.2014	Software optimized <ul style="list-style-type: none"> • The software was attuned to meet the requirements of the Programat App (for iPhones) and various software optimizations implemented. • The software was optimized to meet the requirements of the PrograBase X10 Software. As a result, more functions with regard to PrograBase X10 are available with the Programat CS2. 	No longer available!
V1.0 20.09.2013	Basic version	No longer available!

Version 01.02.2018