

# Programat® CS2

## Software Revision Log

Version	Changes from previous version	Status
<b>V4.0</b> 28.04.2017	<b>Software optimized</b> <ul style="list-style-type: none"> <li>• New program structure with new program parameters</li> <li>• Various software optimizations</li> </ul>	<b>Currently delivered version</b>
<b>V3.0</b> 08.04.2016	<b>Software optimized</b> <ul style="list-style-type: none"> <li>• In-process dehumidification: Function that protects the vacuum system from contamination and its consequential damages</li> <li>• Various software optimizations</li> </ul>	<b>No longer available!</b>
<b>V2.1</b> 09.04.2015	<b>Software optimized</b> <ul style="list-style-type: none"> <li>• Various software optimizations</li> </ul>	<b>No longer available!</b>
<b>V2.0</b> 01.10.2014	<b>Software optimized</b> <ul style="list-style-type: none"> <li>• Various software optimizations</li> </ul>	<b>No longer available!</b>
<b>V1.1</b> 06.06.2014	<b>Software optimized</b> <ul style="list-style-type: none"> <li>• The software was attuned to meet the requirements of the Programat App (for iPhones) and various software optimizations implemented.</li> <li>• The software was optimized to meet the requirements of the PrograBase X10 Software. As a result, more functions with regard to PrograBase X10 are available with the Programat CS2.</li> </ul>	<b>No longer available!</b>
<b>V1.0</b> 20.09.2013	<b>Basic version</b>	<b>No longer available!</b>

Version 28.04.2017