

Programat® EP 3010

Software Revision Log

| Version | Changes from previous version | Status |
|---------------------------|---|------------------------------------|
| V4.0 28.04.2017 | Software optimized <ul style="list-style-type: none"> • Now with FPF-function (Fully automatic Press Function) • New program structure with new program parameters for IPS Style, IPS e.max Press Layering Technique, IPS e.max ZirPress Layering Technique, IPS e.max ZirCAD Staining Technique, IPS e.max CAD Crystall/Glaze, IPS e.max CAD Layering Technique, IPS e.max CAD Layering Technique, IPS d.SIGN, IPS e.max ZirCAD Layering Technique • Optimized Vacuum test program • Various software optimizations | Currently delivered version |
| V3.0 08.04.2016 | Software optimized <ul style="list-style-type: none"> • New programs for IPS Style, IPS Ivocolor and IPS e.max Press MT • In-process dehumidification: Function that protects the vacuum system from contamination and its consequential damages • Various software optimizations | No longer available! |
| V2.1 09.04.2015 | Software optimized <ul style="list-style-type: none"> • Various software optimizations | No longer available! |
| V2.0 01.10.2014 | Software optimized <ul style="list-style-type: none"> • New firing programs for IPS InLine Powder Opaquer (P21, P221, P241) • New press programs for IPS e.max Press Multi • Various software optimizations | No longer available! |
| V1.1 06.06.2014 | Software optimized <ul style="list-style-type: none"> • The software was attuned to meet the requirements of the Programat App (for iPhones) and various software optimization measures were implemented | No longer available! |
| V1.0 10.05.2014 | Basic version | No longer available! |

Version 28.04.2017