


**1 Identification**

- **Product identifier**
- **Trade name:** *Virtual Light Body Regular / Fast Base, Virtual Extra Light Body Regular / Fast Base*
- **Application of the substance / the mixture** Dental impression material
- **Details of the supplier of the safety data sheet**
- **Manufacturer/Supplier:**  
Ivoclar Vivadent Inc.  
175 Pineview Drive, Amherst, N.Y. 14228  
USA  
Tel. +1 800 533 6825  
Fax +1 716 691 2285  
  
Ivoclar Vivadent Inc.  
1-6600 Dixie Road  
Mississauga, Ontario  
L5T 2Y2  
Canada  
Phone: +1 905 670 8499  
Fax: +1 905 670 3102
- **Information department:** Quality Assurance / Regulatory Affairs
- **Emergency telephone number:**  
24 Hour Emergency Assistance:  
Emergency-Call USA - Infotrac: 1-800-535-5053  
Emergency-Call Canada - Canutec: 1-613-996-6666  
  
General SDS Assistance:  
US: 1-800-533-6825  
Canada: 1-800-263-8182

**2 Hazard(s) identification**

- **Classification of the substance or mixture**  
Carc. 1A H350 May cause cancer.  
STOT RE 1 H372 Causes damage to the lung through prolonged or repeated exposure. Route of exposure:  
Inhalation.

---

  - **Label elements**
  - **GHS label elements**  
The product is classified and labeled according to the Globally Harmonized System (GHS).
  - **Hazard pictograms**
- 

GHS08
- **Signal word** Danger
  - **Hazard-determining components of labeling:**  
cristobalite  
Quartz (SiO<sub>2</sub>)
  - **Hazard statements**  
May cause cancer.

# Safety Data Sheet

acc. to OSHA HCS

Printing date 10/23/2019

Version number 8

Reviewed on 04/26/2019

Trade name: **Virtual Light Body Regular / Fast Base, Virtual Extra Light Body Regular / Fast Base**

(Contd. of page 1)

Causes damage to the lung through prolonged or repeated exposure. Route of exposure: Inhalation.

· **Precautionary statements**

Get medical advice/attention if you feel unwell.

Dispose of contents/container in accordance with local/regional/national/international regulations.

· **Additional information:**

Medical devices are exempt from the labeling provisions of GHS according to the Hazard Communication Standard (HCS).

· **Classification system:**

· **NFPA ratings (scale 0 - 4)**



· **HMIS-ratings (scale 0 - 4)**

|            |   |                |
|------------|---|----------------|
| HEALTH     | 0 | Health = *0    |
| FIRE       | 0 | Fire = 0       |
| REACTIVITY | 0 | Reactivity = 0 |

· **Other hazards**

· **Results of PBT and vPvB assessment**

· **PBT:** Not applicable.

· **vPvB:** Not applicable.

### 3 Composition/information on ingredients

· **Chemical characterization: Mixtures**

· **Description:** Mixture of the substances listed below with nonhazardous additions.

· **Dangerous components:**

|                 |                            |           |
|-----------------|----------------------------|-----------|
| CAS: 14464-46-1 | crystalite                 | 10-25%    |
| CAS: 14808-60-7 | Quartz (SiO <sub>2</sub> ) | ≥2.5-<10% |

### 4 First-aid measures

· **Description of first aid measures**

· **After inhalation:**

Supply fresh air and to be sure call for a doctor.

In case of unconsciousness place patient stably in side position for transportation.

· **After skin contact:** Rinse with water.

· **After eye contact:** Rinse opened eye for several minutes under running water. Then consult a doctor.

· **After swallowing:** Seek medical treatment.

· **Information for doctor:**

· **Most important symptoms and effects, both acute and delayed** No further relevant information available.

· **Indication of any immediate medical attention and special treatment needed**

No further relevant information available.

US

(Contd. on page 3)

## Safety Data Sheet acc. to OSHA HCS

Printing date 10/23/2019

Version number 8

Reviewed on 04/26/2019

**Trade name: Virtual Light Body Regular / Fast Base, Virtual Extra Light Body Regular / Fast Base**

(Contd. of page 2)

### 5 Fire-fighting measures

- **Extinguishing media**
- **Suitable extinguishing agents:**  
CO<sub>2</sub>, extinguishing powder or water spray. Fight larger fires with water spray or alcohol resistant foam.
- **Special hazards arising from the substance or mixture** No further relevant information available.
- **Advice for firefighters**
- **Protective equipment:** No special measures required.

### 6 Accidental release measures

- **Personal precautions, protective equipment and emergency procedures** Not required.
- **Environmental precautions:** Do not allow to enter sewers/ surface or ground water.
- **Methods and material for containment and cleaning up:** Pick up mechanically.
- **Reference to other sections**  
See Section 7 for information on safe handling.  
See Section 8 for information on personal protection equipment.  
See Section 13 for disposal information.

### 7 Handling and storage

- **Handling:**
- **Precautions for safe handling**  
Only adequately trained personnel should handle this product.  
For use in dentistry only.
- **Information about protection against explosions and fires:** No special measures required.
- **Conditions for safe storage, including any incompatibilities**
- **Storage:**
- **Requirements to be met by storerooms and receptacles:** Store only in the original receptacle.
- **Information about storage in one common storage facility:** Not required.
- **Further information about storage conditions:**  
Keep receptacle tightly sealed.  
Protect from heat and direct sunlight.
- **Specific end use(s)** No further relevant information available.

### 8 Exposure controls/personal protection

- **Additional information about design of technical systems:** No further data; see item 7.
- **Control parameters**
- **Components with limit values that require monitoring at the workplace:**  
The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.
- **Additional information:** The lists that were valid during the creation were used as basis.
- **Exposure controls**
- **Personal protective equipment:**
- **General protective and hygienic measures:**  
Usual hygienic measures for dental practice and dental laboratories.  
Keep away from foodstuffs, beverages and feed.  
Wash hands before breaks and at the end of work.
- **Breathing equipment:** Not required.

(Contd. on page 4)

## Safety Data Sheet

acc. to OSHA HCS

Printing date 10/23/2019

Version number 8

Reviewed on 04/26/2019

**Trade name: Virtual Light Body Regular / Fast Base, Virtual Extra Light Body Regular / Fast Base**

(Contd. of page 3)

· **Protection of hands:**

The glove material has to be impermeable and resistant to the product/ the substance/ the preparation.  
After use of gloves apply skin-cleaning agents and skin cosmetics.

· **Material of gloves**

The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.

· **Penetration time of glove material**

The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.

· **Eye protection:** Safety glasses

### 9 Physical and chemical properties

· **Information on basic physical and chemical properties**

· **General Information**

· **Appearance:**

|                          |                 |
|--------------------------|-----------------|
| <b>Form:</b>             | Viscous         |
| <b>Color:</b>            | Light blue      |
| · <b>Odor:</b>           | like peppermint |
| · <b>Odor threshold:</b> | Not determined. |

· **pH-value:** Not applicable.

· **Change in condition**

|                                     |               |
|-------------------------------------|---------------|
| <b>Melting point/Melting range:</b> | Undetermined. |
| <b>Boiling point/Boiling range:</b> | Undetermined. |

· **Flash point:** Not applicable.

· **Flammability (solid, gaseous):** Not applicable.

· **Auto igniting:** Product is not selfigniting.

· **Danger of explosion:** Product does not present an explosion hazard.

· **Explosion limits:**

|               |                 |
|---------------|-----------------|
| <b>Lower:</b> | Not determined. |
| <b>Upper:</b> | Not determined. |

· **Vapor pressure:** Not determined.

· **Density:** Not determined.

· **Relative density** Not determined.

· **Vapor density** Not determined.

· **Evaporation rate** Not determined.

· **Solubility in / Miscibility with**

**Water:** Not miscible or difficult to mix.

· **Partition coefficient (n-octanol/water):** Not determined.

· **Viscosity:**

|                   |                 |
|-------------------|-----------------|
| <b>Dynamic:</b>   | Not determined. |
| <b>Kinematic:</b> | Not determined. |

(Contd. on page 5)

## Safety Data Sheet

acc. to OSHA HCS

Printing date 10/23/2019

Version number 8

Reviewed on 04/26/2019

**Trade name: Virtual Light Body Regular / Fast Base, Virtual Extra Light Body Regular / Fast Base**

(Contd. of page 4)

· **Other information** No further relevant information available.

### 10 Stability and reactivity

- **Reactivity** No further relevant information available.
- **Chemical stability** Stable under normal handling and storage conditions.
- **Thermal decomposition / conditions to be avoided:** No decomposition if used according to specifications.
- **Possibility of hazardous reactions** No dangerous reactions known.
- **Conditions to avoid** No further relevant information available.
- **Incompatible materials:** No further relevant information available.
- **Hazardous decomposition products:** None under normal conditions of storage and use.

### 11 Toxicological information

- **Information on toxicological effects**
- **Acute toxicity:**
- **on the skin:** No irritant effect.
- **on the eye:** No irritating effect.
- **Sensitization:** No sensitizing effects known.
- **Additional toxicological information:** No further relevant information available.
- **Carcinogenic categories**

#### · IARC (International Agency for Research on Cancer)

|                 |                            |   |
|-----------------|----------------------------|---|
| CAS: 14464-46-1 | crystalite                 | I |
| CAS: 14808-60-7 | Quartz (SiO <sub>2</sub> ) | I |

#### · NTP (National Toxicology Program)

|                 |                            |   |
|-----------------|----------------------------|---|
| CAS: 14464-46-1 | crystalite                 | K |
| CAS: 14808-60-7 | Quartz (SiO <sub>2</sub> ) | K |

#### · OSHA-Ca (Occupational Safety & Health Administration)

None of the ingredients is listed.

### 12 Ecological information

- **Toxicity**
- **Aquatic toxicity:** No further relevant information available.
- **Persistence and degradability** No further relevant information available.
- **Behavior in environmental systems:**
- **Bioaccumulative potential** No further relevant information available.
- **Mobility in soil** No further relevant information available.
- **Additional ecological information:**
- **General notes:**
- **Water hazard class 1 (Self-assessment):** slightly hazardous for water
- **Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.**
- **Results of PBT and vPvB assessment**
- **PBT:** Not applicable.
- **vPvB:** Not applicable.

(Contd. on page 6)

US

## Safety Data Sheet

acc. to OSHA HCS

Printing date 10/23/2019

Version number 8

Reviewed on 04/26/2019

**Trade name: Virtual Light Body Regular / Fast Base, Virtual Extra Light Body Regular / Fast Base**

· **Other adverse effects** No further relevant information available.

(Contd. of page 5)

### 13 Disposal considerations

- **Waste treatment methods**
- **Recommendation:**  
Take to an approved landfill or a waste incineration plant, under conditions approved by the local authority.
- **Uncleaned packagings:**
- **Recommendation:** Disposal must be made according to official regulations.

### 14 Transport information

- **UN-Number**
- **DOT, ADR/RID/ADN, ADN, IMDG, IATA** Void
- **UN proper shipping name**
- **DOT, ADR/RID/ADN, ADN, IMDG, IATA** Void
- **Transport hazard class(es)**
- **DOT, ADR/RID/ADN, ADN, IMDG, IATA**
- **Class** Void
- **Packing group**
- **DOT, ADR/RID/ADN, IMDG, IATA** Void
- **Environmental hazards:**
- **Marine pollutant:** No
- **Special precautions for user** Not applicable.
- **Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code** Not applicable.
- **Transport/Additional information:** Product is not classified as a dangerous good for transport (ADR, IMDG, IATA).
- **UN "Model Regulation":** Void

### 15 Regulatory information

- **Safety, health and environmental regulations/legislation specific for the substance or mixture**

· **Sara**

· **Section 355 (extremely hazardous substances):**

None of the ingredients is listed.

· **Section 313 (Specific toxic chemical listings):**

None of the ingredients is listed.

· **TSCA (Toxic Substances Control Act):**

CAS: 14464-46-1 cristobalite

CAS: 14808-60-7 Quartz (SiO<sub>2</sub>)

CAS: 9016-45-9 Renex 678

(Contd. on page 7)

US

## Safety Data Sheet

acc. to OSHA HCS

Printing date 10/23/2019

Version number 8

Reviewed on 04/26/2019

**Trade name: Virtual Light Body Regular / Fast Base, Virtual Extra Light Body Regular / Fast Base**

(Contd. of page 6)

· **Proposition 65**

· **Chemicals known to cause cancer:**

|                 |                            |
|-----------------|----------------------------|
| CAS: 14464-46-1 | crystalite                 |
| CAS: 14808-60-7 | Quartz (SiO <sub>2</sub> ) |

· **Chemicals known to cause reproductive toxicity for females:**

None of the ingredients is listed.

· **Chemicals known to cause reproductive toxicity for males:**

None of the ingredients is listed.

· **Chemicals known to cause developmental toxicity:**

None of the ingredients is listed.

· **Carcinogenic categories**

· **EPA (Environmental Protection Agency)**

None of the ingredients is listed.

· **TLV (Threshold Limit Value established by ACGIH)**

|                 |                            |    |
|-----------------|----------------------------|----|
| CAS: 14464-46-1 | crystalite                 | A2 |
| CAS: 14808-60-7 | Quartz (SiO <sub>2</sub> ) | A2 |

· **NIOSH-Ca (National Institute for Occupational Safety and Health)**

|                 |                            |
|-----------------|----------------------------|
| CAS: 14464-46-1 | crystalite                 |
| CAS: 14808-60-7 | Quartz (SiO <sub>2</sub> ) |

· **GHS label elements**

The product is classified and labeled according to the Globally Harmonized System (GHS).

· **Hazard pictograms**



GHS08

· **Signal word** *Danger*

· **Hazard-determining components of labeling:**

crystalite

Quartz (SiO<sub>2</sub>)

· **Hazard statements**

May cause cancer.

Causes damage to the lung through prolonged or repeated exposure. Route of exposure: Inhalation.

· **Precautionary statements**

Get medical advice/attention if you feel unwell.

Dispose of contents/container in accordance with local/regional/national/international regulations.

· **National regulations:**

· **Other regulations, limitations and prohibitive regulations**

The product is a medical device according to the Directive 93/42/EEC.

This product is classified as a medical device under US and Canadian regulations and has been reviewed by the US Food and Drug Administration and Health Canada.

· **Chemical safety assessment:** A Chemical Safety Assessment has not been carried out.

US

(Contd. on page 8)

## Safety Data Sheet

acc. to OSHA HCS

Printing date 10/23/2019

Version number 8

Reviewed on 04/26/2019

**Trade name: Virtual Light Body Regular / Fast Base, Virtual Extra Light Body Regular / Fast Base**

(Contd. of page 7)

### 16 Other information

*This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.*

· **Date of preparation / last revision** 10/23/2019 / 7

· **Abbreviations and acronyms:**

*ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)*

*IMDG: International Maritime Code for Dangerous Goods*

*DOT: US Department of Transportation*

*IATA: International Air Transport Association*

*ACGIH: American Conference of Governmental Industrial Hygienists*

*EINECS: European Inventory of Existing Commercial Chemical Substances*

*ELINCS: European List of Notified Chemical Substances*

*CAS: Chemical Abstracts Service (division of the American Chemical Society)*

*NFPA: National Fire Protection Association (USA)*

*HMS: Hazardous Materials Identification System (USA)*

*PBT: Persistent, Bioaccumulative and Toxic*

*vPvB: very Persistent and very Bioaccumulative*

*NIOSH: National Institute for Occupational Safety*

*OSHA: Occupational Safety & Health*

*TLV: Threshold Limit Value*

*PEL: Permissible Exposure Limit*

*REL: Recommended Exposure Limit*

*Carc. 1A: Carcinogenicity – Category 1A*

*STOT RE 1: Specific target organ toxicity (repeated exposure) – Category 1*

· **\* Data compared to the previous version altered.**

US