

Safety Data Sheet

Section 1: Identification of the material and the Supplier

Product identifier	Multilink Primer A
Product use	Self-curing and self-etching primer
Distributor information	
Company identification	Ivoclar Vivadent (NZ) Ltd
Address	12 Omega Street, Rosedale, Auckland, New Zealand
Telephone Number	0508 486 252
Emergencies within New Zealand	0800 764 766 (National Poison Control Centre)

Section 2: Hazards Identification

Product is classified as hazardous according to the HSNO (Minimum Degrees of Hazard) Regulations 2001

NZ Group Standard & EPA Approval Code Surface Coatings and Colourants (Subsidiary)
HSR002670



Irritant

Hazard Classification	Hazard Code	Hazard Statement
6.4A	H319	Causes serious eye irritation
6.5B	H317	May cause an allergic skin reaction

Prevention Code	Prevention Statement
P103	Read label before use.
P261	Avoid breathing dust/fume/gas/mist/vapours/spray.
P264	Wash hands thoroughly after handling.
P272	Contaminated work clothing should not be allowed out of the workplace.
P280	Wear protective gloves/protective clothing/eye protection/face protection.

Response Code	Response Statement
P302 + P352	IF ON SKIN: Wash with plenty of soap and water.
P305 + P351 + P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P321	Specific treatment see Section 4. First Aid Measures.
P333 + P313	If skin irritation or rash occurs: Get medical advice/attention.
P337 + P313	If eye irritation persists: Get medical advice/attention.
P363	Wash contaminated clothing before reuse.

Storage Code	Storage Statement
-	No storage statements.

Disposal Code	Disposal Statement
P501	Dispose of according to Local Regulations.

Section 3: Composition/Information on Ingredients

Chemical characterization	Mixtures.
Description	Mixture of substances listed below with non-hazardous additions.

Dangerous components:CAS: 3077-12-1
EINECS: 221-359-1

2,2'-[(4-methylphenyl)imino]bisethanol



Xn R22

2.5-10%

Additional information

For the wording of the listed risk phrases refer to section 16.

Section 4: First Aid Measures

Description of first aid measures:

General information	No special measures required.
After inhalation	Supply fresh air; consult doctor in case of complaints.
After skin contact	Rinse with water. Generally the product does not irritate the skin.
After eye contact	Rinse opened eye for several minutes under running water. If symptoms persist, consult a doctor.
After swallowing	Rinse out mouth and then drink plenty of water. If symptoms persist consult doctor.
Information for doctor:	
Most important symptoms and effects, both acute and delayed	No further relevant information available.
Indication of any immediate medical attention and special treatment needed	No further relevant information available.

Section 5: Fire Fighting Measures

Extinguishing media:

Suitable extinguishing agents	CO ₂ , powder or water spray. Fight larger fires with water spray or alcohol resistant foam.
Special hazards arising from the substance or mixture	No further relevant information available.
Advice for firefighters:	
Protective equipment	No special measures required.

Section 6: Accidental Release Measures

Personal precautions, protective equipment and emergency procedures	Not required.
Environmental precautions	Do not allow to enter sewers/surface or ground water.
Methods and material for containment and cleaning up	Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).
Reference to other sections	No dangerous substances are released. See Section 7 for information on safe handling. See Section 8 for information on personal protection equipment. See Section 13 for disposal information.

Section 7: Handling and Storage

Handling:

Precautions for safe handling	Only adequately trained personnel should handle this product. For use in dentistry only.
Information about fire - and explosion protection	No special measures required.
Conditions for safe storage, including any incompatibilities:	
Storage:	
Requirements to be met by storerooms and receptacles	Store only in the original receptacle.

Information about storage in one common storage facility	Not required.
Further information about storage conditions	Store in a cool place. Keep container tightly sealed. Protect from heat and direct sunlight.
Specific end use(s)	No further relevant information available.

Section 8: Exposure Controls/Personal Protection

Additional information about design of technical facilities	No further data; see item 7.
Control parameters:	
Ingredients with limit values that require monitoring at the workplace	The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.
Additional information	The lists valid during the making were used as basis.
Exposure controls:	
Personal protective equipment:	
General protective and hygienic measures	Usual hygienic measures for dental practice.
Respiratory protection	Not required.
Protection of hands	The glove material has to be impermeable and resistant to the product/the substance/the preparation. Due to missing tests no recommendation to the glove material can be given for the product/the preparation/the chemical mixture. Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation.
Material of gloves	The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.
Penetration time of glove material	The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.
Eye protection	Safety glasses.

Section 9: Physical and Chemical Properties

Information on basic physical and chemical properties:	
General Information:	
Appearance:	
Form	Fluid
Colour	Colourless
Odour	Characteristic
Odour threshold	Not determined
pH-value	Not determined
Change in condition:	
Melting point/Melting range	Undetermined
Boiling point/Boiling range	Undetermined
Flash point	Not applicable
Self-igniting	Product is not self igniting
Danger of explosion	Product does not present an explosion hazard
Explosion limits:	
Lower	Not determined

Upper	Not determined
Vapour pressure	Not determined
Density:	Not determined
Relative density	Not determined
Vapour density	Not determined
Evaporation rate	Not determined
Solubility in/Miscibility with water	Soluble
Segregation coefficient (n-octanol/water)	Not determined
Viscosity:	
Dynamic	Not determined
Kinematic	Not determined
Other information	No further relevant information available

Section 10: Stability and Reactivity

Reactivity	No further relevant information available.
Chemical stability	Stable under normal handling and storage conditions.
Thermal decomposition/ conditions to be avoided	No decomposition if used according to specifications.
Possibility of hazardous reactions	No dangerous reactions known.
Conditions to avoid	No further relevant information available.
Incompatible materials	No further relevant information available.
Hazardous decomposition products	None under normal conditions of storage and use.

Section 11: Toxicological Information

Information on toxicological effects:	
Acute toxicity:	
Primary irritant effect:	
On the skin	No irritant effect.
On the eye	No irritating effect.
Sensitization	No sensitizing effects known.
Additional toxicological information:	When used and handled according to specifications, the product does not have any harmful effects to our experience and the information provided to us.

Section 12: Ecological Information

Toxicity:	
Acquatic toxicity	No further relevant information available.
Persistence and degradability	No further relevant information available.
Behaviour in environmental systems:	
Bio accumulative potential	No further relevant information available.
Mobility in soil	No further relevant information available.
Additional ecological information:	
General notes	Water hazard class 1 (German Regulation) (Self-assessment): slightly hazardous for water. Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.

Results of PBT and vPvB assessment:

PBT	Not applicable.
vPvB	Not applicable.
Other adverse effects	No further relevant information available.

Section 13: Disposal Considerations

Waste treatment methods:

Recommendation Take to an approved landfill or a waste incineration plant, under conditions approved by the local authority.

Uncleaned packaging:

Recommendation Disposal must be made according to official regulations.

Section 14: Transport Information

This product is not a Dangerous Good for Road and Rail transport in New Zealand according to NZS 5433 Safe Transport of Dangerous Goods.

Section 15: Regulatory Information

NZ Group Standard & Surface Coatings and Colourants (Subsidiary)
EPA Approval Code HSR002670

HSNO Controls:

Trigger quantities for this substance:

	Trigger Quantity
Approved Handler	Not applicable
Location Certificate	Not applicable
Tracking Trigger Quantities	Not applicable
Signage Trigger Quantities	Not applicable
Emergency Response Plan Trigger Quantities	1,000 L (liquid) / kg (solid)

Section 16: Other Information

Relevant R-phrases: R22 Harmful if swallowed.

This document has been compiled by Ivoclar Vivadent Ltd on behalf of the manufacturer of the product and serves as the manufacturer's Safety Data Sheet ('SDS'). It is based on information concerning the product which has been provided to Ivoclar Vivadent (NZ) Ltd by the manufacturer or obtained from third party sources and is believed to represent the current state of knowledge as to the appropriate safety and handling precautions for the product at the time of issue. Further clarification regarding any aspect of the product should be obtained directly from the manufacturer. While Ivoclar Vivadent (NZ) Ltd has taken all care to include accurate and up-to-date information in this SDS, it does not provide any warranty as to accuracy or completeness. As far as lawfully possible, Ivoclar Vivadent (NZ) Ltd accepts no liability or any loss, injury or damage (including consequential loss) which may be suffered or incurred by any person as a consequence of their reliance on the information contained in this SDS. The information herein is given in good faith, but no warranty, express or implied is made.

Issue date: 26 August 2011

Date Ivoclar Vivadent Ltd prepared: 8 October 2014

Review date: 26 August 2016