

Programat® CS3

Software Revision Log

Version	Changes from previous version	Status
V3.1 01.02.2018	Software optimized <ul style="list-style-type: none"> • WLAN performance has been enhanced • Push function can be tested • Multimedia files can be directly deleted • Various languages have been improved • Various software optimizations 	Currently delivered version
V3.0 28.04.2017	Software optimized <ul style="list-style-type: none"> • New program structure with new program parameters • Various software optimizations 	No longer available!
V2.1 25.09.2016	Software optimized <ul style="list-style-type: none"> • Optimized text boxes for Japanese language • Various software optimizations 	No longer available!
V2.0 08.04.2016	Software optimized <ul style="list-style-type: none"> • DSA - Digital Shade Assistant: Orientation of the reference teeth either "incisal edge to incisal edge" or "incisal edge to tooth neck" possible • In-process dehumidification: Function that protects the vacuum system from contamination and its consequential damages • Various software optimizations 	No longer available!
V1.0 15.10.2015	Basic version	No longer available!

Version 01.02.2018