



# CAD / CAM BLOCK OVERVIEW FOR PrograMill

# IPS Empress® CAD

Leucite glass-ceramic | Leuzit-Glaskeramik


## Multi



C 14 L

## Shades | Farben

A-D



A1 A2 A3 A3.5 B1

Bleach



BL 1 BL 3

## HT (High Translucency)



I 12

A-D



A1 A2 A3 A3.5

B1 B2 B3 C2 D3

## LT (Low Translucency)



C 14

A-D



A1 A2 A3 A3.5

B1 B2 B3 C2 D3

Bleach



BL 1 BL 2

BL 3 BL 4

## Indications | Indikationen

Veneers, inlays, onlays, partial crowns, crowns

Veneers, Inlays, Onlays, Teilkronen, Einzelzahnkronen

# IPS e.max<sup>®</sup> CAD

Lithium disilicate glass-ceramic (LS<sub>2</sub>) | Lithium-Disilikat Glaskeramik (LS<sub>2</sub>)

## HT (High Translucency)




I 12



C 14

### Shades | Farben

A-D



Bleach




## MT (Medium Translucency)



C 14

A-D



Bleach




## LT (Low Translucency)



C 14

A-D



Bleach




C 16



B 32

A-D



Bleach



## MO (Medium Opacity)



C 14



## Impulse



C 14



## Indications | Indikationen

Thin veneers (0,4 mm), veneers, occlusal veneers, inlays, onlays, partial crowns, crowns, minimal invasive crowns (≥1 mm)\*\*, 3-unit-bridges (up to the second premolar as the terminal abutment)

Dünne Veneers (0,4 mm), Veneers, Okklusale Veneers, Inlays, Onlays, Teilkronen, Einzelzahnkronen, Minimalinvasive Kronen (≥1 mm)\*\*, 3-gliedrige Brücken (bis zum zweiten Prämolaren als endständigen Pfeiler)

\*\* Years of clinical experience have shown that the material's high strength in combination with adhesive cementation permit a crown layer thickness of at least 1 mm.

\*\* Die langjährige klinische Erfahrung und die hohe Festigkeit in Kombination mit der adhäsiven Befestigung erlauben eine Kronenschichtstärke von mind. 1 mm.

# IPS e.max<sup>®</sup> ZirCAD

Zirconium oxide ceramic (ZrO<sub>2</sub>) | Zirkoniumoxid-Keramik (ZrO<sub>2</sub>)

## LT (Low Translucency)



C 17



B 45

## Shades | Farben

A-D



A1 A2 A3 B1 B2 C2 D2

Bleach



BL

## Indications | Indikationen

Crowns, 3-unit-bridges

Einzelzahnkronen, 3-gliedrige Brücken

# Telio® CAD

Cross-linked PMMA material | Vernetztes PMMA-Material


## LT (Low Translucency)



B 40 L

## Shades | Farben

A-D



Bleach




---

## Indications | Indikationen

Temporary crowns and bridges  
Temporäre Kronen und Brücken

# Ivoclar Vivadent Block Dimensions Block-Dimensionen



| IPS Empress® CAD            |                 |                     |                |               |
|-----------------------------|-----------------|---------------------|----------------|---------------|
| Abbreviation  <br>Abkürzung | IPS Empress CAD | Size   Grösse in mm |                |               |
|                             |                 | Width   Breite      | Length   Länge | Height   Höhe |
| Inlay                       | I 12            | 10.0                | 12.0           | 15.0          |
| Crown                       | C 14            | 12.0                | 14.0           | 18.0          |
| Crown                       | C 14 L          | 14.4                | 14.4           | 18.0          |

| IPS e.max® CAD              |               |                     |                |               |
|-----------------------------|---------------|---------------------|----------------|---------------|
| Abbreviation  <br>Abkürzung | IPS e.max CAD | Size   Grösse in mm |                |               |
|                             |               | Width   Breite      | Length   Länge | Height   Höhe |
| Inlay                       | I 12          | 10.4                | 12.5           | 15.0          |
| Crown                       | C 14          | 12.4                | 14.5           | 18.0          |
| Crown                       | C 16          | 17.8                | 15.8           | 18.0          |
| Bridge                      | B 32          | 14.5                | 14.5           | 32.0          |

| IPS e.max® ZirCAD           |                  |                     |                |               |
|-----------------------------|------------------|---------------------|----------------|---------------|
| Abbreviation  <br>Abkürzung | IPS e.max ZirCAD | Size   Grösse in mm |                |               |
|                             |                  | Width   Breite      | Length   Länge | Height   Höhe |
| Crown                       | C 17             | 17.0                | 19.0           | 22.0          |
| Bridge                      | B 45             | 17.0                | 19.0           | 45.0          |

| Telio® CAD                  |           |                     |                |               |
|-----------------------------|-----------|---------------------|----------------|---------------|
| Abbreviation  <br>Abkürzung | Telio CAD | Size   Grösse in mm |                |               |
|                             |           | Width   Breite      | Length   Länge | Height   Höhe |
| Bridge                      | B 40 L    | 15.4                | 19.0           | 39.0          |
| Bridge                      | B 55      | 15.5                | 19.0           | 55.0          |

Date information prepared: 2017-09

[www.ivoclarvivadent.com](http://www.ivoclarvivadent.com)

Ivoclar Vivadent AG

Bendererstr. 2 | 9494 Schaan | Liechtenstein | Tel. +423 235 35 35 | Fax +423 235 33 60