



BioUniversal[®] PdF

High-gold universal alloy

for use with special ceramic and composite materials.

Au 71.1	Pt 9.2	Ag 11.7	Cu 4.5	In 1.3	Ir <1.0	Fe <1.0	Ta <1.0	Zn 1.5
-------------------	------------------	-------------------	------------------	------------------	-------------------	-------------------	-------------------	------------------

Advantages

- Palladium-free
- Esthetic, yellow color
- Compatible with special ceramic and composite materials
- Certified biocompatibility
- Wide range of indications

Indication

Onlays, $\frac{3}{4}$ crowns, PFM Crowns, crowns, telescope and conus crowns, posts, short- and long-span bridges

Technical data

Color	yellow
Type	4
Density (g/cm ³)	16.2
Melting range (°C)	910 – 970
Casting temperature (°C)	1080 – 1140
CTE 25 – 500 °C	15.8
Elongation (%)	9.0
Modulus of elasticity (MPa)	85.000
Oxide firing °C / min. / vacuum	700 / 3 / vacuum
Vickers hardness	210
0.2 % Proof stress (MPa)	450



Certificate

Test material: BioUniversal Alloys

Composition in % weight	Au	Pt	Pd	Ag	Cu	In	Ir	Zn	Other
BioUniversal® PKF	75.0	8.6	–	11.9	–	<1.0	<1.0	2.0	Fe <1.0, Ta <1.0, Rh <1.0
BioUniversal® PdF	71.1	9.2	–	11.7	4.5	1.3	<1.0	1.5	Fe <1.0, Ta <1.0
BioUniversal®	59.4	2.0	9.5	25.5	–	–	–	2.0	Fe <1.0, Re <1.0, Ru <1.0, Sn <1.0
BioUniversal® E	–	–	40.0	52.2	–	1.4	–	<1.0	Sn <5.5, Re <1.0, Ru <1.0

Manufacturer

Ivoclar Vivadent Inc., 175 Pineview Drive, Amherst, NY 14228, USA

Corrosion resistance

The test was conducted according to the international regulations of ISO 1562 and ISO 6871–1: static immersion test through analytical determination of the metal ion release after a 7-day immersion.

Test results: The metal ion release after 7 days of immersion was not significant.

Testing facility: Louisiana State University, Dr. Sakar

Cytotoxicity

The Agar Diffusion test determines the biological reactivity of cell culture on test material.

Test results: The test material is considered non-cytotoxic and meets the requirements of the Agar Diffusion test according to ISO 10993–5.

Mutagenicity

An Ames assay was conducted to determine any possible cancer potential.

Test results: No mutagenicity potential was found to exist in these BioUniversal alloys.

Kligman Maximization

This test evaluated the allergenic potential and/or sensitizing capacity of these BioUniversal alloys.

Test results: Based on the standards set by the study protocol, these alloys exhibited no reaction to the challenge (0 % sensitization).

Sensitivity of oral mucosa

Test to determine the contact sensitivity of these BioUniversal alloys at the buccal oral mucosa.

Test results: No reactions were noted in conjunction with these BioUniversal alloys.

Testing facility: Toxikon Corporation, 15 Wiggins Avenue, Bedford, Massachusetts

Amherst, May 2010



Dr. George Tysowsky, D. D. S., M. P. H.
Vice President-Technology