

Astropol[®]

Composite polishing

Lustrous results
in three steps




ivoclar
vivadent[®]
passion vision innovation

Astropol®

Composite polishing

Astropol® is a proven polishing system that provides professional polishing results in three steps. Astropol has been especially designed for polishing composite materials but it can also be used for polishing amalgam. Four different shapes allow a maximum of flexibility to ensure that all restoration areas can be conveniently accessed - for lustrous results.

Proven quality

- Smooth lustrous restorations with a precision finish
- 3 grit sizes for professional results
- Also suitable for polishing amalgam
- 4 shapes for flexible use in the anterior and posterior region

Professional polishing – in 3 steps

- Step 1: Pre-polishing (coarse - blue)
- Step 2: Polishing (fine - green)
- Step 3: High gloss polishing (superfine - pink)

Economical

- Autoclavable and reusable
- Cost-saving

Pre-polishing



Polishing



High gloss polishing



Dr Lukas Enggist, Ivoclar Vivadent, Schaan

Delivery forms:



This product forms part of the "Direct Restoratives" product category. Products from this category are optimally coordinated with each other.

Assortment

Polishing tips, coarse (F):	Polishing tips, fine (P):	Polishing tips, superfine (HP):
Small Flame (3 x)	Small Flame (2 x)	Small Flame (1 x)
Large Flame (3 x)	Large Flame (2 x)	Large Flame (1 x)
Cup (3 x)	Cup (2 x)	Cup (1 x)
Disc (3 x)	Disc (2 x)	Disc (1 x)

Refills

F & P: Each shape is available separately in packages of 6.
HP: Each shape is available separately in packages of 3.

Ivoclar Vivadent AG
Benderstr. 2
9494 Schaan
Liechtenstein
Tel. +423 235 35 35
Fax +423 235 33 60
www.ivoclarvivadent.com

Descriptions and data constitute no warranty of attributes. Printed in Liechtenstein
© Ivoclar Vivadent AG, Schaan Liechtenstein

678864/2015-06-18/EN

ivoclar
vivadent
passion vision innovation