

Programat® CS4

Software Revision Log

Version	Changes from previous version	Status
V2.0 04.02.2019	<p>Software optimized</p> <ul style="list-style-type: none"> • Programs for "IPS e.max ZirCAD MT Multi" added • Program order modified <ul style="list-style-type: none"> ○ PLEASE NOTE: Individual programs on program number 9 and 10 will be deleted in case of an update! • Various software improvements 	Currently delivered version
V1.1 26.07.2017	<p>Software optimized – new programs</p> <ul style="list-style-type: none"> • Parameters of program 3 "IPS e.max ZirCAD Stains/Glaze firing" have been optimized • Various software improvements 	No longer available!
V1.0 03.05.2017	Basic Version of the Programat CS4	No longer available!

Version 04.02.2019